

# Emergency Medical Equipment & Supplies Catalog

IMPROVING THE QUALITY OF PRE-HOSPITAL CARE





# How Can DiaMedical Help You?

You Save Lives, We Provide Equipment & Supplies.



## 01 Consolidate Sourcing to Reduce Costs

As a single source, we maintain vendor relationships so you don't have to. We guarantee trusted products at competitive prices.

## 02 Exclusive EMS Training & Simulation Products

As leaders in healthcare education, we provide products which elevate your simulations and training sessions.

## 03 Over 3,000 Warehoused Medical Products

Enjoy our fastest shipping ever on top emergency medical equipment, supplies and training products.

## 04 Complete Online Account Management

Instantly request quotes, access orders, and view account balances and contract pricing. By utilizing your favorites list you can simplify reordering disposable supplies.

Our Dedicated Team Is Ready to Serve You

(877) 593-6011 | [DiaMedical.com](http://DiaMedical.com) | [Info@DiaMedical.com](mailto:Info@DiaMedical.com)

## CONSOLIDATE PURCHASING

DiaMedical partners with top manufacturers to provide quality and trustworthy products that meet the expectations set forth by the EMS industry. Visit [DiaMedical.com](http://DiaMedical.com) for a complete list of all of the manufacturers we represent.





# TABLE OF CONTENTS

## 4 INFECTION CONTROL

- 5-7 Personal Protective Equipment
- 7-8 Gloves
- 8-9 Disinfectants & Sanitizers
- 9 Sharps Containers

## 10 AIRWAY & OXYGEN DELIVERY

- 11-13 ET Tubes
- 13-17 Airways
- 17-18 Laryngoscopes & Blades
- 19-21 CPR Barriers & Resuscitators
- 21 Oxygen Administration
- 21-23 Suction
- 23-24 Ventilators
- 25 Asthma Management

## 26 TRAUMA & WOUND CARE

- 27-29 Field Ready Loaded Jump Bags™
- 30-31 First Aid & Bleeding Control
- 32 Stop the Bleed Kits
- 33 OB Kits
- 33-34 Tools & Instruments
- 35-36 Gauze & Bandages
- 36 Tapes
- 37 Hot & Cold Packs
- 37-39 Mass Casualty Incident

## 42 PATIENT HANDLING & IMMOBILIZATION

- 43-47 Cots & Stair Chairs
- 48-49 Cot Mattresses
- 49-52 Cot Accessories
- 53-54 Backboards & Transfer Aids
- 55 Backboard Accessories
- 55 Blankets
- 56-57 Head & Neck Immobilizations
- 57 Splints

## 58 MONITORING & DIAGNOSTICS

- 59-61 AEDs
- 61-62 AED Trainers
- 62-64 Defibrillators
- 65 AED Accessories & Cabinets
- 66 Automated Chest Compression Systems
- 67 Electrodes
- 67 Glucose Monitoring
- 68 Blood Pressure Monitors
- 68-69 Stethoscopes
- 69 Thermometers
- 70 Pulse Oximeters
- 71 Diagnostic Bundles

## 72 IV & DRUG DELIVERY

- 73 Vein Finders
- 74 IV Start Kits & Supplies
- 75 IV Administration Sets
- 76 IV Catheters
- 76 Syringes & Needles
- 77 Infusion Pumps & Tubing
- 77 IO Access Systems

## 78 TRAINING & SIMULATION

- 79-83 Loaded Jump Bags™ for Training
- 84-85 SimRig® The Ambulance Trainer
- 86-87 SimRescue™ The Auto Extrication System
- 88-91 Patient Simulator Manikins
- 91-93 Rescue Manikins
- 93-95 Airway Trainers
- 95-96 Cricothyrotomy Trainers
- 96-98 CPR Manikins
- 98-100 IO Trainers
- 101-102 Wound Trainers
- 104-108 Moulage
- 109 IV Skills Trainers
- 109-111 Simulated IV Fluids & Medications

## 112 EMERGENCY BAGS

- 113-115 Field Ready Loaded Jump Bags™
- 116-118 Equipment Bags
- 118-119 Oxygen Bags
- 119 Hard Cases
- 120-121 Supply Modules
- 121 Bag Accessories

A healthcare worker is shown from the chest up, wearing extensive personal protective equipment (PPE). She is wearing a clear face shield with a blue top edge, a white N95-style respirator mask covering her nose and mouth, and purple nitrile gloves. She is also wearing a light blue disposable gown. Her hair is dark and curly. She is holding one of the purple gloves in front of her chest. The background is a white ambulance with red and yellow emergency lights. A blue semi-transparent box is overlaid on the left side of the image, containing the text "INFECTION CONTROL".

# INFECTION CONTROL



## Personal Protective Equipment



## DiaMedical Level 1 General Use 3-Ply Mask

- Lightweight, non-woven polypropylene material
- Elastic loops fit comfortably around ears
- Recommended for use by the general public

SKU	Description	Pkg	Price
COV012034	3-Ply Mask with Ear Loops; Blue	50/bx	\$31.99
COV012035	3-Ply Mask with Ear Loops; Blue	2000/cs	\$424.99



COV012034



## DiaMedical Level 3 Procedure 4-Ply Mask

- Meets ASTM F2100-19 Level 3 standards
- Latex-free, lightweight, non-woven polypropylene material
- Fluid-resistant and disposable

SKU	Description	Pkg	Price
COV012060	4-Ply Mask with Ear Loops; Blue	100/bx	\$42.99
COV012061	4-Ply Mask with Ear Loops; Blue	2000/cs	\$800.99



COV012060



## PoweCom® KN95 Masks

- Multi-layer filtration system of non-woven soft and breathable fibers
- 46% non-woven polypropylene, 28% melt-blown fabric, 26% ES hot air cotton
- Latex-free ear loops

SKU	Description	Pkg	Price
COV012019	N95 Flat Fold Mask	20/pk	\$53.99
COV012020	N95 Flat Fold Mask	200/cs	\$427.99



COV012019



## Disposable Face Shields

- Optically clear, anti-fog polyester
- Full length face shield
- Dimensions: 12.5" W x 8.5" H

SKU	Description	Pkg	Price
INF091054	Face Shield with Foam and Elastic Band	100/bx	\$188.99
INF091054-CS400	Face Shield with Foam and Elastic Band	400/cs	\$747.99



INF091054



### Dynarex® Safety Glasses

- Meets ANSI Z87.1-2010 standard
- Fits over most prescription eyeglasses
- Durable and sturdy construction

SKU	Description	Pkg	Price
INF100973-BLK	Black Frame with Clear/Uncoated Lens	1 ea	\$4.19



INF100973-BLK



### Medline Yellow Fluid Resistant Isolation Gowns

- Lightweight, fluid-repelling polypropylene material
- Available with classic neck and waist ties
- One size fits most

SKU	Description	Pkg	Price
SK020217	Regular/Large Size; Yellow	50/bx	\$60.99



SK020217

### AAMI Level 2 Isolation Gown

- Meets AAMI Level 2 standards
- 3-layer design of polypropylene non-woven fabric
- AATCC 42 - Water Resistance: Spray Impact Penetration Test
- AATCC 127 - Water Resistance: Hydrostatic Pressure Test
- Universal size

SKU	Description	Pkg	Price
COV012028	Neck Ties and Elastic Wrist Cuff; White	100/cs	\$213.99



COV012028



### DiaMedical Personal Protection Kit

**Kit includes:**

- One (1) Yellow Isolation Gown
- One (1) Pair of Shoe Covers
- One (1) Bouffant Cap
- One (1) Face Mask with Ear Loops
- One (1) Pair Non-Powdered Gloves

SKU	Description	Pkg	Price
IC028601	Disposable, Single Use PPE Kit	1 ea	\$4.99
IC028601-PK5	Disposable, Single Use PPE Kit	5/pk	\$24.99
IC028601-PK25	Disposable, Single Use PPE Kit	25/pk	\$99.95
IC028601-CS100	Disposable, Single Use PPE Kit	100/cs	\$485.00



IC028601

+ customer favorite

### Pocket Nurse® Wire Glove Box Holder

- Holds standard size glove boxes
- Mount on wall or under counter
- Mounting hardware included
- Vertical or horizontal position

SKU	Description	Pkg	Price
INF100914-SNGL	Pocket Nurse® Wire Glove Box Holder - Single	1 ea	\$12.49
INF100914-DBL	Pocket Nurse® Wire Glove Box Holder - Double	1 ea	\$12.49
INF100914-TRPL	Pocket Nurse® Wire Glove Box Holder - Triple	1 ea	\$17.99



INF100914-TRPL



## Clinton White Steel Glove Box Holder

- Powder coated wire steel construction
- All units mount horizontal or vertical
- Mounting hardware included
- Universal size fits most glove boxes

SKU	Description	Pkg	Price
IC042310	Holds One (1) Glove Box	1 ea	\$44.95
IC042311	Holds Two (2) Glove Boxes	1 ea	\$51.95
IC042312	Holds Three (3) Glove Boxes	1 ea	\$57.95
IC042313	Holds Four (4) Glove Boxes	1 ea	\$75.95



IC042310

## Gloves

### Nitrile Powder-Free Exam Gloves

- Glove Color: Blue
- Material: Nitrile
- Latex-Free, Powder-Free, Non-Sterile
- Textured Fingertips

SKU	Description	Pkg	Price
COV012040-SM	Small, Blue	100/bx	\$16.99
COV012040-MD	Medium, Blue	100/bx	\$16.99
COV012040-LG	Large, Blue	100/bx	\$16.99
COV012040-XL	X-Large, Blue	100/bx	\$16.99



COV012040-SM



### Pulse® Nitrile Exam Gloves

- Glove Color: Aqua Blue
- Material: Nitrile
- Latex-Free, Powder-Free, Non-Sterile
- Textured Fingertips

SKU	Description	Pkg	Price
INF101080-SM	Small, Aqua Blue	100/bx	\$51.99
INF101080-MED	Medium, Aqua Blue	100/bx	\$51.99
INF101080-LG	Large, Aqua Blue	100/bx	\$51.99
INF101080-XLG	X-Large, Aqua Blue	100/bx	\$51.99



INF101080-SM



### NitriDerm® Ultra Black Nitrile Exam Gloves

- Glove Color: Black
- Material: Nitrile
- Latex-Free, Powder-Free, Non-Sterile
- Textured
- Fentanyl Tested: Yes

SKU	Description	Pkg	Price
INF101081-SM	Small, Black	100/bx	\$31.99
INF101081-MED	Medium, Black	100/bx	\$31.99
INF101081-LG	Large, Black	100/bx	\$31.99
INF101081-XLG	X-Large, Black	100/bx	\$31.99



INF101081-SM

### Pocket Nurse® Nitrile Gloves

- Contains no natural rubber latex proteins
- Material: Nitrile
- Latex-Free, Powder-Free, Non-Sterile
- Ambidextrous

SKU	Description	Pkg	Price
INF100906-S	Pocket Nurse® Nitrile Gloves - Small	100/bx	\$14.99
INF100906-M	Pocket Nurse® Nitrile Gloves - Medium	100/bx	\$14.99
INF100906-L	Pocket Nurse® Nitrile Gloves - Large	100/bx	\$14.99
INF100906-XL	Pocket Nurse® Nitrile Gloves - X-Large	100/bx	\$14.99



INF100906-L



### NitriDerm® Ultra Blue Nitrile Exam Gloves

- Glove Color: Light Blue
- Material: Nitrile
- Latex-Free, Powder-Free, Non-Sterile
- Textured
- Fentanyl Tested: Yes

SKU	Description	Pkg	Price
INF101082-SM	Small, Light Blue	100/bx	\$25.99
INF101082-MED	Medium, Light Blue	100/bx	\$25.99
INF101082-LG	Large, Light Blue	100/bx	\$25.99
INF101082-XLG	X-Large, Light Blue	100/bx	\$25.99



INF101082



### Dynarex® Safe-Touch Blue Nitrile Exam Gloves

- Glove Color: Blue
- Material: Nitrile
- Latex-Free, Powder-Free, Non-Sterile

SKU	Description	Pkg	Price
INF101073-SM	Small, Blue	100/bx	\$24.99
INF101073-MED	Medium, Blue	100/bx	\$24.99
INF101073-LG	Large, Blue	100/bx	\$24.99
INF101073-XLG	X-Large, Blue	100/bx	\$24.99



INF101073

## Disinfectants & Sanitizers

### Sanidry Disinfectant Wipes

- Kills 99.99% of bacteria and viruses
- For use on most hard, non-porous surfaces and equipment
- Wipe dimensions: 5.9" x 6.7"



SKU	Description	Pkg	Price
COV012057	300 Wipes per Bucket	1 ea	\$18.99
COV012059	300 Wipes per Bucket	4/cs	\$69.99



### Purell® Advanced Hand Sanitizer Gel

- Kills 99.99% of most common germs that may cause illness
- Contains 70% ethyl alcohol
- Fragrance-free



SKU	Description	Pkg	Price
COV012038	64 oz. Pump Bottle	1 ea	\$34.95
COV012039	64 oz. Pump Bottle	4/cs	\$124.95



### CaviCide™ Surface Disinfectant Spray

- Effective against TB and MRSA in three (3) minutes
- Effective against HIV-1, HCV and HBV in two (2) minutes
- For use on most hard, non-porous surfaces



SKU	Description	Pkg	Price
INF100912-1G	1 Gallon	1 ea	\$45.25
INF100912-24oz	24 oz.; Spray Bottle	1 ea	\$19.99

### CleanCide® Disinfectant Wipes

- Virucidal, germicidal, fungicidal and tuberculocidal efficacy
- Effective against antibiotic resistant organisms
- For use on most hard, non-porous surfaces
- Wipe dimensions: 7" x 7"



SKU	Description	Pkg	Price
COV012026	Tub: 400 Wipes	1 ea	\$41.99
COV012027	Tub: 400 Wipes	4/cs	\$117.99

## Sharps Containers



### Covidien SharpSafety™ Transportable Sharps Container

- Simple, one-handed operation opens and closes disposal opening
- Dimensions: 8.75" H x 4.5" W x 2.5" D

SKU	Description	Pkg	Price
INF100971	1 qt; Red	1 ea	\$8.31



INF100971



### Curaplex® Sharps Solo

- Easy to use, one-handed operation
- Locking lid
- Puncture-resistant lid

SKU	Description	Pkg	Price
EM090225	Compact Size	1 ea	\$4.95
EM090226	Compact Size	24/cs	\$98.95



EM090225



### Dynarex® Tabletop Sharps Container

- Available in multiple sizes/fill capacities
- Transparent lid allows for easy visual monitoring to avoid overfilling
- Created with puncture-resistant plastic
- Locking lid

SKU	Description	Pkg	Price
NF101087-1GAL	1 Gallon Slide Lid	1 ea	\$8.29
INF101087-2GAL	2 Gallon Slide Lid	1 ea	\$9.49
INF101087-3GAL	3 Gallon Slide Lid	1 ea	\$18.99
INF101087-4GAL	4 Gallon Slide Lid	1 ea	\$24.99





**AIRWAY &  
OXYGEN  
DELIVERY**



**ET Tubes**



**Dynarex® Uncuffed Endotracheal Tubes**

- Uncuffed
- High-Volume, Low-Pressure
- Thermosensitive tubes are kink-resistant
- Features smooth Murphy eye holes to help protect vocal cords
- Kink-Resistant, Sterile, Latex-Free, Individually Wrapped

SKU	Description	Pkg	Price
RES091245-2.5MM	2.5 mm, Uncuffed	1 ea	\$2.49
RES091245-3.0MM	3.0 mm, Uncuffed	1 ea	\$2.49
RES091245-4.0MM	4.0 mm, Uncuffed	1 ea	\$2.49
RES091245-4.5MM	4.5 mm, Uncuffed	1 ea	\$2.49



RES091245

**Teleflex®**

**Teleflex® Pediatric Cuffed Endotracheal Tubes**

- Cuffed
- High-Volume, Low-Pressure
- Thermosensitive tubes are kink-resistant
- Features smooth Murphy eye holes to help protect vocal cords
- Kink-Resistant, Sterile, Latex-Free, Individually Wrapped

SKU	Description	Pkg	Price
RES091233-3.0MM	3.00 mm, Cuffed	1 ea	\$12.99
RES091233-4.0MM	4.00 mm, Cuffed	1 ea	\$12.99
RES091233-5.0MM	5.00 mm, Cuffed	1 ea	\$12.99
RES091239-6.0MM	6.00 mm, Cuffed	1 ea	\$4.19
RES091239-6.5MM	6.50 mm, Cuffed	1 ea	\$4.19
RES091239-7.0MM	7.00 mm, Cuffed	1 ea	\$4.19
RES091239-7.5MM	7.50 mm, Cuffed	1 ea	\$4.19
RES091239-8.0MM	8.00 mm, Cuffed	1 ea	\$4.19



RES091233-3.0MM



### Dynarex® Cuffed Endotracheal Tube With Preloaded Stylet

- Cuffed
- Smoothly crafted Murphy eye holes to help protect vocal cords
- Sterile
- Latex-free
- Pre-assembled and ready to use

SKU	Description	Pkg	Price
RES091246-5.0MM	5.0 mm; Cuffed with Stylet	1 ea	\$8.19
RES091246-5.5MM	5.5 mm; Cuffed with Stylet	1 ea	\$8.19
RES091246-6.0MM	6.0 mm; Cuffed with Stylet	1 ea	\$8.19
RES091246-6.5MM	6.5 mm; Cuffed with Stylet	1 ea	\$8.19
RES091246-7.0MM	7.0 mm; Cuffed with Stylet	1 ea	\$8.19
RES091246-7.5MM	7.5 mm; Cuffed with Stylet	1 ea	\$8.19



RES091246



### Dynarex® Uncuffed Endotracheal Tube With Preloaded Stylet

- Uncuffed
- Smoothly crafted Murphy eye holes to help protect vocal cords
- Pre-assembled and ready to use
- Sterile, Latex-Free, Individually Wrapped

SKU	Description	Pkg	Price
RES011057	2.0 mm; Uncuffed with Stylet	10/bx	\$15.99
RES011058	2.5 mm; Uncuffed with Stylet	10/bx	\$15.99
RES011059	3.0 mm; Uncuffed with Stylet	10/bx	\$15.99
RES011060	3.5 mm; Uncuffed with Stylet	10/bx	\$15.99
RES011061	4.0 mm; Uncuffed with Stylet	10/bx	\$15.99
RES011062	4.5 mm; Uncuffed with Stylet	10/bx	\$15.99
RES011063	5.0 mm; Uncuffed with Stylet	10/bx	\$15.99




RES011057



### Dynarex® Resp-O2™ Endotracheal Tube Introducer Bougies

- Designed to guide and support endotracheal tubes during difficult intubations and “blind” intubations
- Crafted from high-quality plastic materials for optimal durability
- Features clearly visible depth markings for improved safety
- Sterile, Latex-Free, Individually Wrapped

SKU	Description	Pkg	Price
RES012261	10 FR, Pediatric	10/bx	\$72.95
RS120209	15 FR, Adult	10/bx	\$79.95



### Laerdal® Thomas™ Tube ET Holder

- Hook-and-loop padded strap provides secure and comfortable fit, even when wet
- Large opening allow for view of patient's mouth and lip color
- Easy access for suctioning and other oral care needs
- Pliable bit block prevents patient from biting or pinching the tub and restricting airflow
- Non-Sterile, Latex-Free, Disposable

SKU	Description	Pkg	Price
RES012231-PED	Pediatric	1 ea	\$139.99
RES012231-ADLT	Adult	1 ea	\$139.99



RES012231

## Airways



### ADC® Berman Oral Airway Kit

- Kit includes six (6) sizes: 40 mm, 60 mm, 80 mm, 90 mm, 100 mm, 110 mm
- Smooth polyethylene airway with reinforced bite block
- Color coded for easy identification
- Non-Sterile, Latex-Free, Disposable

SKU	Description	Pkg	Price
RES011014	Disposable Berman Airway Kit; Boxed	6/bx	\$10.41



RES011014



### Dynarex® Resp-O2™ Berman Oral Airways

- Vented Septum
- Designed with firm, polyethylene plastic and built-in bite block
- Non-Sterile, Latex-Free, Individually Wrapped, Disposable

SKU	Description	Pkg	Price
RES091258-40MM	40 mm, Latex-Free	1 ea	\$1.69
RES091258-60MM	50 mm, Latex-Free	1 ea	\$1.69
RES091258-80MM	60 mm, Latex-Free	1 ea	\$1.69
RES091258-90MM	70 mm, Latex-Free	1 ea	\$1.69
RES091258-100MM	80 mm, Latex-Free	1 ea	\$1.69



RES091258-40MM





### Medline Disposable Guedel Airways

- Designed with smooth, curved edges for improved comfort
- Reinforced bite block provides optimal stability
- Latex-Free, Individually Wrapped, Disposable

SKU	Description	Pkg	Price
SK020102	50 mm, Infant; Blue	10/bx	\$8.99
SK020103	60 mm, Small Child; Black	10/bx	\$8.99
SK020104	70 mm, Child; White	10/bx	\$8.99
SK020105	80 mm, Adult Small; Green	10/bx	\$8.99
SK020106	90 mm, Adult Medium; Yellow	10/bx	\$8.99
SK020107	100 mm, Adult Large; Red	10/bx	\$8.99
EM090345	110 mm, Adult X-Large; Orange	10/bx	\$9.99
EM090346	120 mm, Adult XX-Large; Lavender	10/bx	\$10.99



### ADC® Guedel Oral Airway Kit

- Kit includes eight (8) sizes: 40 mm, 50 mm, 60 mm, 70 mm, 80 mm, 90 mm, 100 mm, 110 mm
- Smooth polyethylene airway with reinforced bite block
- Color coded for easy identification
- Non-Sterile, Latex-Free, Disposable

SKU	Description	Pkg	Price
EM090179	Disposable Guedel Airway Kit; Boxed	8/bx	\$24.99




EM090179



### Dynarex Nasopharyngeal Airway Kit

- Kit includes nasopharyngeal airways and lubricating jelly packets
- 6-Piece kit includes sizes: 22 FR, 24 FR, 26 FR, 28 FR, 30 FR, 32 FR
- 9-Piece kit includes sizes: 20 FR, 22 FR, 24 FR, 26 FR, 28 FR, 30 FR, 32 FR, 34 FR, 36 FR
- Soft, flexible tubing resists collapse and kinking
- Rounded tip and fixed flange design ensure safe, secure placement
- Sterile, Latex-Free

SKU	Description	Pkg	Price
RES011073	6 Nasopharyngeal Airways; 6 Lubricating Jelly Packets	6/bx	\$30.99
EM090180	9 Nasopharyngeal Airways; 9 Lubricating Jelly Packets	9/bx	\$33.99



RES011073

**Ambu® King LTS-D™ Disposable Laryngeal Tube Kit**

- Kit includes: King LTS-D disposable laryngeal tube, one (1) lubricant packet and one (1) syringe
- Seals in the esophagus and oropharynx to provide positive pressure inflation
- Drain tube for gastric and suction catheters
- Latex-Free, Single Port Inflation

SKU	Description	Pkg	Price
RS120112	Size 3.0	5/cs	\$299.99
RS120113	Size 4.0	5/cs	\$299.99
RS120114	Size 5.0	5/cs	\$299.99



**Teleflex® LMA® Unique™ Laryngeal Mask Airway with Cuff Pilot™ Technology**

- Features aperture bars, integrated cuff pressure monitoring and silicone cuff
- Cuff Pilot™ Technology integrates cuff pressure indicator that monitors cuff pressure
- Aperture bars are designed to prevent the blockage of airflow by the epiglottis
- Silicone cuff is soft, flexible and contours to the anatomy
- Sterile, Individually Wrapped, Disposable

SKU	Description	Pkg	Price
RS120153	Size 1.0, Neonatal/Infant; ≤ 5 kg	10/cs	\$268.99
SK020121	Size 1.5, Infant; 5-10 kg	10/cs	\$268.99
SK020122	Size 2.0, Infant/Small Child; 10-20 kg	10/cs	\$268.99
SK020123	Size 2.5, Small Child; 20-30 kg	10/cs	\$268.99
SK020124	Size 3.0, Large Child; 30-50 kg	10/cs	\$268.99
SK020125	Size 4.0, Small Adult; 50-70 kg	10/cs	\$268.99
SK020291	Size 5.0, Large Adult; 70-100 kg	10/cs	\$268.99



SK020122

**Intersurgical® i-gel Supraglottic Airways**

- Shape is contoured to perilaryngeal anatomy eliminating the need for cuff inflation
- Color coded for easy identification
- Individually Wrapped, Disposable

SKU	Description	Pkg	Price
RES012245	Size 1.0, Neonatal; Pink	10/cs	\$239.99
RS120103	Size 1.5, Infant; Blue	10/cs	\$239.99
RES012246	Size 2.0, Small Child; Gray	10/cs	\$239.99
RS120104	Size 2.5, Large Child; White	10/cs	\$239.99
RES012241	Size 3.0, Small Adult; Yellow	25/cs	\$595.99
RES012242	Size 4.0, Medium Adult; Green	25/cs	\$595.99
RES012243	Size 5.0, Large Adult; Orange	25/cs	\$595.99





### North American Rescue® CricKit® Tactical Cricothyotomy Kit

**Kit includes:**

- One (1) Nylon Bag
- One (1) Alcohol Packet
- One (1) Povidone Iodine Packet
- One (1) 2-Pack Gauze Pad - 4" x 4"
- One (1) Protected #10 Scalpel
- One (1) NAR Tracheal Hook
- One (1) Syringe - 10cc
- One (1) Cric. Tube with Obturator, Cuffed - 6.0 mm
- One (1) Tube Securing Strap
- One (1) Plastic Vacuum Bag

Item requires Medical Device Authorization for purchase

SKU	Description	Pkg	Price
RES012225	6.0 mm; Emergency Cricothyotomy Kit	1 ea	\$129.95



RES012225



### Curaplex® Cricothyotomy Field Kit

**Kit includes:**

- One (1) Murphy Endotracheal Tube, Cuffed - 6.0 mm
- One (1) Luer Lock Syringe - 12 mL
- One (1) Face Mask
- One (1) Safety Scalpel - #10
- One (1) Sterile Gauze Pad - 4" x 4"
- One (1) Endotracheal Tube Holder
- One (1) PVP Prep Pad
- One (1) Tracheal Dilator - 5.5"
- One (1) Stainless Steel Trach Hook

Item requires Medical Device Authorization for purchase

SKU	Description	Pkg	Price
EM090182	6.0 mm; Emergency Surgical Airway Kit	1 ea	\$53.99



EM090182



### Dynarex® DynaLube™ Lubricating Jelly

- Designed for medical purposes to lubricate orifices and assists in the insertion of devices for diagnostic/therapeutic purposes
- Provides lubrication for hinged instruments
- Color: Clear
- Sterile, Water-Soluble, Alcohol-Free, Non-Staining

SKU	Description	Pkg	Price
MED081055	2.7 g Packets	144/bx	\$9.79
PTC111526-2OZ	2.0 oz Tube	1 ea	\$2.59
PTC111526-4OZ	4.0 oz Tube	1 ea	\$2.99



PTC111526-4OZ





## Dynarex® Non-Sterile Tongue Depressors

- High-quality, splinter-free wood
- Non-sterile, disposable

SKU	Description	Pkg	Price
PTC111527-5.5IN	Non-Sterile Depressor, 5.5" Junior	500/bx	\$11.99
PTC111527-6IN	Non-Sterile Depressor, 6" Adult	500/bx	\$13.49
PTC111528-5.5IN	Non-Sterile Depressor, 5.5" Junior	100/bx	\$7.79
PTC111528-6IN	Non-Sterile Depressor, 6" Adult	500/bx	\$8.59



## Laryngoscopes & Blades



## MDPro EZvision Laryngoscope

- Equipped with state-of-the-art optics and illumination technology, the MDPro EZvision Laryngoscope ensures optimal visualization of the vocal cords and surrounding structures
- The high-resolution camera captures crisp, high-definition images, allowing healthcare providers to assess airway anatomy with precision and accuracy
- High-resolution camera for clear and precise visualization
- Versatile blade sizes for patients of all ages
- Compatible with standard video systems for seamless integration

SKU	Description	Pkg	Price
RES091290	EZVision Laryngoscope	Each	\$1,795.00
RES091291-2	Macintosh Blades - Size 2, Pediatric	10/box	\$115.00
RES091291-3	Macintosh Blades - Size 3, Adult	10/box	\$115.00
RES091291-4	Macintosh Blades - Size 4, Large Adult	10/box	\$115.00



## Disposable Fiberoptic Laryngoscope Handle

- Standard holding fixture for fiberoptic Miller and Macintosh blades

SKU	Description	Pkg	Price
RES091261	Disposable Laryngoscope Handle	1 ea	\$26.99
RES091229-2	Miller Blade #2	1 ea	\$8.39
RES091229-3	Miller Blade #3	1 ea	\$8.39
RES091229-4	Miller Blade #4	1 ea	\$8.39
RES091228-2	Macintosh Blade #2	1 ea	\$8.89
RES091228-3	Macintosh Blade #3	1 ea	\$8.89
RES091228-4	Macintosh Blade #4	1 ea	\$8.89



RES091261

### Zulco® Macintosh Laryngoscope Kit

**Set includes:**

- Battery-operated handle
- Sizes 2, 3 and 4 standard Macintosh blades
- Zippered carry case

SKU	Description	Pkg	Price
RES091217	Macintosh; Standard	1 ea	\$85.99



RES091217

### Zulco® Miller Laryngoscope Kit in Case

**Set includes:**

- Battery-operated handle
- Sizes 1, 2, 3 and 4 standard Miller blades
- Zippered carry case

SKU	Description	Pkg	Price
RES091218	Miller; Standard	1 ea	\$85.99



RES091218

simlabsolutions

### SimLabSolutions Rapid Sequence Intubation Training Bundle

**Training bundle includes:**

- Simulated Propofol/Diprivan
- Simulated Etomidate
- Simulated Midazolam
- Simulated Fentanyl
- Simulated Ketamine
- Simulated Succinylcholine
- Simulated Rocuronium
- Simulated Vecuronium
- Simulated Glycopyrrolate
- Oral Airways
- Disposable Laryngoscope Kit
- Cuffed Endotracheal Tubes
- Intubation Stylet
- ET Tube Holder
- Colorimetric CO2 Indicator
- Luer Lock Syringe, 10 mL
- Adult BVM Resuscitator
- Yankauer Suction with Tubing
- Open Suction Catheter Kit

SKU	Description	Pkg	Price
KB016003	RSI Training Bundles; For Training Use Only	1 ea	\$199.95



KB016003

For Training Use Only

## Maintain Competencies with Manikins & Trainers

Keep your staff up to date on emergency care and ACLS with low to high-fidelity manikins, CPR management trainers, airway simulators and more.

[SHOP AT DIAMEDICAL.COM](http://DIAMEDICAL.COM)



**CPR Barriers & Resuscitators**



**Dynarex® Resp-O2 CPR Face Shield with One Way Valve**

- Lightweight face shield to prevent mouth to mouth contact during CPR
- One way valve eliminates blow back to caregiver
- Completely seals to the faces of infants, children and adults



SKU	Description	Pkg	Price
RES091193	CPR Face Shield	1 ea	\$5.19



**Dynarex® Resp-O2™ CPR Rescue Mask Kit**

- Helps protect rescuers in an emergency situation
- Completely seals to the faces of infants, children and adults



SKU	Description	Pkg	Price
RES012115	CPR Rescue Mask Kit	25/bx	\$176.99



**Pocket Nurse® CPR Face Shield in Key Ring Case**

- Integrated filter protects rescuer from exhaled air
- Individually packed in nylon key ring pouches



SKU	Description	Pkg	Price
RES091151-BLK	Black	1 ea	\$3.79



**Dynarex® CPR Rescue Mask**

- Helps protect rescuers in an emergency situation
- Completely seals to the faces of infants, children and adults
- Latex-free



SKU	Description	Pkg	Price
RS120214	Mask and Gloves	25/bx	\$171.99



**Dynarex® Resp-O2™ Manual Resuscitator**

- 360° swivel port and intake valve to ensure one-way air flow
- Textured design for non-slip grip
- Translucent face mask
- Not made with natural rubber latex



SKU	Description	Pkg	Price
RES081007-E	Infant	1 ea	\$39.99
RES081004-E	Pediatric	1 ea	\$39.99
RES081006-E	Adult	1 ea	\$39.99



**Laerdal® The BAG II® Disposable Resuscitator**

- Includes mask, universal Star Lumen® Tubing and reservoir bag
- Elbow between mask and bag swivels 360°
- Latex-free, disposable



SKU	Description	Pkg	Price
RS044204	Infant #1	12/pk	\$239.00
RS044205	Child #3	12/pk	\$239.00
RS044206	Adult #5	12/pk	\$239.00



**MedSource Labs CPR Board**

- Designed for CPR training and administration
- Anatomically contoured to help lift and position the patient for effective CPR



SKU	Description	Pkg	Price
EM090422	CPR Board; Blue	1 ea	\$64.95



**Iron Duck® Bariboard Bariatric CPR Board**

- Raised surface is designed for bariatric patients
- Dual pane supports both the torso and spine
- Weight capacity: 275-500 lbs





SKU	Description	Pkg	Price
EM019306	Bariatric CPR Board; Blue	1 ea	\$160.60





# Repairs & Preventive Maintenance

DiaMedical offers Depot and On-Site Repairs for Movable Medical Equipment. Reduce wasted and unnecessary spend by partnering with our service network to provide cost-effective service for your entire facility. Our technicians work with facilities of all sizes to help optimize efficiency and increase profitability.



## Medtronic

### Nellcor™ Adult Colorimetric Easy Cap II CO2 Detector

- Designed to help verify endotracheal tube placement
- Quickly responds to exhaled CO2 with a color change from purple to yellow
- Breath-to-breath response for up to two (2) hours
- Latex-Free



SKU	Description	Pkg	Price
RES031047	Pediatric; 1-15 kg	1 ea	\$25.99



### Westmed CO2 Easy™ CO2 Detectors

- Carbon dioxide detector with large 2" visual indicator
- Anti-fog lens reduces condensation from exhaled gas
- Purple-to-yellow change indicates presence of exhaled CO2 and tracheal intubation
- Latex-Free




SKU	Description	Pkg	Price
RS120147	Pediatric; 1-15 kg	10/bx	\$209.95
RS120146	Adult; ≥ 15 kg	10/bx	\$159.95

## Oxygen Administration



### Dynarex® Resp-O2™ Nasal Oxygen Cannulas with Standard Tips

- Designed with non-flared nasal prongs
- Crush-resistant, StarLumen oxygen supply tubing
- Latex-Free, Individually Wrapped, Disposable



SKU	Description	Pkg	Price
RES031029	Adult; 7' Standard Connector	50/cs	\$34.95
RS120210	Pediatric; 7' Standard Connector	50/cs	\$39.95



### Dynarex® Resp-O2™ Oxygen Supply Tubing

- Crush-resistant, StarLumen oxygen supply tubing
- Latex-Free, Individually Wrapped, Disposable



SKU	Description	Pkg	Price
RES031033	7' Standard Connector	50/cs	\$19.95
RES031080	25' Standard Connector	25/cs	\$39.95



### Dynarex® Resp-O2™ Non-Rebreather Masks

- Transparent vinyl mask features smooth edges, an adjustable nose clip and elastic strap
- Includes a vinyl reservoir bag and 7' of crush-resistant, StarLumen oxygen supply tubing
- Non-Sterile, Latex-Free, Individually Packaged

SKU	Description	Pkg	Price
RES081007	Infant; Standard	6/cs	\$109.95
RES081004	Pediatric; Elongated	6/cs	\$109.95
RES081006	Adult; Elongated	6/cs	\$109.95



RES081006

## Suction



### Dynarex® Resp-O2™ Medium Concentration Oxygen Mask

- Transparent vinyl with smooth edges, latex-free
- Adjustable nose clip and elastic strap
- Pre-attached 7 ft oxygen supply tubing

SKU	Description	Pkg	Price
RES091118	Pediatric	1 ea	\$1.79
RES091168	Adult	1 ea	\$2.79



RES091118



### Dynarex® Resp-O2™ Single Suction Catheter Kit

- Kit includes pop-up basin and one (1) pair of exam gloves
- Soft, pliable catheter with beveled tip provides optimal patient comfort
- Vented cap with control valve allows for greater procedural control
- Catheter tip is designed with two (2) staggered eyelet openings for maximum suctioning efficiency
- Catheter features easy-to-read depth-number markings
- Sterile, Latex-Free, Individually Wrapped, Disposable



RES071015

SKU	Description	Pkg	Price
RES071011	6 FR, Pediatric; Light Green	50/cs	\$53.99
RES071012	8 FR, Pediatric; Light Blue	50/cs	\$53.99
RES071013	10 FR, Pediatric/Adult; Black	50/cs	\$53.99
RES071015	12 FR, Adult; White	50/cs	\$53.99
RES071017	14 FR, Adult; Dark Green	50/cs	\$53.99
RES071018	16 FR, Adult; Orange	50/cs	\$53.99
RS120212	18 FR, Adult; Red	50/cs	\$53.99



### Dynarex® Resp-O2™ Yankauer Suction Handle

- Blue-tinted, transparent material with slip-resistant handle
- Ribbed connector for secure attachment to suction tubing
- Sterile, Latex-Free, Individually Wrapped, Disposable



SKU	Description	Pkg	Price
RES071001	Bulb Tip, Vented	50/cs	\$34.99
RES071070	Bulb Tip, Non-Vented	50/cs	\$34.99
RES071002	Straight Tip, Vented	50/cs	\$34.99
RES071071	Straight Tip, Non-Vented	50/cs	\$34.99



### Dynarex® Resp-O2™ Yankauer Suction Handle with Tubing

- Blue-tinted, transparent material with slip-resistant handle
- Non-vented bulb tip to minimize tissue trauma
- Sterile, Latex-Free, Individually Wrapped, Disposable



SKU	Description	Pkg	Price
SK020222	Bulb Tip, Non-Vented with 1/4" x 6' Tubing	50/cs	\$106.99



### Laerdal® V-VAC™ Starter Kit

- Manual suction unit features large bore intake to resist clogging
- Vacuum range of 170-380 mmHg and peak flow rate of 70 LPM
- Kit includes handle, two (2) replacement cartridges and double male connector



SKU	Description	Pkg	Price
EM090399	Manual Suction Unit Kit	1 ea	\$179.00



### Laerdal® Compact Suction Unit® 4 (LCSU® 4)

- Small and lightweight unit features on-hand-grip design
- Vacuum range of 50-550 mmHg with an internal integrated filter system
- 45-minute battery run time
- Includes suction unit, 300 mL disposable canister, patient tubing, carry bag, AC/DC adapter charger and battery



SKU	Description	Pkg	Price
RS014202	300 mL Compact Suction Unit	1 ea	\$739.00



## Laerdal® Suction Unit (LSU)

- Portable, powerful and reliable suction unit designed for optimal preparation in emergency situations
- Vacuum range of 80-550 mmHg and flow rate of >25 LPM
- 30-minute battery run time
- Includes suction unit, disposable canister, patient tubing, vacuum tubing, US power cord and battery

SKU	Description	Pkg	Price
RS014201	LSU with Disposable Bemis Canister	1 ea	\$1,199



RS014201



## Drive DeVilbiss® 7305 Series Vacu-Aide® Suction Unit

- Small, lightweight unit weighs only 3.8 lbs without battery
- Vacuum range of 80-550 mmHg and flow rate of 27 LPM
- Includes suction unit, 800cc disposable canister, tubing, AC/DC power cord

SKU	Description	Pkg	Price
EM090187	Unit Only	1 ea	\$244.95
EM090186	Unit with Battery	1 ea	\$374.95
EM090188	Unit with Battery & External Filter	1 ea	\$369.95



EM090188



## SCOR S-SCORT® III Portable Suction Unit

- Ultra-reliable and durable suction device weighs 7 lbs
- Two (2) position regulator adjusts for adult and tracheal/pediatric procedures
- Vacuum rate of 525 mmHg or 120 mmHg (+/- 15%) and flow rate of >30 LPM
- 45-minute battery run time
- Includes suction unit, 1200cc canister, patient tubing, catheter, charger and DC power cord

SKU	Description	Pkg	Price
EM090189	Suction Unit	1 ea	\$829.95



EM090189

## Ventilators



### Vyaire CareFusion LTV™ Series 1200 Ventilator

- Offers invasive and non-invasive modes of ventilation in pressure control, pressure support, volume control and spontaneous breath types
- Choose infant, pediatric or adult patient (≥ 5 kg) pre-sets for quick initial ventilation configuration
- Turbine technology eliminates the need for a high-pressure air source

SKU	Description	Pkg	Price
RS120194	New; Portable Ventilator	1 ea	\$16,995.00
RS020503	Reconditioned; Portable Ventilator	1 ea	\$7,995.00

#### PARTS & ACCESSORIES

EM090190	Vyaire CareFusion SprintPack™ Lithium-Ion Power System for LTV™	1 ea	\$1,559.95
RS120176	Teleflex® Dual-Limb Adult Ventilator Circuits	10/bx	\$96.95



RS120194





### Vyaire CareFusion ReVel™ Ventilator

- High-performance ventilation for pediatric or adult patients weighing ≥ 5 kg
- Includes integrated pulse oximeter with heart rate monitoring
- Offers invasive and non-invasive modes of ventilation
- User interface is based on the LTV series platform, easing use and reducing training time

SKU	Description	Pkg	Price
RS120196	New; Transport Ventilator	1 ea	\$24,995.00
RS120196-R	Reconditioned; Transport Ventilator	1 ea	\$16,995.00

#### PARTS & ACCESSORIES

RS120185	Vyaire Patient Circuit without PEEP with Elbow, 22mm SPU	10/bx	\$209.95
RS120187	Vyaire Patient Circuit without PEEP with Elbow, 15mm SPU	10/bx	\$214.95



RS120196



### Smiths Medical Pneupac VR1® Automatic Resuscitator

- Gas-powered, hand-held control unit
- Auto/manual control allows the operator to choose automatic mode or to provide ventilation manually
- Single knob for tidal volume/frequency control enables quick set-up
- Flow range 11-32 LPM; Pressure relief valve 40 cm H<sub>2</sub>O standard

SKU	Description	Pkg	Price
RES012016	VR1 Ventilator/Resuscitator	1 ea	\$3,274.95
RES012017	VR1 Ventilator/Resuscitator with Air Mix	1 ea	\$3,474.95

#### PARTS & ACCESSORIES

RES012018	Smiths Medical Pneupac Single Limb Circuit without PEEP Valve	10/bx	\$241.99
RES012021	Smiths Medical Pneupac Single Limb Patient Circuit with PEEP Valve	10/bx	\$289.99
RES012020	Smiths Medical Pneupac CPAP Circuit with Medium Adult Mask	6/bx	\$359.99
RES012019	Smiths Medical Pneupac CPAP Circuit with Large Adult Mask	6/bx	\$359.99



RES012016



### Dynarex® Resp-O2™ Compressor Nebulizer

- Compact design with rounded corners and built-in nebulizer holder
- Includes nebulizer, air tube, mouthpiece with T-adapter, adult mask, pediatric mask and five (5) disposable filters

SKU	Description	Pkg	Price
MED053253	Compact Nebulizer	1 ea	\$35.99



RES041004

### Hudson Nebulizer Kits & Masks

- Nebulizer cup features anti-spill design which allows patient to receive treatment at an angle
- Nebulizer Cup Capacity: 6cc
- Connection Type: Standard
- Latex-free

SKU	Description	Pkg	Price
RES091107	Nebulizer Kit T-Piece/Mouthpiece	1 ea	\$3.29
RES091237-ADLT	Adult Mask	1 ea	\$5.39
RES091237-PED	Pediatric Mask	1 ea	\$5.39



RES091107

### Hudson RCI® AsthmaMD™ Peak Flow Meter

- 60-810 LPM range for use with adults and children
- Red, yellow and green markers for use with asthma action plans

SKU	Description	Pkg	Price
RES091110	Peak Flow Meter; 60-800 mL	1 ea	\$29.99



RES091110





**TRAUMA &  
WOUND  
CARE**



# Field Ready Loaded BLS Jump Bag™

Save valuable time with the Field Ready Loaded BLS Jump Bag™. This jump bag comes loaded with first aid supplies and diagnostic equipment required to perform basic life support.



SKU	Description	Pkg	Price
LC129912-FR	DiaMedical Field Ready Loaded BLS Jump Bag™	1 ea	\$325.99
LC129922-FR	DiaMedical Field Ready Loaded BLS Jump Bag™ Refill Kit	1 ea	\$269.99
LC129912	SimLabSolutions Loaded BLS Jump Bag™; For Training Use Only	1 ea	\$325.99



### PPE/Infection Control

- Two (2) Pairs Sterile Exam Gloves, Vinyl
- Nine (9) Alcohol Gel Hand Wipes



### Diagnostic

- One (1) Glucose Meter
- Fifty (50) Test Strips
- Fifty (50) Safety Lancets
- One (1) Aneroid BP Cuff, Adult
- One (1) Stethoscope
- One (1) Penlight with Pupil Gauge
- One (1) Digital Thermometer
- One (1) Pulse Oximeter, Fingertip



### Respiratory

- One (1) BVM, Adult
- One (1) Disposable CPR Face Shield
- One (1) Each Disposable Guedel Airway - 50 mm, 60 mm, 70 mm, 80 mm, 90 mm & 100 mm
- One (1) Each Rusch Nasal Airway - 16 FR, 26 FR & 28 FR
- Six (6) Sterile Lubricating Jelly Packets
- One (1) Bite Stick



### Wound Care/Dressings

- One Hundred (100) Plastic Adhesive Strips - 0.75" x 3"
- One (1) Elastic Bandage - 3" x 5 yd
- Two (2) Sterile Trauma Dressings - 12" x 30"
- Two (2) Sterile Burn Dressings
- Two (2) Triangular Bandages
- Five (5) Sterile ABD Pads - 5" x 9"
- Two (2) Sterile Transparent Dressings - 4" x 4"
- Fifteen (15) Sterile Tongue Depressors
- Two (2) Sterile Sof-Form Gauze Rolls - 1"
- Three (3) Sterile Sof-Form Gauze Rolls - 3"
- Three (3) Sterile Gauze - 4" x 4"
- Four (4) Sterile Eye Pads
- Two (2) Transparent Adhesive Tape Rolls - 1"



### Additional Supplies

- One (1) Splinter Forceps
- One (1) Utility Scissors - 7.5"
- One (1) Kelly Forceps, Straight - 5.5"
- One (1) Instant Cold Pack - 5" x 9"
- One (1) Obstetrical Kit
- One (1) Emergency Survival Blanket
- Twenty-Five (25) Alcohol Prep Pads
- One (1) Each Cervical Collar - Adult & Pediatric





# Field Ready Loaded ALS Jump Bag™

Stay life support ready with the Field Ready Loaded ALS Jump Bag™! This loaded jump bags arrives stocked with airway supplies, diagnostic equipment, IV supplies and first aid supplies required for advanced life support.



SKU	Description	Pkg	Price
LC129911-FR	DiaMedical Field Ready Loaded ALS Jump Bag™	1 ea	\$745.99
LC129921-FR	DiaMedical Field Ready Loaded ALS Jump Bag™ Refill Kit	1 ea	\$595.99
LC129911	SimLabSolutions Loaded ALS Jump Bag™; For Training Use Only	1 ea	\$745.99

**PPE/Infection Control**

- Two (2) Pairs Sterile Exam Gloves, Vinyl
- Nine (9) Alcohol Gel Hand Wipes
- One (1) Procedure Mask with Face Shield

**Diagnostic**

- One (1) Glucose Meter
- Fifty (50) Test Strips
- Fifty (50) Safety Lancets
- One (1) Aneroid BP Cuff, Adult
- One (1) Stethoscope
- One (1) Penlight with Pupil Gauge
- One (1) Digital Thermometer
- One (1) Pulse Oximeter, Fingertip

**Infusion/Injection**

- Twenty-Five (25) Alcohol Prep Pads
- Three (3) Latex-Free Tourniquets, Rolled - 1" x 18"
- One (1) IV Armboard, Adult - 3" x 9"
- Two (2) Each Luer Lock Syringes - 3 mL, 5 mL & 10 mL
- Two (2) Tuberculin Safety Syringes
- Four (4) Each Safety Hypodermic Needles - 25G & 18G
- Four (4) IV Start Kits with ChloroPrep®
- One (1) IV Extension Set - 7"
- One (1) 3-Way Stopcock
- Two (2) Primary IV Sets - Microdrip & Macro drip
- Twenty (20) Alcohol Swabsticks
- One (1) Each Safety IV Catheters - 24G, 22G, 18G & 14G

**Additional Supplies**

- Two (2) Transparent Adhesive Tape Rolls - 1"
- One (1) Splinter Forceps
- One (1) Utility Scissors - 7.5"
- One (1) Kelly Forceps, Straight - 5.5"
- One (1) Instant Cold Pack - 5" x 9"
- One (1) Obstetrical Kit
- One (1) Emergency Survival Blanket
- One (1) Each Cervical Collars - Adult & Pediatric

**Respiratory**

- One (1) BVM, Adult
- One (1) Intubation Bag
- One (1) Emergency Cricothyrotomy Kit
- One (1) Each Nasal Cannulas - Adult & Pediatric
- One (1) Each Non-Rebreather Oxygen Masks - Adult & Pediatric
- Six (6) Sterile Lubricating Jelly Packets
- One (1) Bite Stick
- Five (5) Sterile Tongue Depressors
- One (1) Each Disposable Guedel Airway - 50 mm, 60 mm, 70 mm, 80 mm, 90 mm & 100 mm
- One (1) Each Rusch Nasal Airway - 16 FR, 26 FR, 28 FR & 30 FR
- One (1) Disposable Laryngoscope Handle
- One (1) Each Disposable Laryngoscope Blades, Miller - #0, #1 & #2
- One (1) Each Disposable Laryngoscope Blades, Macintosh - #2, #3 & #4
- One (1) Each Cuffed Endotracheal Tubes - 3.0 mm, 4.0 mm, 5.0 mm, 6.0 mm, 7.0 mm, 8.0 mm & 9.0 mm
- One (1) Each Intubation Stylets - 10 FR & 14 FR
- One (1) CO2 Indicator, Colorimetric
- One (1) ET Tube Holder
- One (1) Magill Catheter Forceps, Adult - 9.75"

**Wound Care/Dressings**

- One Hundred (100) Plastic Adhesive Strips - 0.75" x 3"
- Two (2) Sterile Trauma Dressings - 12" x 30"
- Two (2) Sterile Sof-Form Gauze Rolls - 1"
- Three (3) Sterile Sof-Form Gauze Rolls - 3"
- Three (3) Sterile Gauze - 4" x 4"
- Five (5) Sterile Gauze - 2" x 2"
- Two (2) Sterile Burn Dressings
- Five (5) Sterile ABD Pads - 5" x 9"
- Four (4) Sterile Eye Pads
- Two (2) Sterile Transparent Dressings - 4" x 4"
- One (1) Elastic Bandage - 3" x 5 yd
- Two (2) Triangular Bandages



# Field Ready Loaded Pediatric Jump Bag™

Pediatric emergencies are stressful enough, let DiaMedical make sure your jump bag is ready! The Field Ready Loaded Pediatric Jump Bag™ arrives with pediatric emergency equipment and supplies organized in Broselow Color-Coded modules for simple identification.



SKU	Description	Pkg	Price
LC129913-FR	DiaMedical Field Ready Loaded Pediatric Jump Bag™	1 ea	\$755.99
LC129923-FR	DiaMedical Field Ready Loaded Pediatric Jump Bag™ Refill Kit	1 ea	\$605.99
LC129913	SimLabSolutions Loaded Pediatric Jump Bag™; For Training Use Only	1 ea	\$755.99



## Color-Coded Weight Based Modules

- One (1) Suction Catheter Kit
- One (1) Cuffed Murphy Endotracheal Tube
- One (1) Disposable Laryngoscope Blade, Miller and/or Macintosh
- One (1) Intubating Stylet
- One (1) Disposable Guedel Airway
- One (1) Rusch Nasal Airway
- One (1) Safety IV Catheter
- One (1) Butterfly Infusion Set
- Two (2) IV Start Kits
- Ten (10) Alcohol Prep Pads
- Three (3) Sterile Lubricating Jelly Packets



## Additional Supplies

- One (1) Each Disposable IV Armboards - Neonatal, Infant & Adult
- Two (2) Each Luer Lock Syringes - 10 mL & 3 mL
- Four (4) Pediatric/Neonatal ECG Electrodes
- One (1) Each BVMs - Child & Infant
- One (1) Each Aneroid BP Cuffs - Pediatric & Infant
- One (1) Pediatric Stethoscope
- Two (2) Sterile Bulb Syringes - 60 mL
- One (1) Obstetrical Kit
- One (1) PediaQuick Tape
- One (1) Each Primary IV Sets - Microdrip & Macro drip
- One (1) IV Extension Set - 7"
- One (1) Disposable Laryngoscope Handle
- One (1) Whole Blood Tube Plastic Vacutainer
- One (1) 3-Way Stopcock
- One (1) Magill Catheter Forceps, Child
- One (1) Each Nasal Cannulas - Pediatric & Infant
- One (1) Non-Rebreather Oxygen Mask - Pediatric
- One (1) CO2 Indicator, Colorimetric

## TUESDAY TEACHINGS

### Bring Our Products Into Your Practice

Learn why our EMS Loaded Jump Bags™ are an essential tool for life-saving procedures & training scenarios!

WATCH NOW



First Aid & Bleeding Control

**Graham-Field First Aid Kit**

- Ideal for workplaces with up to twenty-five (25) people
- Kit comes in a water-resistant plastic case with dividers and a carry handle
- Carrying case may be wall-mounted



PTC111644

SKU	Description	Pkg	Price
PTC111644	First Aid Kit for 25 People	1 kit	\$55.99



**North American Rescue® S.T.O.R.M. Operator IFAK**

Kit Includes:

- One (1) Bag (STORM)
- Two (2) Pairs of Black Talon® Nitrile Trauma Gloves
- One (1) Nasopharyngeal Airway 28 FR with Lubricant
- One (1) HyFin® Vent Chest Seal Twin Pack
- One (1) ARS® Needle Decompression Kit - 14G x 3.25"
- One (1) C-A-T® (Combat Application Tourniquet®)
- One (1) ETD™ 6" Emergency Trauma Dressing
- Two (2) S-Rolled Gauze™ - 4.5" x 4.1 yd
- One (1) Polycarbonate Eye Shield
- One (1) Combat Casualty Card

Item requires Medical Device Authorization for purchase



EM090201

SKU	Description	Pkg	Price
EM090201	Black	1 ea	\$165.95



**North American Rescue® C-A-T (Combat Application Tourniquet), Gen 7**

- Performance compression band minimizes pressure loss over time after application
- Slack indicator wedge provides visual confirmation that you have removed all slack from the band
- Rugged buckle is designed to remove excess slack in a single fluid motion



MED171008

SKU	Description	Pkg	Price
MED171008	Black	1 ea	\$31.59
EM090280	Training Blue; For Training Use Only	1 ea	\$31.59



EM090280

For Training Use Only



### TacMed Solutions™ SOF® Tourniquet

- Performance compression band minimizes pressure loss over time after application
- Slack indicator wedge provides visual confirmation that you have removed all slack from the band
- Rugged buckle is designed to remove excess slack in a single fluid motion

SKU	Description	Pkg	Price
MED121041	Black	1 ea	\$34.00
MED121042	Rescue Orange	1 ea	\$34.00
EM090416	Training Blue; For Training Use Only	1 ea	\$34.00



For Training Use Only



### RevMedx™ TX Tourniquet

- Self-locking, non-slip ratcheting design with mechanical release
- High-density polymer ladder maintains function in presence of sand and dirt
- UV-resistant MILSPEC webbing
- Recommended by CoTCCC

SKU	Description	Pkg	Price
EM090197	TX2; 2" x 39"	1 ea	\$41.99
EM090198	TX3; 3" x 60"	1 ea	\$42.99



### SAM® XT Tourniquet

- Designed for rapid application by normalizing the number of windlass turns needed
- TRUFORCE™ Buckle technology auto-locks at a predetermined amount of circumferential force
- Recommended by CoTCCC

SKU	Description	Pkg	Price
EM090199	High Visibility Orange	1 ea	\$39.99
EM090417	High Visibility Blue	1 ea	\$39.99
EM090418	Tactical Black	1 ea	\$34.99



### Safeguard Medical™ SWAT-T™ Tourniquet

- Unique, multipurpose tourniquet can be utilized in three (3) applications: tourniquet, pressure dressing and elastic bandage
- Named after a description of usage - Stretcher, Wrap and Tuck

SKU	Description	Pkg	Price
EM090200	Rescue Orange	1 ea	\$18.85
EM090419	Tactical Black	1 ea	\$18.85
EM090420	Training Blue; For Training Use Only	1 ea	\$18.85



For Training Use Only



Stop the Bleed Kits



# North American Rescue® Public Access Bleeding Control

- Designed to provide the general public with essential equipment to stop life threatening bleeding
- Includes easy-to-use tools empowering non-medical professionals to take action if an everyday emergency occurs
- Intended to treat one (1) individual

SKU	Description	Pkg	Price
EM090153	Basic	1 ea	\$74.95
EM090154	Intermediate	1 ea	\$86.95
EM090156	Advanced	1 ea	\$139.95



Contents of **ADVANCED** Individual Bleeding Control Kit



EM090153



EM090154



EM090156

### Basic Kit

- One (1) C-A-T® Tourniquet
- One (1) Emergency Trauma Dressing - 6"
- Two (2) Wound Packing Gauze
- Two (2) Pairs Nitrile Gloves, Large
- One (1) Trauma Shears - 7.25"
- One (1) Survival Blanket
- One (1) Permanent Marker
- One (1) Instruction Card
- One (1) Red Nylon Bag with PVC Bleeding Control Patch

### Intermediate Kit

- One (1) HyFin® Vent Compact Chest Seal Twin Pack
- One (1) C-A-T® Tourniquet
- One (1) Emergency Trauma Dressing - 6"
- Two (2) Wound Packing Gauze
- Two (2) Pairs Nitrile Gloves, Large
- One (1) Trauma Shears - 7.25"
- One (1) Survival Blanket
- One (1) Permanent Marker
- One (1) Instruction Card
- One (1) Red Nylon Bag with PVC Bleeding Control Patch

### Advanced Kit

- One (1) Hemostatic Gauze
- One (1) HyFin® Vent Chest Seal, Twin Pack
- One (1) C-A-T® Tourniquet
- One (1) Emergency Trauma Dressing - 6"
- Two (2) Wound Packing Gauze
- Two (2) Pairs Nitrile Gloves, Large
- One (1) Trauma Shears - 7.25"
- One (1) Survival Blanket
- One (1) Permanent Marker
- One (1) Instruction Card
- One (1) Red Nylon Bag with PVC Bleeding Control Patch

OB Kits



**Dynarex Obstetrical Kit**

**Kit includes:**

- One (1) Pair Sterile Nitrile Exam Gloves
- One (1) Drape Sheet
- Six (6) Sterile Gauze Sponges
- One (1) Plastic Bag to Hold Placenta
- One (1) Underpad
- One (1) Receiving Blanket
- One (1) Sterile Bulb Syringe
- Two (2) OB Towelettes
- One (1) Sterile Disposable Scalpel
- Four (4) Disposable Towels
- Two (2) Sterile Umbilical Clamps
- Two (2) Twist Ties
- One (1) Sterile OB Pad
- Two (2) Nylon Tie-Offs
- Two (2) Alcohol Prep Pads, Medium



SK020238

SKU	Description	Pkg	Price
SK020238-E	Disposable OB Kit; Bagged	1 ea	\$10.49



**Medline OB Kit**

- Provides supplies required for pre-hospital vaginal deliveries
- Sterile, Disposable, Latex-Free



EM090202

SKU	Description	Pkg	Price
EM090202	Disposable OB Kit; Bagged	1 ea	\$18.99
EM090203	Disposable OB Kit; Bagged	20/cs	\$353.99

Tools & Instruments



**North American Rescue® Trauma Shears**

- Premium stainless steel blade with enhanced serrated lower edge
- Clothing-lift lip and rounded tip
- Soft-grip, non-slip contoured handle
- Includes an O2 wrench



SKU	Description	Pkg	Price
EM090421	6.25"; Black	1 ea	\$12.99
EM090204	7.25"; Black	1 ea	\$12.99



**North American Rescue® Responder Shears**

- Premium-grade trauma shears
- Designed for safety and durability



SKU	Description	Pkg	Price
EM090205	7.25"	1 ea	\$3.99



### ADC® Medicut™ EMT Shears

- Stainless steel blade with serrated edge
- Safety bandage tip and contoured polypropylene handles
- Blades tempered to a Rockwell hardness of C56
- Autoclavable
- Additional colors available online:



SKU	Description	Pkg	Price
SK020266-TCT	7.25"; Tactical Black	1 ea	\$5.69
SK020266-BLK	7.25"; Black	1 ea	\$5.69



### ADC® Magill Catheter Forceps

- Prema Traditional Magill Catheter Forceps offset design with serrated tips
- Forged from corrosion resistant 410 surgical stainless steel
- Designed for safety and durability
- Autoclavable, Floor Grade



SKU	Description	Pkg	Price
SK020259-BX	Child, 8"	12/bx	\$156.95
SK020260-BX	Adult, 9.75"	12/bx	\$156.95



### EMI Window Punch with Pocket Clip

- Designed to shatter tempered auto window glass
- Adjustable spring loaded for light or heavy punch
- Lightweight aluminum barrel with knurled finish and replaceable, steel tip
- Steel tip is hardened to Rockwell hardness of C60



SKU	Description	Pkg	Price
EM090207	0.625" x 5"	1 ea	\$16.95



### ADC® Plain Splinter Forceps

- Traditional Splinter forceps design with serrated fine tips and non slip grooved handle design
- Forged from corrosion resistant 410 surgical stainless steel
- Autoclavable, Floor Grade



SKU	Description	Pkg	Price
SK020264	3.5"	1 ea	\$1.99
SK020263	4.5"	1 ea	\$2.95



### Medline Ring Cutter

- Specialized emergency tool designed to cut rings from fingers
- Surgical grade, stainless steel ring cutter



SKU	Description	Pkg	Price
EM090206	Ring Cutter with Blade	1 ea	\$19.99



### EMI Lifesaver Plus™ Seat Belt Cutter with O2 Wrench

- Created to cut through restraining belts, clothes and seat belts
- Easy to grip, high-impact plastic with surgical grade, replaceable stainless steel blade



SKU	Description	Pkg	Price
EM090208	6.75" x 1.25"	1 ea	\$11.95

Gauze & Bandages



**Dynarex® Non-Stick Gauze Pads**

- Highly absorbent, easy to remove gauze pad will not stick to wounds
- Individually packaged and packed in convenient dispenser boxes
- Sterile, Non-Stick, Latex-Free



SKU	Description	Pkg	Price
PTC111185-2IN	2" x 3"	100/bx	\$7.79
PTC111185-3IN	3" x 4"	100/bx	\$11.49
PTC111261	3" x 8"	1 ea	\$0.49



**Dynarex® Surgical Gauze Sponge**

- Highly absorbent, 100% cotton gauze, non-woven
- Non-Sterile pads ideal for general cleaning



SKU	Description	Pkg	Price
MED900008	2" x 2"; 8-Ply	200/bx	\$1.99
MED900011	4" x 4"; 8-Ply	200/bx	\$5.79



**Dynarex® Conforming Stretch Gauze Bandage**

- Knit construction provides compression without constriction
- Non-Sterile, Latex-Free, Non-Linting, Ravel-Free



SKU	Description	Pkg	Price
PTC111176-2IN	2" x 4.1 yd	12/bx	\$2.59
PTC111176-3IN	3" x 4.1 yd	12/bx	\$2.89
PTC111176-4IN	4" x 4.1 yd	12/bx	\$3.59
PTC111176-6IN	6" x 4.1 yd	6/bx	\$3.59



**Dynarex® Sensi-Wrap™ Cohesive Bandage**

- Lightweight, self-adherent bandage is easily torn by hand
- Offers consistent support and compression
- Non-Sterile, Latex-Free



SKU	Description	Pkg	Price
PTC111533-1IN	1" x 5 yd	1 ea	\$1.09
PTC111533-2IN	2" x 5 yd	1 ea	\$1.99
PTC111533-3IN	3" x 5 yd	1 ea	\$2.69
PTC111533-4IN	4" x 5 yd	1 ea	\$3.49
PTC111533-6IN	6" x 5 yd	1 ea	\$4.49



**Dynarex® Elastic Bandage with Self-Closure**

- Poly-wrapped elastic bandage features self-closure
- Offers consistent stretch and maximum compression
- Non-Sterile, Latex-Free, Individually Wrapped



SKU	Description	Pkg	Price
MED131083	2" x 5"; Tan	50/bx	\$30.95
MED131086	3" x 5"; Tan	50/bx	\$51.99
MED131089	4" x 5"; Tan	50/bx	\$65.99
MED131092	6" x 5"; Tan	50/bx	\$80.99



**Dynarex® Triangular Bandage**

- 100% cotton triangular design offers maximum adaptability for wound care
- Each bandage includes two (2) safety pins
- Non-Sterile, Latex-Free, Non-Woven, Individually Wrapped



SKU	Description	Pkg	Price
EM090209	36" x 36" x 51"	12/cs	\$9.99





### Dynarex® Sterile Combine ABD Pads

- Multi-layer dressing with non-adherent outer layer
- Sterile, Latex-Free, Non-Woven

SKU	Description	Pkg	Price
PTC111156	5" x 9"	1 ea	\$0.79
PTC111121	8" x 7.5"	12/bx	\$4.49



PTC111121

## Tapes



### Dynarex® Paper Surgical Tape

- Natural material offers maximum breathability
- Non-Sterile, Latex-Free

SKU	Description	Pkg	Price
MED053498	0.5" x 10 yd	24/bx	\$9.59
MED011040	1" x 10 yd	12/bx	\$9.59
MED011042	2" x 10 yd	6/bx	\$9.59
MED053497	3" x 10 yd	4/bx	\$9.59



### Dynarex® Transparent Surgical Tape

- Ideal for securing IV lines
- Waterproof, Non-Sterile, Latex-Free

SKU	Description	Pkg	Price
PTC111592-0.5IN	0.5" x 10 yd	24/bx	\$19.49
PTC111592-1IN	1" x 10 yd	12/bx	\$16.99
PTC111592-2IN	2" x 10 yd	6/bx	\$16.99



### Dynarex® Cloth Surgical Tape

- Strong, lightweight construction provides durable skin adhesion
- Non-Sterile, Latex-Free

SKU	Description	Pkg	Price
PTC111194-0.5IN	0.5" x 10 yd	24/bx	\$18.49
PTC111194-1IN	1" x 10 yd	12/bx	\$18.49
PTC111157-2IN	2" x 10 yd	6/bx	\$18.49
PTC111194-3IN	3" x 10 yd	4/bx	\$18.49



MED900013

Hot & Cold Packs



**Dynarex® Instant Hot Packs**

- Offers instant heat therapy for minor injuries
- Latex-Free, Single-Use



SKU	Description	Pkg	Price
PTC111583	5" x 9"	1 ea	\$3.39



**Dynarex® Instant Cold Packs**

- Refrigerator-free cooling relief for minor sprains, burns and ailments
- Latex-Free, Single-Use



SKU	Description	Pkg	Price
MED090102	4" x 5"	24/cs	\$17.99
MED090103	5" x 9"	24/cs	\$62.99



**Dynarex® Reusable Hot and Cold Gel Packs**

- Pack can be either frozen or heated (microwaved)
- Retains temperature for up to twenty (20) minutes
- Reusable, Non-Toxic

SKU	Description	Pkg	Price
PTC111540-SM	4" x 6"	1 ea	\$2.29
PTC111540-MED	5" x 11"	1 ea	\$3.09
PTC111540-LG	6" x 9"	1 ea	\$3.19



Mass Casualty Incident



**North American Rescue® QuikClot® Combat Gauze LE™**

- Z-folded gauze promotes clotting within minutes of application
- Eight (8) tear points for quick application
- Heat-Free, Sterile



EM090210

SKU	Description	Pkg	Price
EM090210	3" x 4 yd	1 ea	\$42.89



**North American Rescue® QuikClot® Bleeding Control Z-fold Dressing**

- Advanced hemostatic agent works with the body's normal clotting process
- Clots up to 5-times faster than blood on its own
- Non-Allergenic, Sterile



EM090211

SKU	Description	Pkg	Price
EM090211	3" x 4 yd	1 ea	\$44.89



### North American Rescue® Combat Gauze™ Trainer

- Gauze training aid for use during instructional scenarios
- Ideal for treatment scenarios and wound packing exercises
- No hemostatic properties

SKU	Description	Pkg	Price
EM090212	3" x 4 yd; For Training Use Only	1 ea	\$10.99



EM090212

For Training Use Only



### PerSys Medical® Emergency Bandage

- Non-pneumatic trauma bandage
- Applies between 30-40+ lbs of pressure to the wound
- Non-adherent pad with built-in closure bar
- Sterile

SKU	Description	Pkg	Price
EM090213	4"	100/cs	\$815.00
EM090214	6"	100/cs	\$850.00



EM090214



### North American Rescue® HyFin® Vent Chest Seal

- Designed for the prevention, management and treatment of an open and/or tension pneumothorax potentially caused by a penetrating chest trauma
- 3-vent channels designed to prevent airflow into the chest cavity during inhalation while allowing air and blood to escape
- Latex-Free

SKU	Description	Pkg	Price
EM090215	HyFin® Vent Chest Seal	2/pk	\$19.99
EM090216	HyFin® Trainer; For Training Use Only	2/pk	\$39.99



EM090215



EM090216

For Training Use Only



### North American Rescue® ARS® for Needle Decompression

- Easy-open container for quick access
- Capless flash chamber for immediate confirmation of needle placement
- Convenient, compact size

SKU	Description	Pkg	Price
EM090217	14G x 3.25"	1 ea	\$10.95
EM090218	10G x 3.25"	1 ea	\$16.95



EM090218



### Dynarex® Post Mortem Bag Kit

- Non-porous bag with envelope-style double zipper
- Shroud sheet, toe tags, and precaution label
- Weight capacity: 400 lbs
- Dimensions: 90" x 36"

SKU	Description	Pkg	Price
PTC111421	36" x 96"; Black	1 ea	\$19.49




PTC111421




### Medline Medical Emergency Triage Tags

- Made from high-density water- and tear-resistant synthetic paper
- Each tag is coded with a unique sequential barcode serial number
- Standard graphic symbols and a four-color tear-off triaging system make the triage tag, a universal tool

SKU	Description	Pkg	Price
EM090220	Triage Tags	50/pk	\$134.22



EM090220



### Curaplex® Small Complete Triage Kit

Kit includes:

- One (1) Vest, COMMANDER, MCI Management, Royal Blue
- One (1) Vest, TRIAGE, MCI Management, Yellow
- One (1) Vest, TREATMENT, MCI Management, Red
- One (1) Vest, TRANSPORT, MCI Management, Green
- One (1) Vest, STAGING, MCI Management, Orange
- One (1) Posse BX Plain PB37C
- Twenty-Five (25) Triage Tags
- One (1) Tape, Triage Red Conterra
- One (1) Tape, Triage Yellow Conterra
- One (1) Tape, Triage Green Conterra
- One (1) Tape, Triage Black/White Conterra
- One (1) San Marcos Cert Red Bag

SKU	Description	Pkg	Price
EM090221	Triage Kit	1 ea	\$599.99



EM090221

### Safety Flag Mesh Safety Vests ANSI Class 2

- Fluorescent green, open weave mesh vest
- Includes two (2) vertical and one (1) horizontal stripes on front and back
- Vest has two (2) chest pockets and one (1) outside pocket

SKU	Description	Pkg	Price
EM090222	Large/X-Large	1 ea	\$35.79
EM090223	XX-Large/XXX-Large	1 ea	\$35.79
EM090224	XXXX-Large	1 ea	\$35.79



EM090222



# Custom Student Supply Kits

Choose from **over 500 items** to pack your kit!

STEP

01

**Pick Your Bag**

STEP

02

**Select Your Products**

Tote Bag



Oxygen Delivery



Infection Control Supplies



IV Start Supplies



Wound Care Supplies




Duffel Bag



Simulated Medications & IV Fluids




CPR Masks



Suction



ID & Organization



Backpack



Pocket Supplies



Diagnostics



Syringes & Safety Needles



Injection Supplies



\* Contact your rep for a complete list of bag options



STEP

03

**Personalize Your Kit**

STEP

04

**Contact Your  
DiaMedical Account  
Manager to Design  
Your Custom Kit**

[\(877\) 593-6011](tel:(877)593-6011)

[Info@DiaMedical.com](mailto:Info@DiaMedical.com)

IV Tubing Sets



Medication Administration Supplies



Respiratory Supplies



Print Your Logo



Embroider Your Logo







# PATIENT HANDLING & IMMOBILIZATION



stryker®

# Stryker® 6500 Power-PRO™ XT Powered Ambulance Cot

 Expertly Reconditioned by DiaMedical

- Retractable head section allows cot to be shortened in any height position for maximum versatility
- Innovative battery-powered hydraulic system that raises and lowers the patient with the touch of a button
- Easy to use manual back-up system provides complete cot operation in the event of power loss
- Includes new MedMattress Premium Bolster Cot Mattress, new restraints and two (2) batteries with one (1) charger
- Dimensions: 81" L x 23" W x 14"-41.5" H
- Weight capacity: 700 lbs
- One (1) year warranty



NEW  
Battery



NEW  
MedMattress  
Replacement Mattress



NEW  
Restraints

SKU	Description	Pkg	Price
BS047359	Reconditioned, Ambulance Cot	1 ea	\$12,350

#### PARTS & ACCESSORIES

CWS-7319B	MedMattress Premium Bolster Cot Mattress	1 ea	\$297.99
EM090411	Stryker® XPS® Retrofit Kit	1 ea	\$3,995.00
BS130109	Stryker® Defibrillator Platform	1 ea	\$1,084.99
BS130107	Stryker® Base Soft Storage Net	1 ea	\$309.95
BS130108	Stryker® Head End Flat Storage	1 ea	\$359.99

## Expertly Reconditioned by DiaMedical



#### Cleaned & Degreased

Entire Cot is Thoroughly Inspected for Any Extensive Damage to Welds and Pinch Points



#### NEW Battery

OEM Compatible Replacement Battery Providing 800+ Charge Cycles



#### Wheels

Multipoint Inspection of Wheels and Bearings



#### NEW MedMattress Replacement Mattress

Premium Mattress Upgrade is Ultra-Durable and Easy to Clean



#### NEW Restraints

OEM Certified Complete Safety Straps and Restraints





**stryker**

**Stryker® MX PRO Ambulance Cot**

- Features positive action height adjustment, seven (7) height positions and a floor-mounted safety hook
- One-hand release, fold down side rails and breakaway head section with safety bar
- Includes new MedMattress Premium Bolster Cot Mattress and new Restraints
- Dimensions: 80.5" L x 23" W x 13.5"-37.5" H
- Weight capacity: 500 lbs
- One (1) year warranty

SKU	Description	Pkg	Price
BS047407	Reconditioned, Ambulance Cot	1 ea	\$3,550.00

**PARTS & ACCESSORIES**

CWS-7319B	MedMattress Premium Bolster Cot Mattress	1 ea	\$297.99
BS130105	Stryker® MX-PRO Foot End Lead Pull Handle Kit	1 ea	\$234.00
BS130104	Stryker® MX-PRO Base Storage Net	1 ea	\$309.99
BS130106	Stryker® MX-PRO Head End Pocketed Backrest Storage Pouch	1 ea	\$399.99



BS047407

Expertly Reconditioned by DiaMedical

**stryker**

**Stryker® Performance-PRO XT Ambulance Cot**

- Features 3-position adjustable load height, fold down side rails and pneumatic gas-assist backrest
- Lightest-in-class, all-aluminum X-frame with high-visibility powder coated frame
- Includes new MedMattress Premium Bolster Cot Mattress and new Restraints
- Dimensions: 80" L x 23" W x 14"-42" H
- Weight capacity: 700 lbs
- One (1) year warranty

SKU	Description	Pkg	Price
BS047410	Reconditioned, Ambulance Cot	1 ea	\$5,050.00

**PARTS & ACCESSORIES**

CWS-7319B	MedMattress Premium Bolster Cot Mattress	1 ea	\$297.99
EM090411	Stryker® XPS® Retrofit Kit	1 ea	\$3,995.00
BS130109	Stryker® Defibrillator Platform	1 ea	\$1,084.99
BS130107	Stryker® Base Soft Storage Net	1 ea	\$309.95
BS130108	Stryker® Head End Flat Storage	1 ea	\$359.99



BS047410

Expertly Reconditioned by DiaMedical



**PROVIDE REALISM & PRECISION TO EMS TRAINING**

SimRig® allows students to train within the confined space of a real ambulance, enabling them to enter the workforce with improved efficiency and accuracy.

**SEE PAGES 84-85 FOR MORE INFORMATION**





### FERNO® 35 PROFlexx® Ambulance Cot

- Features 360° turning radius, an EZ-Pull release handle and rotate down side arms
- Longer, fixed front leg in the X-Frame and carefully engineered geometry create a low center of gravity and minimal flex
- 5-position drop frame provides better leverage for lifting and loading
- Locked on mattress ensures easy patient transfer and security during transport
- Includes bolster mattress and restraints
- Dimensions: 79" L x 24" W x 14"-38" H
- Weight capacity: 700 lbs
- One (1) year warranty



BS059907

SKU	Description	Pkg	Price
BS059907	New, Ambulance Cot	1 ea	\$9,769.00

#### PARTS & ACCESSORIES

BS047449	FERNO® Oxygen Cylinder Holder	1 ea	\$332.99
BS059909	FERNO® Dialysis Kickstand	1 ea	\$761.99
BS059910	FERNO® Cot Mounted IV Pole	1 ea	\$503.99



### FERNO® 35XST PROFlexx® StatTrac Ambulance Cot

- Features the SAE-compliant StatTrac cot fastening system, an EZ-Pull release handle and a low center of gravity
- 35XST EZ-Pull™ release system is designed to reduce wrist strain and allows first responders to get in proper positioning before assuming the weight of the patient and cot
- 5-position drop frame provides better leverage for lifting and loading
- Pull-out side lift handles provide additional grip points for assistance with heavy patients
- Dimensions: 79" L x 24" W x 34" H
- Weight capacity: 700 lbs
- Two (2) year warranty



BS052832-R-UA

SKU	Description	Pkg	Price
BS052832-R-UA	New, Ambulance Cot	1 ea	\$10,445.00



### FERNO® POWER X1 Powered Ambulance Cot

- Industry-leading battery technology which provides high-powered and reliable performance on every call
- Programmable load height adjusts to match the height of your ambulance floor
- Features four (4) 6" x 2" oversized swivel wheels with front-end locks for added maneuverability
- Back-rest mounted equipment hook holds equipment to keep hands free
- Pneumatic backrest adjusts to a comfortable position for each patient
- Shock frame allows you to raise patient's feet for treatment of OTC shock
- Includes bolster mattress and FERNO® X-restraints
- Dimensions: 80" L x 24" W x 44" H
- Weight capacity: 700 lbs
- Two (2) year warranty



BS047444

SKU	Description	Pkg	Price
BS047444	New, Ambulance Cot	1 ea	\$19,989.00

#### PARTS & ACCESSORIES

BS047445	FERNO® O2 Bottle Holder Kit	1 ea	\$239.99
BS047446	FERNO® X-Restraint Kit	1 ea	\$293.99
BS047447	FERNO® Power X1 IV Pole	1 ea	\$529.99

**stryker**

**Stryker® 6252 Stair-Pro Stair Chair**

- Innovative Stair-TREAD system allows caregivers to transport patients down stairs without lifting
- Features locking rear life handles and extendable foot and lift handles to support ergonomic lifting technique and enhanced up stair mobility
- Upper control handle can be extended for optimal operatory control
- Compact folded design with easy-to-fold functionality allows for convenient storage
- Weight capacity: 500 lbs
- One (1) year warranty

SKU	Description	Pkg	Price
BS047398	Reconditioned, Stair Chair	1 ea	\$3,400

**PARTS & ACCESSORIES**

BS130113	Stair Chair Extended Foot Rest	1 ea	\$424.99
BS130114	Stair Chair Vinyl Head Support	1 ea	\$89.99



BS047398

Expertly Reconditioned by DiaMedical

**FERNO**

**FERNO® Transcend™ Stair Chair with PowerTraxx®**

- PowerTraxx® feature allows the chair and track system to carry the weight which increases patient and operator safety
- Includes ENDURACharge 36V battery which runs up to forty (40) flights of stairs per charge
- 5-position lift bar and handles offers operators of various heights maximum control and helps reduce the strain of lifting
- Built-in electronic speed control limits track speed during decent for enhance patient and operator safety
- 1" track-to-ground clearance provides smooth rolling over carpet and rough surfaces
- Positive locking hinge mechanism ensures chair is locked in place whether deployed for use or folded for carrying and storage
- Weight capacity: 550 lbs
- One (1) year warranty

SKU	Description	Pkg	Price
EM090165	New, Battery Powered Stair Chair	1 ea	\$9,549.00



EM090165

**FERNO**

**FERNO® Transcend Stair Chair**

- 5-position lift bar and handles offers operators of various heights maximum control and helps reduce the strain of lifting
- All metal construction increases durability and extends the chairs lifespan
- Extending footrest and lower track angle increases patient comfortability
- 1" track-to-ground clearance provides smooth rolling over carpet and rough surfaces
- Positive locking hinge mechanism ensures chair is locked in place whether deployed for use or folded for carrying and storage
- Weight capacity: 550 lbs
- One (1) year warranty

SKU	Description	Pkg	Price
EM090166	New, Manual Stair Chair	1 ea	\$4,584.00



EM090166



**Spencer 4BELL STAIR+**

- Foldable transport and evacuation chair designed for smooth patient transport and operator safety
- Features a reinforced, lightweight frame with telescopic front handles, integrated belts and slides
- Pivoting front wheels increases chair maneuverability and stability
- Soft backrest and semi-rigid seat enhances patient comfort
- Weight capacity: 396 lbs (180kg)

SKU	Description	Pkg	Price
EM090438	New, Stair Chair	1 ea	\$4,395.00



EM090438

**MOBI 5C PRO Stair Stretcher**

- Features non-slip front and rear handles and ABS seating
- Compact, foldable chair design make storage easy
- Made from high-strength aluminum alloy
- Weight capacity: 450 lbs
- One (1) year warranty

SKU	Description	Pkg	Price
BS047435	New, Manual Stair Chair	1 ea	\$319.00



BS047435



**GLOVES & INFECTION CONTROL**

**WOUND CARE SUPPLIES**

**OXYGEN & RESPIRATORY ACCESSORIES**

**QUICK SHIP  
DISPOSABLE MEDICAL PRODUCTS**

In efforts to keep your facility stocked and prepared at all times, DiaMedical warehouses over 3,000 medical products available to ship immediately.

(877) 593-6011 | [DiaMedical.com](http://DiaMedical.com)





Cot Mattresses



# MedMattress Premium Bolster Cot Mattress

- Designed as a replacement bolster cot mattress for non-XPS Stryker Emergency Cots
- Durable heavy duty black vinyl resistant to cracking and easy to clean/sanitize
- Sealed seams for increased infection control such as mitigating fluid penetration
- Features side bolsters to minimize lateral movement and enhance patient comfort
- Dimensions: 73" L x 18" W x 4" H with Bolsters



+ customer favorite

SKU	Description	Pkg	Price
CWS-7319B	Replacement Mattress for non-XPS Stryker Cots	1 ea	\$297.99



RF Welded Sealed Seams



Heavy Duty Vinyl Cover



Side Bolsters

## MedMattress Premium Bolster Cot Mattress Features



**RF Welded Sealed Seams**

Increase Infection Control by Mitigating Fluid Penetration



**Heavy Duty Vinyl Cover**

Easy to Clean, Durable Heavy Duty Black Vinyl is Crack-Resistant



**Side Bolsters**

Designed to Cradle the Patient and Minimize Lateral Movement



**Dimensions**

18" W x 73" L x 4" H with Bolsters



\*Cot mattress designed & crafted by MedMattress, universal replacement mattress compatible with Stryker® EMS Stretchers



### Medline Critical Response™ Mattress

- Heavy-duty vinyl is highly resistant to tearing
- Mattress designed to increase patient comfort while reducing transfer gap
- Made for Stryker XPS stretchers only, must be used in conjunction with Stryker XPS side rails

SKU	Description	Pkg	Price
EM090235	Replacement Mattress for XPS Stryker® Cots	1 ea	\$297.99



EM090235



### FERNO® PROFlexx Series Replacement Mattress

- Black vinyl-coated nylon cover features heat-sealed seams
- Abrasion-Resistant, Fire-Resistant, Bacteriostatic

SKU	Description	Pkg	Price
EM090236	Replacement Mattress for FERNO® PROFlexx Series Cots	1 ea	\$498.50



EM090236



### FERNO® Power X1 Replacement Mattress

- Designed specifically for the POWER X1 Ambulance Cot
- Internally bolstered, externally sealed
- Dimensions: 73" L x 18" W x 5" H"

SKU	Description	Pkg	Price
EM090237	Replacement Mattress for FERNO® Power X1 Cot	1 ea	\$420.00



EM090237

## Cot Accessories



### Stryker® Power-LOAD® Powered Cot Fastener

- Lifts or lowers the cot into and out of the ambulance, eliminating spinal loads that can result in cumulative trauma injuries
- Minimizes patient drops by supporting the cot until the wheels are on the ground
- Features a manual back-up system for operation in the event of power loss
- Includes floor plate, head end lock LED indicators and installation hardware
- Compatible with Stryker Power-PRO XT, Stryker Power-PRO IT and Stryker Performance-PRO XT
- Meets SAE J3027 crash test safety standards
- Ninety (90) day warranty

SKU	Description	Pkg	Price
BS047413	Reconditioned, Powered Cot Fastener	1 ea	\$23,381.00

#### Features

- Lifting arms
- Head end lock LED indicators
- Control panel
- Cot release handles
- Foot end lock LED indicator
- Linear transfer system
- Inductive charging
- Operation guide
- Includes floor plate and installation hardware

#### Compatibility

The Power-LOAD is compatible with the Power-PRO XT and Power-PRO IT. This system meets dynamic crash test standards for maximized occupant safety and will automatically charge the Power-PRO XT and Power-PRO IT SMRT battery.



BS047413



### Stryker® Performance-LOAD® Manual Cot Fastener

- Includes head end lock and foot-end latch to secure both ends of the cot while in transport position
- Linear transfer system guides and centers the cot during loading and unloading
- Compatible with Stryker Power-PRO 2, Stryker Power-PRO XT, Stryker Performance-PRO XT and Stryker Power-PRO IT
- Meets SAE J3027 crash test safety standards

SKU	Description	Pkg	Price
EM090238	Reconditioned, Manual Cot Fastener	1 ea	\$3,795.00



EM090238



### FERNO® iNLINE® Cot Fastening System

- Unique in-line design prevents cot drift and eliminates the need to steer when loading and unloading
- One touch tab provides easy release from fastener
- Designed to securely retain the iNX® Integrated Patient Transport & Loading System
- Built-in safety hook ensures that cot remains secured while operators raise or lower the undercarriage during loading or unloading
- Fits FERNO® Stat Trac® and Stryker Power-LOAD® floor plates for interchangeability between vehicles
- Meets SAE J3027 crash test safety standards

SKU	Description	Pkg	Price
EM090363	New, 86" Cot Fastener	1 ea	\$5,438.00
EM090364	New, 95" Cot Fastener	1 ea	\$5,438.00
EM090365	New, 102" Cot Fastener	1 ea	\$5,438.00



EM090364



### FERNO® PRO F1™ Universal Cot Fastener

- Cross-compatible design accommodates the FERNO® POWER X1, Stryker Power-PRO™ XT and Stryker Performance-PRO™XT ambulance cots
- One-touch foot-end release button makes unfastening cots simple and quick
- 110° approach angle when combined with the POWER X1 Ambulance Cot
- Optional Integrated Charging System (ICS®) provides in-fastener charging for the FERNO® POWER X1 cot
- Meets SAE J3027 crash test safety standards
- Two (2) year warranty

SKU	Description	Pkg	Price
EM090242	New, Manual Universal Cot Fastener	1 ea	\$4,896.00
EM090366	New, Universal Cot Fastener with Integrated Charging System	1 ea	\$7,854.00



EM090366



### FERNO® Stat Trac Cot Fastening System

- Compatible with the FERNO® POWERFlexx ST and FERNO® 35XST PROFlexx Ambulance Cot
- Head-end and foot-end, dual-locking safety mechanism ensures a secure and comfortable ride for the patient
- Single-lever release and a safety stop provide security while loading and unloading
- Straight in-track design eliminates the need to steer the cot when loading and unloading
- A simple, single-hand release unlocks the cot from the Stat Trac

SKU	Description	Pkg	Price
EM090238	New, 86" Cot Fastener	1 ea	\$2,577.00
EM090239	New, 91" Cot Fastener	1 ea	\$3,047.00
EM090240	New, 96" Cot Fastener	1 ea	\$3,047.00
EM090241	New, 103" Cot Fastener	1 ea	\$3,047.00




EM090239



### FERNO® Model 175 Universal Floor Mount

- Designed to limit the movement of a FERNO® ambulance cot inside the compartment of a ground-based ambulance
- Compatible with FERNO® POWER X1 and FERNO® 35X ambulance cot
- Heavy duty stainless steel bracketry and rail fastener construction

SKU	Description	Pkg	Price
SRM00100	Antler and Rail Fastening System, Floor Mount	1 ea	\$1,413.00



EM090366

### Replacement Battery Kit for Stryker® Power-PRO™ Cots

- Not compatible with Stryker® SMRT Power Systems
- Voltage: 24v
- Capacity: 3300mAh / 3.3AH
- Chemistry: Ni-CD
- Bundle includes two (2) batteries and one (1) charger

SKU	Description	Pkg	Price
DW0246B	Replacement Battery & Charger Bundle; Non-OEM	1 kit	\$354.99
DW0242	Replacement Battery; Non-OEM	1 ea	\$129.95
DW0246	Replacement Charger; Non-OEM	1 ea	\$119.95



DW0246B



### EMS Solutions Emergency Child Restraint System

- Features three (3) color coded, full adjustable harnesses in sizes small, medium and large
- Load Capacity: Small: 4-18 lbs; Medium: 13-46 lbs; Large: 42-110 lbs
- No stitching or hook and loop fasteners ensure maximum infection control
- Exceeds NASEMSO standards

SKU	Description	Pkg	Price
EM090253	Child Restraint System	1 ea	\$858.99



EM090253



### XDcuff® Starter Pack

- Compatible with Stryker® and FERNO® emergency cots
- Pack includes four (4) permanent clip anchors, one (1) pair of reusable ankle restraints and one (1) pair of reusable wrist restraints
- XDcuff® pre-connects the anchors (clips) to the cot allowing medics to quickly anchor limb restraints

SKU	Description	Pkg	Price
EM090244	Clip-In; Plastic	1 kit	\$134.95
EM090245	Clip-In; Metal	1 kit	\$154.95
EM090246	Secure-In; Plastic	1 kit	\$154.95
EM090247	Secure-In; Metal	1 kit	\$212.95



EM090244

+ customer favorite





### FERNO® NeoMate™ Infant Restraint System

- Fully-adjustable, 5-point harness system transforms ambulance cot to safely transport infants
- Halo Pad keeps the infant's head stable during transport
- Rolls compactly for convenient storage and easy deployment
- Vinyl construction is nontoxic and easy-to-clean
- Three (3) restraint straps easily attach to any cot
- Load limit: 5-14 lbs



SKU	Description	Pkg	Price
EM090255	Infant Restraint System	1 ea	\$482.99



### FERNO® Pedi Mate® + Pediatric Restraint System

- Fully-adjustable, 5-point harness system securely holds patients, providing safe restraint for transport
- Rolls compactly for convenient storage and easy deployment
- Vinyl construction is nontoxic and easy-to-clean
- Load limit: 10-100 lbs



SKU	Description	Pkg	Price
EM090254	Pediatric Restraint System	1 ea	\$482.99

### SnugFit® Fitted Barrier Sheet

- Non-woven, poly barrier
- Designed to prevent slippage during patient movement
- Help prevent equipment contamination



SKU	Description	Pkg	Price
EM019412	30" x 84"; Blue	50/cs	\$199.99



### Curaplex® Fitted Stretcher Sheet for Stryker XPS Mattress

- Designed to fit the Stryker® XPS Mattress
- Made from spunbound polypropylene
- Disposable, fluid-resistant
- Lift weight capacity: 400 lbs



SKU	Description	Pkg	Price
EM090243	36" x 90"; Blue	30/cs	\$102.95



### DMS Stair Chair Restraint

- 2-point adjustable 5' strap features looped ends
- Buckle: Plastic and metal push button
- Material: Nylon
- Color: Black



SKU	Description	Pkg	Price
EM090251	2-Point Nylon Strap with Metal Buckle & Looped End; 5'	1 ea	\$24.95



### DMS Nylon Stretcher/Cot Straps

- 2" x 5'
- Buckle: Metal push button
- Material: Nylon
- Color: Black



SKU	Description	Pkg	Price
EM090248	Nylon Strap with Metal Push Button Buckle & Loop Ends	1 ea	\$19.99



### DMS Stair Chair Ankle Restraint

- 2-point adjustable 4.5' strap features looped ends
- Buckle: Plastic and metal push button
- Material: Polypropylene
- Color: Black



SKU	Description	Pkg	Price
EM090252	1-Piece Polypropylene Strap with Plastic Buckle; 4.5'	1 ea	\$14.95



### DMS 4-Point Shoulder Harness Restraint System

- Made up of two (2) 60" adjustable shoulder straps and one (1) 5' double adjusting torso strap
- Buckle attaches into L-links to secure on patient
- Buckle: Metal push button
- Material: Nylon
- Color: Black




SKU	Description	Pkg	Price
EM090250	4-Point Nylon Harness with Metal Buckle & Looped End	1 ea	\$109.99

Backboards & Transfer Aids



Iron Duck® Baseboard

- Rated for the lowest amount of deflection in the 400-500 lbs weight class
- Made of a roto-molded HDPE shell over non-toxic foam for durability and cleaning
- Fourteen (14) over-sized hand holds can accommodate gloved hands
- X-Ray and CT scan translucent, MRI scan compatible applications
- Dimensions: 72" L x 16" W x 1.75" H
- Load capacity: 500 lbs
- Float capacity: 250 lbs
- Additional colors available online



SKU	Description	Pkg	Price
EM019301-YEL	Yellow	1 ea	\$156.45



Iron Duck® Ultra Short Board

- Rugged alternative to standard length spineboards
- Ideal for close-space extraction such as motor vehicle accidents
- Six (6) over-sized hand holds accommodate gloved hands
- High-density polyethylene shell and filled with eco-friendly foam for durability and easy cleaning
- Available with or without eight (8) board pins
- Dimensions: 31.25" L x 15.75" W x 1" H
- Weight: 5 lbs
- Additional colors available online



SKU	Description	Pkg	Price
EM019302-ORG	Standard; Orange	1 ea	\$65.95
EM019310-ORG	Pinned; Orange	1 ea	\$78.95



FERNO® Millennia Plastic Backboard

- Easy-to-clean plastic shell with foam-filled construction is lightweight yet durable
- Special mounting slots for the FERNO® Model 445 Head Immobilizer
- X-ray-translucent
- Minimal deflection for improved stabilization
- Dimensions: 72" L x 16" or 18" W
- Load capacity: 350 lbs
- Weight: 12 lbs



SKU	Description	Pkg	Price
EM012801-ORG	16"; Orange	1 ea	\$450.99
EM012801-YEL	16"; Yellow	1 ea	\$450.99
EM012802-ORG	18"; Orange	1 ea	\$450.99



Junkin Plastic Backboard

- Rugged single piece lightweight polyethylene spineboard with twelve (12) large hand holes for easy handling
- Built in runners
- Low profile and X-ray translucent
- Arrives with four (4) straps
- Dimensions: 72" L x 16" W
- Load capacity: 400 lbs
- Weight: 15 lbs



SKU	Description	Pkg	Price
EM019401	Bright Yellow	1 ea	\$272.99

Reeves® Flexible Stretcher

- Flexible design allows for maneuvering around corners or on staircases
- Features five (5) removable hardwood slats, six (6) ergonomically placed handles and three (3) patient straps
- Dimensions: 78" L x 28" W
- Load capacity: 700 lbs
- Weight: 15 lbs



SKU	Description	Pkg	Price
EM090256-ORG	Orange	1 ea	\$529.99
EM090256-GRY	Gray	1 ea	\$529.99
EM090256-BLK	Black	1 ea	\$529.99

Graham Medical MegaMover® Portable Transport Unit

- Durable, cost effective patient transport device
- Strong, disposable and compact alternative for transporting patients through confined spaces or restricted areas
- Fourteen (14) sewn-in carrying handles
- Dimensions: 80" L x 40" W
- Load capacity: 1,000 lbs



SKU	Description	Pkg	Price
EM019411	White/Black	10/cs	\$279.50



### FERNO® Basket Stretcher

- High-density, polyethylene shell is supported by a permanently attached, heavy-duty, aluminum frame
- Shell is chemical, UV, rust and corrosion-resistant, ensuring long life
- Molded-in runners and handholds make the stretcher easy to handle
- Dimensions: 80" L x 40" W
- Load capacity: 600 lbs
- Weight: 23 lbs

SKU	Description	Pkg	Price
EM090257	Orange/Black	1 ea	\$1,391.00



EM090257



### Spencer Shell Basket Stretcher

- Orange high-density polyethylene stretcher is joined to an aluminum frame and has a uniform thickness for strength and support
- Includes a removable mattress, applied with Velcro
- Features twelve (12) built-in handles and three (3) nylon patient straps
- Dimensions: 84.25" L x 25.6" W
- Load capacity: 617 lbs
- Weight: 27.5 lbs

SKU	Description	Pkg	Price
EM090147	Basket Stretcher; Orange	1 ea	\$799.00



EM090147



### Junkin Basket Stretcher

- High-density polyethylene shell, supported by a stainless steel outer rail, permanently attached with stainless steel semi-tubular rivets
- Features molded runners, fully exposed outer rail, non-absorbent foam pad secured to stretcher and four (4) patient restraint straps
- Break-away stretcher separates into two (2) halves and utilizes positive-locking quick-release pins to secure the halves together
- Dimensions: 84.5" L X 24" W X 7.5" H
- Load capacity: 1,200 lbs
- Weight: 31/32 lbs

SKU	Description	Pkg	Price
EM090258	Plastic Stretcher; Yellow	1 ea	\$1,015.99
EM090259	Break-Away Plastic Stretcher; Yellow	1 ea	\$1,326.00



EM090259



### FERNO® Scoop EXL Stretcher

- Interlocking halves allow operators to "scoop" patients up, eliminating log-roll maneuvers
- Head section is recessed for proper cervical alignment
- Radiolucent
- Dimensions: 79" L x 17" W x 3" D
- Load capacity: 500 lbs
- Weight: 18 lbs

SKU	Description	Pkg	Price
EM090260	Scoop™ EXL with 430P Restraint Straps	1 ea	\$1,422.00
EM090261	Scoop™ EXL with Pins & 436 Speed Clip Restraints with Metal Buckle	1 ea	\$1,544.00



EM090260

Backboard Accessories



**Kemp USA 10-Point Patient Restraint Spineboard Straps**

- Color-coded 10-point straps are made of heavy duty nylon and fasten with a hook and loop
- Holds patient securely and comfortably to spineboard



SKU	Description	Pkg	Price
EM090263	Color-Coded Spineboard Strap	1 ea	\$49.95



**Kemp USA Spineboard Straps**

- Made of high quality nylon with extremely durable plastic clips
- Designed for fast, firm immobilization on all types of spineboards



SKU	Description	Pkg	Price
EM090264	One (1) Piece with Plastic Buckle; Black	1 ea	\$9.00
EM090265	Two (2) Piece with Plastic Buckle; Black	1 ea	\$8.50



**Kemp USA Spineboard Mounting Bracket**

- Made of high quality nylon with extremely durable plastic clips
- Designed for fast, firm immobilization on all types of spineboards



SKU	Description	Pkg	Price
EM090266	Universal Wall Mount for Backboards	1 ea	\$55.50

EM090266

Blankets



**Dynarex® Survival Rescue Blanket**

- Mylar foil material assists in body heat retention in emergency situations
- Non-Sterile, Latex-Free, Individually Wrapped, Disposable

SKU	Description	Pkg	Price
MED160104	84" L x 52" W	1 ea	\$1.69



MED160104



**Dynarex® Emergency Highway Blankets**

- High-quality, water-resistant blankets designed to provide warmth and coverage to help prevent shock in emergency situations
- Latex-Free, Disposable

SKU	Description	Pkg	Price
MED160110	80" L x 54" W; Economy	50/cs	\$160.99



MED160110



Head & Neck Immobilizations



**FERNO® Model 125 K.E.D.**

- Immobilizes the torso, head, and neck enabling prompt extrication of patients from confined spaces
- Two (2) part wrap-around design provides horizontal flexibility and vertical rigidity
- Adapts to a hip and pelvic splint by simple inversion
- May be used for children and pregnant women
- Color-coded sewn-in securing straps and snaplock buckles



SKU	Description	Pkg	Price
EM012808	Kendrick Extrication Device	1 ea	\$199.99



**Laerdal Stifneck® Extrication Collar**

- Easy access for pulse checks, advanced airway procedures and visualization through oversized trachea hole
- Stabilize spinal patients correctly with four (4) sizing options (adult) or three (3) sizing options (pediatric) that all lock quickly and safely into place ensuring symmetrical adjustments
- Radiolucent, MRI and CT-scan compatible
- Latex-Free



SKU	Description	Pkg	Price
EM090267	Stifneck® Select™; Adult	1 ea	\$15.99
EM090268	Stifneck® Pedi-Select Collar; Pediatric	1 ea	\$15.99



**Laerdal Stifneck® No-Neck Extrication Collar**

- Easy access for pulse checks, advanced airway procedures and visualization through oversized trachea hole
- Radiolucent, MRI and CT-scan compatible
- Latex-Free



SKU	Description	Pkg	Price
EM090269	Stifneck® No-Neck Collar; Baby	1 ea	\$11.99
EM090270	Stifneck® No-Neck Collar; Adult	1 ea	\$11.99
EM090271	Stifneck®; Pediatric	1 ea	\$11.99
EM090272	Stifneck® Regular; Adult	1 ea	\$11.99

**Philadelphia Patriot Collar**

- Quick multi-height adjustability for a one-size-fits-all design
- Features a large trachea opening for quick access to carotid pulse monitoring and emergency tracheotomies



SKU	Description	Pkg	Price
EM090434	Adult	1 ea	\$15.95
EM090433	Pediatric	1 ea	\$15.95



**Junkin Head Immobilizer**

- Includes one (1) base plate, two (2) head blocks and two (2) head/chin straps
- Vinyl coated, buoyant and waterproof
- Can be left in place during X-Ray



SKU	Description	Pkg	Price
EM019402	Head Immobilizer; Yellow	1 ea	\$139.99



**Iron Duck® Sof Loc Head Immobilizer**

- Includes one (1) base, two (2) head blocks, one (1) head strap and one (1) chin strap
- Vinyl coated head blocks provide maximum support and minimum head movement




SKU	Description	Pkg	Price
EM019307	Head Immobilizer; Blue	1 ea	\$163.86



### Kemp USA Head Immobilizer Set

- Includes one (1) base plate, two (2) head block pillows and two (2) head/chin straps
- Features two (2) large circular ear holes to monitor bleeding and fluid drainage
- Can be worn during an X-Ray, MRI and CT scan



EM090432-ORG

SKU	Description	Pkg	Price
EM090432-BLK	Head Immobilizer; Black	1 ea	\$94.50
EM090432-RED	Head Immobilizer; Red	1 ea	\$94.50
EM090432-ORG	Head Immobilizer; Orange	1 ea	\$94.50
EM090432-RBL	Head Immobilizer; Royal Blue	1 ea	\$94.50

## Splints




### SAM® Splint

- Created from a thin core of aluminum alloy, sandwiched between two layers of closed-cell foam
- Pliable design becomes extremely strong when bent into shape
- Lightweight, radiolucent, waterproof



### SAM® Pelvic Sling II

- 1-piece design safely and effectively reduces and stabilizes open-book pelvic ring fractures
- SAM® AUTOSTOP buckle acts to avoid over or under tightening
- Radiolucent and MRI safe



SKU	Description	Pkg	Price
EM090281	36" L x 4.25" W, Rolled; Orange/Blue	1 ea	\$25.99
EM090282	18" L x 4.25" W, Flat Folded; Orange/Blue	1 ea	\$14.99

SKU	Description	Pkg	Price
EM090273	Small, 27" - 47"	1 ea	\$98.99
EM090274	Medium, 32" - 50"	1 ea	\$98.99
EM090275	Large, 36" - 54"	1 ea	\$98.99




### FernoTrac™ Traction Splint Combination Kit

- Includes one (1) adult traction splint with case and one (1) pediatric traction splint with case
- Low-profile ischial pad permits splint to be applied with minimum movement of the patient's leg
- Simple adjustable wrap-around ankle hitch designed to allow EMTs to check pedal pulses



SKU	Description	Pkg	Price
EM012809	Traction Splint Kit	1 ea	\$1,051.00



### Curaplex® Traction Splint

- Realigns a fractured femur to the proper position
- Features telescoping design
- Includes carrying case



SKU	Description	Pkg	Price
EM090276	Adult	1 ea	\$270.99
EM090277	Pediatric	1 ea	\$324.99



### Sager® S-301 Form III Single Lower Limb Splint

- Apply dynamic traction and document exactly how much traction you have applied to your patient
- Designed for use with adults or children



SKU	Description	Pkg	Price
EM090278	Splint with Soft Case	1 ea	\$779.95



### DMS Padded Board Splint Set

- Sealed in reusable orange vinyl bags
- Features 0.5" of foam on one side for added patient comfort
- Set includes two (2) 15", two (2) 36 and two (2) 54" splints



SKU	Description	Pkg	Price
EM090279	Splint Set with Carrying Case	1 kit	\$64.95





**MONITORING**  
& DIAGNOSTICS

**AEDs**

**PHILIPS**

**Philips HeartStart FRx AED**

- Arrives ready to use with SMART Pads II and battery already installed
- FRx AED can withstand rough handling, extreme temperatures and dusty or wet surfaces
- Tested to withstand loads up to 1,100 lbs and drops from 4'
- FRx provides real-time voice guidance and runs automatic self-diagnostic check
- One (1) year warranty



DF026203

SKU	Description	Pkg	Price
DF026203	New, Automated External Defibrillator	1 ea	\$2,254.00

**PARTS & ACCESSORIES**

DF026203EL	HeartStart FRx AED SMART Pads II	1 ea	\$62.00
DF026203PD	HeartStart FRx AED Infant/Child Key	1 ea	\$111.95
DF026203YC	Philips HeartStart FRx Waterproof Hard Case	1 ea	\$269.00

**PHILIPS**

**Philips HeartStart OnSite AED**

- Arrives ready to use with Adult SMART Pads Cartridge and battery already installed
- Voice directions guide you through every step and instructions are repeated as needed
- SMART pads on the patient's chest to monitor the heart rhythm and chest compressions
- Eight (8) year warranty



DF026201

SKU	Description	Pkg	Price
DF026201	New, Automated External Defibrillator	1 ea	\$1,584.00

**PARTS & ACCESSORIES**

DF152040	HeartStart OnSite Adult SMART Pads	1 ea	\$109.95
DF152041	HeartStart OnSite Infant/Child SMART Pads	1 ea	\$115.95
DF152042	Philips HeartStart OnSite/FRX Replacement Battery	1 ea	\$189.00



**Cardiac Science® Powerheart® G5 AED**

- Arrives with adult pads and battery
- Allows rescuers to adjust prompts to match their skill level as well as change shock and CPR protocols to match clinical preferences
- Settings can easily be updated as local protocols and resuscitation guidelines change
- Shock level is customized to deliver the appropriate energy dose to the patient
- Defibrillation energy automatically escalates within 10 seconds or less
- Dual-language
- Eight (8) year warranty



DF028345-AN

SKU	Description	Pkg	Price
DF028345-SA	New, Semi-Automatic	1 ea	\$2,618.56
DF028345-AN	New, Fully Automatic	1 ea	\$2,618.56

**PARTS & ACCESSORIES**

DF152006	Adult Intellisense Defibrillation Pads	1 ea	\$81.00
DF152007	Pediatric Intellisense Defibrillation Pads	1 ea	\$110.00



## ZOLL

### ZOLL® AED Plus Defibrillator

- Features step-by-step text and voice prompts to help rescuers
- Real CPR Help® technology includes a compression depth indicator that provides real-time CPR feedback
- ZOLL® CPR-D-padz have a unique one-piece design to ensure easy and accurate placement
- Can withstand extreme conditions including extreme temperature, altitude, vibration and shock
- Five (5) year warranty

SKU	Description	Pkg	Price
DF028302	New, Semi-Automatic	1 ea	\$2,921.00
DF152079	New, Fully Automatic	1 ea	\$2,921.00
PARTS & ACCESSORIES			
DF152078	ZOLL® Replacement CPR-D-Padz	1 ea	\$223.00
DF152016	ZOLL® Pedi-Padz II	1 ea	\$150.00
DF152077	ZOLL® AED Plus Replacement Lithium Batteries	10/pk	\$97.00



DF152079

## ZOLL

### ZOLL® AED 3® Defibrillator

- Features a full-color display with vivid rescue images, a CPR cycle timer and a large color bar gauge that shows CPR compression depth
- Real CPR Help® technology includes a compression depth indicator that provides real-time CPR feedback
- RapidShock™ analysis allows ZOLL® AED 3® to deliver a necessary shock in five (5) seconds
- CPR Uni-padz® electrodes provide a solution for both adult and pediatric victims
- Intuitive touchscreen allows quick access to AED information and enables configuration
- Five (5) year warranty

SKU	Description	Pkg	Price
DF152052	New, Semi-Automatic	1 ea	\$2,144.00
DF152053	New, Fully Automatic	1 ea	\$4,101.00
DF152054	BLS for Professionals and First Responders	1 ea	\$3,195.00
PARTS & ACCESSORIES			
EM090283	ZOLL® CPR Uni-padz® III	1 ea	\$234.00
EM090284	ZOLL® AED 3 Non-Rechargeable Battery	1 ea	\$195.00



DF152052

## ♥ defibtech

### Defibtech Lifeline

- Lightweight, compact design is simple for untrained responders to use
- Lifeline AED has only two (2) buttons - green to turn it on and red to deliver a shock if needed
- Calm, clear voice prompts lead rescuers through all steps of defibrillation
- Active Status light shows that the unit is ready to perform any time, every time

SKU	Description	Pkg	Price
DF029701	New, Semi-Automatic	1 ea	\$1,495.50
EM090285	New, Fully Automatic	1 ea	\$1,595.50
PARTS & ACCESSORIES			
EM090286	Defibtech Lifeline Adult Defibrillation Electrode Pads	1 ea	\$68.00
EM090287	Defibtech Lifeline Pediatric Defibrillation Electrode Pads	1 ea	\$117.89
DCF-200	Defibtech Lifeline AED Standard Battery Pack	1 ea	\$189.95



DF029701

### Defibtech Lifeline View

- The first and only AED with an interactive display and integrated full color video to demonstrate chest compressions and defibrillation
- Rescuers can easily switch from compression-only CPR to compressions with rescue breathing
- Suitable for use in high vibration settings including jet aircraft, helicopters and boats
- One-touch status screen shows readiness of both electrode pads and battery
- Lifeline ECG allows the user to display patient's ECG data and event information while in AED mode using defibrillation pads

SKU	Description	Pkg	Price
DF152055	New, Semi-Automatic	1 ea	\$1,880.00
EM090289	New, Semi-Automatic with ECG Display	1 ea	\$2,595.00



DF152055

**PARTS & ACCESSORIES**

EM090290	Adult Electrodes for Defibtech Lifeline VIEW AED	1 ea	\$75.17
EM090291	Pediatric Electrodes for Defibtech Lifeline VIEW AED	1 ea	\$123.50

### AED Trainers

#### ZOLL

##### ZOLL® AED Plus Trainer2 AED Trainer

- Simulates how the ZOLL® AED Plus Defibrillator will function during a SCA event
- Prepares students to save a life with four (4) different fixed rescue scenarios or manual control with a wireless remote controller
- Trainer provides audible real-time feedback about compression depth and rate
- One (1) year warranty

SKU	Description	Pkg	Price
DF018301	AED Trainer; For Training Purposes Only	1 ea	\$446.00



DF018301

For Training Use Only

#### PHILIPS

##### Philips HeartStart OnSite AED Trainer

- Simulates how the Philips HeartStart OnSite AED will function during a SCA event
- Voice directions guide you through every step and instructions are repeated as needed
- HeartStart Trainer includes nylon carrying case, one set of Adult Training Pads and one External Manikin Adapter Strip for use with standard manikins

SKU	Description	Pkg	Price
EM090435	AED Trainer; For Training Purposes Only	1 ea	\$469.00



EM090435

For Training Use Only

**PARTS & ACCESSORIES**

EM090436	HeartStart OnSite AED Replacement Adult SMART Training	1 ea	\$45.00
----------	--	------	---------



### Defibtech Standalone Trainer AED

- Simulates how the Defibtech Lifeline AED will function during a SCA event
- Voice directions guide you through every step and instructions are repeated as needed
- Includes six (6) different rescue scenarios or manual control with a wireless remote controller
- One (1) year warranty

SKU	Description	Pkg	Price
DF019701	AED Trainer; For Training Purposes Only	1 ea	\$499.99

#### PARTS & ACCESSORIES

DF152037	Adult Training Pad Set	1 ea	\$31.83
DF152038	Pediatric Training Pad Set	1 ea	\$49.95
DF152039	Battery Charger	1 ea	\$84.95



DF019701

For Training Use Only



### Cardiac Science Full-Size AED Trainer

- Dual mode AED trainers with infrared remote control for teaching both Semi-Automatic or Fully-Automatic defibrillation
- Students can practice using both adult and pediatric AED trainer
- Can be paused at any time and picks up the simulation where it left off
- Simulates how the Cardiac Science Powerheart G3 AED will function during a SCA event

SKU	Description	Pkg	Price
EM090292	AED Trainer; For Training Purposes Only	1 ea	\$454.00

#### PARTS & ACCESSORIES

EM090293	TRAINING Electrode Pads, Adult	1 ea	\$28.00
EM090294	TRAINING Electrode Pads, Pediatric	1 ea	\$83.00



EM090292

For Training Use Only

## Defibrillators

### ZOLL

#### ZOLL® X Series Defibrillator

- Designed to meet the specifications of military and EMS
- Equipped with Real CPR Help® giving the user real-time assistance and feedback in the event of a sudden cardiac emergency
- All-inclusive unit with unmatched mobility
- Weighs less than 12 lbs making it lighter than other full-featured devices
- Large, bright screen allows for simultaneous viewing of four traces, all 12 ECG leads, or static and dynamic 12-lead views on a "split screen"
- CPR Dashboard provides real-time numeric depth and rate values and has a CPR timer that switches to an idle timer when compressions are stopped
- Features integrated WiFi, Bluetooth® and USB cellular modem capabilities
- Medics can now directly send X Series files to your system with the push of a button from RescueNet CaseReview
- One (1) year warranty

SKU	Description	Pkg	Price
DF152090	Reconditioned; 12-Lead with AED, Biphasic Pacing, SpO2, SPCO, NIIBP & EtCo2	1 ea	\$27,850.00



DF152090

## Physio-Control LIFEPAK® 15 Defibrillator

- Continuous 12-lead ECG monitoring
- ST-Segment Trending feature
- Compact, portable and easy-to-use construction
- Features six (6) main operating modes: AED, manual, archive, setup, service and demo
- AED mode for automated ECG analysis and prompted treatment protocol for patients in cardiac arrest
- Manual mode for performing manual defibrillation, synchronized cardioversion, noninvasive pacing and ECG vital signs
- Archive mode for accessing stored patient information
- Setup mode for changing default settings of operating functions
- Service mode for authorized personnel to perform diagnostic tests and calibration
- Demo mode for simulated waveforms and trend graphs for demonstration and training purposes
- One (1) year warranty



DF028341

SKU	Description	Pkg	Price
DF028341	Reconditioned; 12-Lead with AED, Pacing, SPO2, SPCO, NIBP & ETCO2 (Version 2)	1 ea	\$23,995.00
DF028332	Reconditioned; 12-Lead with AED, Pacing, SPO2, SPCO, NIBP, ETCO2 & Bluetooth (Version 4)	1 ea	\$26,995.00



# QUICK SHIP DISPOSABLE MEDICAL PRODUCTS

In efforts to keep your facility stocked and prepared at all times, DiaMedical warehouses over 3,000 medical products available to ship immediately.

(877) 593-6011 | [DiaMedical.com](http://DiaMedical.com)







# SimVS EMS Virtual Diagnostic Platform

SimVS EMS Virtual Diagnostic Platform is a responsive simulation platform incorporating monitor, defibrillator and transport ventilator functions to promote student development of critical thinking skills in an ultrarealistic simulation setting.

SimVS enables instructors to manipulate patient vital signs and diagnostic results, empowering students to use their skills to react accordingly and treat the patient in a realistic manner. Additionally, SimVS allows instructors to convert any manikin or standardized patient into a high-fidelity simulator!

*\*No annual fees! Lifetime software upgrades*



For Training Use Only

SKU	Description	Pkg	Price
DG073218	SimVS EMS; For Training Purposes Only	1 ea	\$6,995.00
DG073214	SimVS EMS Plus; For Training Purposes Only	1 ea	\$7,995.00



Transport Ventilator



Patient Monitor



Zoll X or Lifepak 15 Defibrillator Monitor



Instructor Interface

## SimVS EMS Features



### Instructor Tablet

Controls the physiological displays of monitors empowering educators to influence student care decisions

### Realistic Student Monitors

Up to five (5) unique monitor interfaces including ultrarealistic Zoll X or Lifepak 15 Defibrillator, GE Patient Monitor, and telemetry transport ventilator

### Carrying Case

Durable bag that is configured to function identically to a defibrillator bag students will find throughout their training

### Scenario Booklet

Comprehensive simulations specifically designed to be used alongside the SimVS EMS Platform

### SimVS ECG Cables

Simulated ECG cables connect to SimVS bag to reinforce learning objectives required for starting an 12-lead ECG

### Training Accessories

Includes essential diagnostic equipment students need to provide a comprehensive patient exam

**AED Accessories & Cabinets**



**DiaMedical Indoor AED Cabinet with Alarm**

- Fits all AED brands
- Designed to be installed in an indoor location
- Dimensions: 15" H x 15" W x 8" D



SKU	Description	Pkg	Price
DF152050	Indoor AED Cabinet	1 ea	\$139.99

**Outdoor AED Cabinet with Alarm**

- Two (2) Gauze, 2" x 2"
- One (1) Suresite™ Window, 2.375" x 2.75"
- One (1) Alcohol Prep Pad
- One (1) CURAD Tape, 0.75" x 18", roll
- One (1) Tourniquet, Latex-Free



SKU	Description	Pkg	Price
EM090295	Outdoor AED Cabinet	1 ea	\$290.99

**Responder® Flexible AED Wall Sign**

- Designed to be mounted flat or can be converted to 3D
- Dimensions: 14.75" H x 9" W



SKU	Description	Pkg	Price
EM090296	AED Wall Sign	1 ea	\$14.95

**AED Check™ Tag**

- Made from weather-resistant plastic
- Includes reminder of four (4) important AED check components
- Space to record AED serial number, model number, electrode expiration date and battery installation date



SKU	Description	Pkg	Price
EM090297	AED Inspection Tag	1 ea	\$1.99
EM090298	AED Inspection Tag	5/pk	\$5.95

**Thermal Paper for Physio-Control Lifepak 15**

- For use with certain Physio-Control Lifepak 15 Defibrillators



SKU	Description	Pkg	Price
EM090299	100 mm; Plain White Thermal Paper	1 ea	\$6.99

**Thermal Paper for ZOLL® R & R/X Series**

- For use with ZOLL® R & R/X Series Monitor/Defibrillator



SKU	Description	Pkg	Price
EM090300	80 mm; Plain White Thermal Paper	6/pk	\$32.00
EM090301	80 mm; Thermal Paper with Grid	6/pk	\$32.00

**Symbio 12 Lead Arrhythmia Simulator**

- Simulate cardioversion by activating convert feature
- Indicator illuminates when simulators shocked by a defib set to 50J or more
- Connectors: V1, V2, V3, V4, V5, V6, RA, LA, LL, RL and PADS
- Paced beat is displayed and indicator blinks off when simulator is paced by an external pacer
- Use capture control to vary the pacing capture level
- Three (3) year warranty

SKU	Description	Pkg	Price
DF046221	Philips; For Training Use Only	1 ea	\$854.99
DF046311	Physio-Control; For Training Use Only	1 ea	\$854.99
DF048302	ZOLL®; For Training Use Only	1 ea	\$862.99



DF046311



## Automated Chest Compression Systems



### Defibtech Lifeline ARM Package

Lifeline ARM Automated Chest Compression (ACC) device performs mechanical chest compressions when effective manual CPR is not possible. In an emergency situation, the Lifeline ARM when used as an adjunct to manual CPR, provides extended CPR with minimal interruptions when fatigue, insufficient personnel or patient transport (e.g., ambulance, intra-hospital) may prohibit the delivery of effective/consistent compressions to the victim. Its award winning design provides portability, speedy deployment and delivers precise CPR with unmatched operating times.

- Single piece design provides increased structural integrity during compressions
- Simplified Control Panel requires just two steps to initiate mechanical CPR
- “Run Continuous” Setting performs chest compressions only without pauses for ventilation
- “Run with Breaths” Setting utilizes programmed pauses for rescue breaths
- Long runtime battery pack life with fast recharge time
- External AC power adapter also recharges the battery pack during use
- Single piece design enhances its usability and provides an incredibly strong platform for performing mechanical CPR
- Comes packed in a red canvas carrying case designed for backpack portability
- Manufactured to withstand small particulates and water spray
- Meets military standards for shock and vibration testing

SKU	Description	Pkg	Price
DF028329	Hands-Free CPR Device	1 ea	\$13,970.00

#### What's Included?

- Complete Lifeline Arm Hands-Free CPR Device
- Battery Pack
- Power Module
- Stabilization Strap
- Backboard
- Red canvas carrying case



DF028329

+ customer favorite



### Physio-Control LUCAS 3 Chest Compression System

LUCAS 3 Chest Compression System delivers high-performance, continuous chest compressions with less strain, micromanagement and risk for the caregiver. The LUCAS chest compression system, provides benefits both to the cardiac arrest patient and the resuscitation team. The LUCAS device extends the reach of care by maintaining chest compressions during transport to advanced lifesaving therapies, including ECMO or PCI in the cath lab. It can improve provider safety, avoiding fatigue over long durations and reducing transport risks by allowing caregivers to sit belted. And it can increase the efficiency of caregiving teams, calming the scene and providing an extra pair of hands so you can focus on the patient and underlying conditions.

- Deploys easily and quickly with a low-profile back plate for fast, smooth patient applications
- Ability to configure settings to meet specific EMS protocols
- Adjustable ventilation alerts, pause length and count
- Easily access and share post-event data
- Easy to store, clean and carry
- Ability to integrate with CODE-STAT 11 data review software
- A median of only seven (7) seconds of interruption demonstrated when transitioning from manual to mechanical CPR in clinical use
- Connects wirelessly to the LIFENET System
- Timer to remind rhythm and pulse checks
- Sets auto-lowering and pressure pad parameters to your preferences
- Straps secure patient arms and device during transport
- Spacious support structure accommodates larger patients

SKU	Description	Pkg	Price
DF152084	Hands-Free CPR Device	1 ea	\$20,015.00

#### What's Included?

- Complete LUCAS 3 Chest Compression System
- Battery Pack
- Carrying Case
- Patient Straps
- Back Plate



DF152084

**Electrodes**



**Medline MedGel®  
Pediatric Foam ECG Electrodes**

- Smaller footprint for pediatric patients
- High percentage of gel to electrode pad surface area allows for gentle use and removal
- Specially formulated high performance gel provides clear tracings
- Flexible leadwire attachment with perforated strip backing
- Latex-Free



**Medline MedGel®  
General Monitoring Foam Electrodes**

- Specially formulated high performance gel provides clear tracings
- Electrodes stay securely on the skin, yet remove cleanly with top-quality adhesive
- Flexible leadwire attachment with perforated strip backing
- Easy removal with full-edge lift tab
- Latex-Free, Fluid-Resistant Foam



SKU	Description	Pkg	Price
EM090302	Pediatric ECG Electrode	1000/ct	\$171.99

SKU	Description	Pkg	Price
DG028232	General Electrode	50/pkg	\$7.39



**Dynarex® Adult ECG Monitoring Electrodes**

- For ECG/EKG monitoring or diagnosis with the connected equipment as medical sensors
- Standard snap connector
- ECG Snap Electrodes are for adult use only
- Foam dimensions: 50 mm x 55 mm

SKU	Description	Pkg	Price
DG028459	Adult ECG Monitoring Electrodes	4/pkg	\$2.62
MED053023	Adult ECG Monitoring Electrodes	50/pkg	\$13.49



DG028459

**Glucose Monitoring**



**Arkray Assure®  
Platinum Blood Glucose Meter**

- LCD backlit display increases visibility in low-light situations
- Test strip ejector was engineered for touch-free test strip disposal
- 0.5 µL sample size, results in seven (7) seconds
- qcProGuard® Control solution test notification audio alert provides added quality control
- Features hypoglycemic alert for improved patient safety



**Arkray Assure® Lance Lockout  
Safety Lancets**

- Hidden-needle construction
- Automatic retraction after use
- Secured grip and handling
- Turn and pull protective cap



SKU	Description	Pkg	Price
MED160029	Glucose Meter	1 ea	\$48.99

SKU	Description	Pkg	Price
MED053566	21G	100/bx	\$20.99
DG028268	25G	100/bx	\$26.99
DG028268	30G	100/bx	\$26.99

**PARTS & ACCESSORIES**

MED160030	Blood Glucose Test Strips	50/bx	\$24.99
MED160032	Control Solutions, Normal and High	1 set	\$16.49



**Blood Pressure Monitors**



**ADC® Prosplyg™ 760 Pocket Aneroid Sphyg**

- Black enamel 300mmHg manometer provides dependability and exceeds industry standards
- Features Adcuff™ nylon cuff with ADC's proprietary Size Guide™ marking system which prevents mis-cuffing
- Includes leatherette carrying case with durable nylon zipper
- Latex-Free
- Additional colors available online



SKU	Description	Pkg	Price
DG010672	Infant; Black	1 ea	\$50.99
DG010673	Child; Black	1 ea	\$50.99
DG010670	Small Adult; Black	1 ea	\$50.99
DG010602	Adult; Black	1 ea	\$50.99
DG010659	Large Adult; Black	1 ea	\$61.99
DG010671	Thigh; Black	1 ea	\$68.99



**Pocket Nurse® BP Cuff**

- Nylon cuff with size-range finder for quick convenient fitting
- Latex-free PVC tubing, bladder and bulb
- Navy blue cuff with black zipper carrying case



SKU	Description	Pkg	Price
DG028489-ADLT	Adult	1 ea	\$21.99
DG028489-CHLD	Child	1 ea	\$21.99



**ADC® Advantage™ 6021N Automatic Digital BP Monitor**

- Quickly measures systolic, diastolic pressure and pulse
- Detects irregular heartbeat
- Smart Logic Technology automatically determines ideal inflation levels
- Internal memory recalls last ninety-nine (99) readings with date/time stamps for two (2) users
- Five (5) year warranty



SKU	Description	Pkg	Price
DG010611	Adult	1 ea	\$83.99

**Stethoscopes**



**Pocket Nurse® Dual-Head Teaching Stethoscope**

- Angled headset with soft ear tips
- Extended 53" tubing length
- Latex-free tubing resists stains
- 3-year warranty

SKU	Description	Pkg	Price
DG028383-BLKROY	Black/Royal Blue	1 ea	\$101.99
DG028383-BLKBURG	Black/Burgundy	1 ea	\$101.99





### ADC® Adscope® 603 Clinician Stethoscope

- Features Adsoft™ Plus silicone ear tips and a dual sided chestpiece with non-chill material for patient comfort
- Extra large bell (1 3/8" diameter) for unsurpassed low-frequency response
- Flexible 22" PVC tubing and stainless-steel binaurals fixed at 15° angle
- Limited lifetime warranty
- Additional colors available online



SKU	Description	Pkg	Price
DG090602-BLK	Black	1 ea	\$73.67



### 3M™ Littmann® Classic III™ Stethoscope

- Designed to be use with adult and pediatric patients
- Tunable, dual-sided stainless steel chestpiece with open or closed bell
- Non-chill bell sleeve for greater patient comfort
- Soft-sealing ear-tips provide optimal comfort and excellent sound occlusion
- Five (5) year warranty
- Additional colors available online



SKU	Description	Pkg	Price
DG250106	Black	1 ea	\$169.49

## Thermometers



### DiaMedical Diagnostics Oral-Axillary Digital Thermometers

- Measures temperatures from 90°-107.6° F, ±0.2°F
- Accurate readings in approximately thirty (30) seconds
- LCD display is easy to read with audible beep when measurement is complete



SKU	Description	Pkg	Price
COV012030	Digital Thermometer	10/pk	\$26.99



### DiaMedical Diagnostics Infrared Non-Contact Thermometer

- Measures temperatures from 89.6°-109.4°
- Accurate readings in less than two (2) seconds
- 10-reading memory



SKU	Description	Pkg	Price
COV012003	Non-Contact Thermometer	1 ea	\$52.99



### Amsino® Infrared Non-Contact Digital Thermometer

- Measurement range: 89.6° to 109° F (32-43° C)
- Temperature units: Celsius/Fahrenheit
- Measurement distance: 0-2" (0-5cm)
- 30-Second automatic shutoff function
- (2) AAA batteries (included)



SKU	Description	Pkg	Price
DG028306	Non-Contact Thermometer	1 ea	\$82.99

### Tempa-DOT Disposable Oral Thermometers

- Single-use oral thermometers prevent cross contamination and spread of infection
- Color changing dots indicate the patient's temperature
- Non-sterile and latex-free



SKU	Description	Pkg	Price
DG028476	Disposable Thermometers	100/bx	\$29.99

**Pulse Oximeters**



**Proactive Medical Protekt® Finger Pulse Oximeter**

- Displays digital SpO2 and pulse rate while using bar graph technology to monitor heart rate
- Portable non-invasive spot check oxygen saturation and heart rate
- Automatic power off function preserves battery life



SKU	Description	Pkg	Price
COV012004	Finger Pulse Oximeter	1 ea	\$46.50



**Pocket Nurse® Fingertip Pulse Oximeter**

- Reads SpO2 and pulse rate
- Easy to read OLED display
- Six (6) display modes
- Low voltage indicator



SKU	Description	Pkg	Price
RES091095	Finger Pulse Oximeter	1 ea	\$25.99



**ADC® Adimals® 2150 Fingertip Pulse Oximeter**

- Features child-friendly design on instrument face
- Displays digital SpO2 and pulse rate
- Includes six (6) changeable display modes and ten (10) brightness settings
- Accuracy: SpO2 ±2-3% / PR ±5 BPM
- Two (2) year warranty



SKU	Description	Pkg	Price
RS030601	Pediatric Finger Pulse Oximeter	1 ea	\$72.99

**IMDK® Fingertip Pulse Oximeter**

- Two (2) parameters: SpO2 and pulse rate
- Bright multi-directional OLED display
- Programmable audio alarms
- Automatic power off



SKU	Description	Pkg	Price
DG028123	Fingertip Pulse Oximeter	1 ea	\$47.99



**DiaMedical Diagnostics DM100 Handheld Pulse Oximeter**


- Features a large LCD display with SpO2 values and waveforms along with pulse rate and real time clock
- Latest 10 minutes trend graph and table of SpO2 and pulse rate can be reviewed on screen
- Data can be transferred to PC for storage, review and printing
- Two (2) year warranty

SKU	Description	Pkg	Price
DG250114	Handheld Pulse Oximeter	1 ea	\$362.99



DG250114

+ customer favorite



**CREATE AN ONLINE ACCOUNT ▶▶**

Utilize your DiaMedical Online Account to request quotes, manage orders, check tracking updates, and access invoices all in one location!



Diagnostic Bundles



**DiaMedical EMS Diagnostic Bundle #1**

Kit includes:

- 3M™ Littmann® Lightweight II S.E. Stethoscope
- ADC® Prosplyg™ 775 Pocket Aneroid Sphyg, Navy
- ADC® Metalite™ Reusable Penlight
- ADC Medicut™ EMT Shears - 7.25"
- ADC® Stethoscope Hip Clip Holster, Black
- 3M™ Littmann® Stethoscope Identification Tag, Black
- Additional colors available online



SKU	Description	Pkg	Price
EM090440-BL	Bundle #1, Black	1 ea	\$176.99
EM090440-CB	Bundle #1, Caribbean Blue	1 ea	\$176.99
EM090440-PRLP	Bundle #1, Pearl Pink	1 ea	\$176.99



EM090440-PRLP

+ customer favorite



**DiaMedical EMS Diagnostic Bundle #2**

Kit includes:

- 3M™ Littmann® Classic III™ Stethoscope
- ADC® Prosplyg™ 760 Pocket Aneroid Sphyg, Black
- ADC® Metalite™ Reusable Penlight
- ADC Medicut™ EMT Shears - 7.25"
- ADC® Stethoscope Hip Clip Holster, Black
- 3M™ Littmann® Stethoscope Identification Tag, Black
- Additional colors available online



SKU	Description	Pkg	Price
EM090441-BLK	Bundle #2, Black	1 ea	\$198.99
EM090441-LEL	Bundle #2, Lemon Lime	1 ea	\$198.99
EM090441-CRB	Bundle #2, Caribbean Blue	1 ea	\$198.99
EM090441-PRP	Bundle #2, Pearl Pink	1 ea	\$198.99



EM090441-CRB



**DiaMedical EMS Diagnostic Bundle #3**

Kit includes:

- 3M™ Littmann® Classic III™ Stethoscope
- ADC® Diagnostix™ 703 Palm Aneroid Sphyg, Black
- ADC® Metalite™ Reusable Penlight
- ADC Medicut™ EMT Shears - 7.25"
- ADC® Stethoscope Hip Clip Holster, Black
- 3M™ Littmann® Stethoscope Identification Tag, Black
- Additional colors available online



SKU	Description	Pkg	Price
EM090442-BLK	Bundle #3, Black	1 ea	\$234.95
EM090442-LEL	Bundle #3, Lemon Lime	1 ea	\$234.95
EM090442-CRB	Bundle #3, Caribbean Blue	1 ea	\$234.95
EM090442-PRP	Bundle #3, Pearl Pink	1 ea	\$234.95



EM090442-LEL





**IV & DRUG**  
DELIVERY

Vein Finders



**Veinlite® PED12®**

- Simplifies and improves venous access in children of all ages (0-17 years)
- Uses side-transillumination to make the veins more visible and accessible to quickly and accurately locate suitable veins for IV access in pediatric patients
- Switchable light colors for vein imaging, with 3 levels of adjustable brightness control
- Includes two (2) AAA batteries, twenty-five (25) disposable plastic covers and hardshell case

SKU	Description	Pkg	Price
DG250146	Pediatric Vein Finder	1 ea	\$299.00

PARTS & ACCESSORIES

DG170013	Disposable Plastic Covers	50/bx	\$18.00
----------	---------------------------	-------	---------



**Veinlite® EMS Pro**

- Side-transillumination technology reveals the structure of deep and superficial veins
- C-shaped opening provides traction and prevents the vein from rolling during venipuncture procedures
- Integrated white exam light for patient assessment in the field
- Includes two (2) AA batteries, fifty (50) disposable plastic covers, light shield, pediatric adapter and carrying case

SKU	Description	Pkg	Price
DG170006	Vein Finder	1 ea	\$379.00

PARTS & ACCESSORIES

DG170014	Disposable Plastic Covers	50/bx	\$24.00
----------	---------------------------	-------	---------



**Veinlite® LED+®**

- Side-transillumination technology reveals underlying venous structures, including veins, valves and bifurcations up to 6 mm deep
- Optimum high-contrast viewing area
- Effective on all skin tones and patients of all ages
- Includes rechargeable Li-ion battery with USB charging cable, fifty (50) disposable plastic covers, light shield, pediatric adapter and carrying case

SKU	Description	Pkg	Price
DG170005	Vein Finder	1 ea	\$549.00

PARTS & ACCESSORIES

DG170011	Disposable Plastic Covers	50/bx	\$24.00
----------	---------------------------	-------	---------



**Veinlite® LEDX®**

- Varicose vein mapping and access device designed specifically for sclerotherapy applications
- High contrast imaging and the largest vein visualization area, allowing for fast, efficient and safe sclerotherapy sessions
- C-shaped opening enables direct access for sclerosant injection
- Effective on all skin tones and patients of all ages
- Includes rechargeable Li-ion battery with USB charger, fifty (50) disposable plastic covers, light shield, size adapter and carrying case

SKU	Description	Pkg	Price
DG170004	Vein Finder	1 ea	\$679.00

PARTS & ACCESSORIES

DG170008	Disposable Plastic Covers	50/bx	\$24.00
----------	---------------------------	-------	---------

IV Start Kits & Supplies



Medline IV Start Kit with Chloraprep

- Two (2) Gauze - 2" x 2"
- One (1) Suresite™ Window - 2.375" x 2.75"
- One (1) Chloraprep - 1 mL
- One (1) CURAD Tape, Roll - 0.75" x 18"
- One (1) Tourniquet, Latex-Free



SKU	Description	Pkg	Price
SK020229	IV Start Kit with Chloraprep	100/cs	\$216.99



Dynarex® IV Start Kit with Vinyl Gloves

- One (1) Latex-Free Tourniquet
- One (1) Pair of Vinyl Gloves
- One (1) 0.75" x 18" Transparent Tape Roll
- One (1) Transparent IV Dressing
- One (1) Alcohol Prep Pad
- Two (2) 2" x 2" Gauze Sponges



SKU	Description	Pkg	Price
IV080435	IV Start Kit with Vinyl Gloves	1 ea	\$4.69



Medline IV Start Kit with Prevantics

- Two (2) Gauze - 2" x 2"
- One (1) Suresite™ Window - 2.375" x 2.75"
- One (1) Prevantics Prep Pad - 1 mL
- One (1) CURAD Tape, Roll - 0.75" x 18"
- One (1) Tourniquet, Latex-Free



SKU	Description	Pkg	Price
IV080434	IV Start Kit with Prevantics	1 ea	\$3.89



Dynarex® Alcohol Prep Pads

- Saturated with 70% isopropyl alcohol and gamma sterilized
- Non-woven pads are soft and absorbent
- Sterile, Single-Use



SKU	Description	Pkg	Price
SK020207-BX	Large	100/bx	\$4.19



Dynarex® Alcohol Swabsticks

- Each 4" swabstick is saturated with 70% isopropyl alcohol
- Soft sponge tip with plastic shaft
- Sterile, Single-Use



SKU	Description	Pkg	Price
MED053044	Alcohol Swabstick; 3/pk	25/bx	\$8.29



Pocket Nurse® Latex-Free Tourniquet

- Rolled and banded for clinician convenience
- Latex-Free, Single-Use



SKU	Description	Pkg	Price
PTC111223	Tourniquet, 1" x 18"; Blue	1 ea	\$0.89



Dynarex® View Guard Transparent Dressings

- Sterile, Latex-Free
- Provides a breathable barrier that lets air in while keeping moisture vapors and fluid out
- Use around IV sites, for wound care, and for surgical site protection
- Conforms to difficult body contours

SKU	Description	Pkg	Price
MED053190	2.38" x 2.75"	100/bx	\$39.99
MED053190-4IN	4" x 4.75"	50/bx	\$40.99



MED053190

IV Administration Sets



**Dynarex® IV Administration Sets**

- All IV sets are gravity sets; not compatible with infusion pumps
- Latex-Free, Non-DEHP Tubing, Individually Wrapped



SKU	Description	Pkg	Price
IV080011	83", 10 gtt/mL; One (1) Needless Port, One (1) Pre-Pierced Port	50/bx	\$115.95
IV080012	83", 60 gtt/mL; One (1) Needless Port, One (1) Pre-Pierced Port	50/bx	\$106.99



**Dynarex® IV Extension Set**

- Needle-free Luer connection for IVP
- For use with gravity administration sets only; not compatible with infusion pumps
- Sterile, Latex-Free, Individually Wrapped



SKU	Description	Pkg	Price
IV080455-6IN	6" IV Extension; Luer Lock	1 ea	\$4.59
IV080455-7IN	7" IV Extension; Luer Lock	1 ea	\$4.59
IV080455-8IN	8" IV Extension; Luer Lock	1 ea	\$4.59

IV080455



**ICU Medical 34" Secondary IV Set**

- Sterile, 34" Length
- Convertible piercing pin and sight chamber
- Secure lock mate
- Priming capacity: 7 mL




SKU	Description	Pkg	Price
IV080360	34" Secondary Macrobore Set	1 ea	\$3.69

IV080360



**Amsino AMSafe® Extension Set with Flow Controller**

- Priming Volume: 3.5 mL; Length: 20"
- Pre-Pierced Injection Y-Site: 1 for IVP
- Intended for use with blunt cannulas



SKU	Description	Pkg	Price
EM090319	20"; One (1) Pre-Pierced Port, Flow Controller, Slide Clamp	50/cs	\$279.95

EM090319



IV Catheters

**B. BRAUN**

**B. Braun Introcath Straight Safety IV Needles**



SKU	Description	Pkg	Price
MED053270	18G x 1.25"	50/bx	\$172.99
MED053268	20G x 1"	50/bx	\$157.99
MED053217	20G x 1.25"	50/bx	\$172.99
MED053269	22G x 1"	50/bx	\$172.99
MED041007	24G x 0.75"	50/bx	\$188.99

**BD**

**BD Insyte™ Autoguard™ Safety IV Needles**



SKU	Description	Pkg	Price
SK020174	14G x 1.75"	50/bx	\$161.99
SK020173	16G x 1.25"	50/bx	\$161.99
MED041016	18G x 1.25"	50/bx	\$191.57
MED041103	20G x 1"	50/bx	\$180.24
SK020170	22G x 1"	50/bx	\$183.15
SK020169	24G x .75"	50/bx	\$155.99

Syringes & Needles

**BD**

**BD Luer-Lock™ Tip Syringes without Needles**



SKU	Description	Pkg	Price
MED053236	1 mL	100/bx	\$80.90
MED053237	3 mL	200/bx	\$35.99
MED051540	5 mL	125/bx	\$33.99
MED053239	20 mL	48/bx	\$38.99
MED053203	30 mL	56/bx	\$51.99
MED051579	50 mL	40/bx	\$58.99

**BD**


**BD SafetyGlide™ Sliding Safety Hypodermic Needle**



SKU	Description	Pkg	Price
MED053222	21G x 1"	100/bx	\$33.99
MED053223	21G x 1.5"	100/bx	\$33.99
MED053283	23G x 1"	100/bx	\$25.99
MED053197	23G x 1.5"	100/bx	\$39.99
MED053282	25G x 0.625"	100/bx	\$25.99
MED053284	25G x 1"	100/bx	\$33.99

**BD**

**BD Eclipse™ Hypodermic Safety Needles**



SKU	Description	Pkg	Price
MED051042	18G x 1.5"	100/bx	\$46.99
MED053453	22G x 1.5"	100/bx	\$43.99
MED053452	23G x 1"	100/bx	\$46.99
MED053449	23G x 1.5"	100/bx	\$43.99
MED053450	25G x 0.625"	100/bx	\$44.99
MED053451	25G x 1"	100/bx	\$43.99
MED053584	25G x 1.5"	100/bx	\$43.99

**Monoject™ Angel™ Winged 12" Infusion Set**



SKU	Description	Pkg	Price
MED053101	21G x 0.75"	50/bx	\$81.99
MED053102	23G x 0.75"	50/bx	\$81.99

Infusion Pumps & Tubing

**Baxter**

**Baxter® Sigma Spectrum IV Pump**

- Expertly Reconditioned Baxter Sigma Spectrum Pump can infuse IV, subQ, intrathecal and epidural fluids and medications with an integrated drug library
- Lightweight, compact design with a bright, easy to read color display and quick drug library search
- Includes Primary, Secondary, Multi-Step and Cyclic TPN therapy modes
- Compatible with all standard Baxter tubing sets
- Infusions can be easily switched to and from gravity and pump; No separate infusion pump tubing is required
- One (1) year warranty

\*Medication formularies may vary. Facility is responsible for updating drug library on pump.



IV011410

SKU	Description	Pkg	Price
IV011410	Reconditioned IV Pump; For Training Use Only	1 ea	\$2,357.00
IV011410-PR	Reconditioned IV Pump; Patient Ready*	1 ea	\$2,357.00



**Baxter**

**Baxter® ClearLink™ Continu-Flo Solution Set**

- MFR #: 2C8519
- Luer Ports: Two (2) ClearLink™ for IVP and IVPB

SKU	Description	Pkg	Price
IV080327	112", 10 gtt/mL; Two (2) Needle-Free Ports	1 ea	\$8.09

**LifeFlow**

**LifeFlow® PLUS Blood & Fluid Infuser**

- Hand-operated, rapid infusion device delivers fluid and blood products quickly and effectively
- Deliver 500 mL of fluid in less than 2 minutes through a 20G IV catheter
- Deliver one (1) unit of whole blood in two (2) minutes
- Handle and tubing packaged together for quick and easy set-up



Item requires Medical Device Authorization for purchase

SKU	Description	Pkg	Price
EM090344	Rapid Infusion Device	1 ea	\$459.95

IO Access Systems

**SAM MEDICAL**

**SAM Medical® IO Access System**

- Manually operated intraosseous (IO) access system
- Intuitively designed for current IO users
- Features color coding and depth indicators matching industry standards



SKU	Description	Pkg	Price
EM090309	IO Driver	1 ea	\$184.99

PARTS & ACCESSORIES

EM090310	15 mm IO Needle; Pink	1 ea	\$129.99
EM090311	15 mm IO Needle; Pink	5/bx	\$599.99
EM090312	25 mm IO Needle; Blue	1 ea	\$129.99
EM090313	25 mm IO Needle; Blue	5/bx	\$599.99
EM090314	45 mm IO Needle; Yellow	1 ea	\$129.99
EM090315	45 mm IO Needle; Yellow	5/bx	\$599.99

**SAFEGUARD MEDICAL**

**Safeguard Medical™ NIO® Intraosseous Device**

- Pre-set needle penetration depth based on patient age
- Pocket-sized, lightweight device is pre-assembled and ready to use
- Self-contained with no exposed needle
- Sterile, disposable, single-use



SKU	Description	Pkg	Price
EM090318	18G Stepped Needle®, Infant; Pink	1 ea	\$129.00
EM090317	18G, Pediatric; Red	1 ea	\$129.00
EM090316	15G, Adult; Blue	1 ea	\$129.00



**TRAINING &  
SIMULATION**



simlabsolutions

# SimLabSolutions Loaded BLS Jump Bag™

Simulate basic life support skills with the SimLabSolutions Loaded BLS Jump Bag™. The Loaded BLS Jump Bag™ arrives simulation ready with the supplies required to teach first aid and basic life support skills.



For Training Use Only

SKU	Description	Pkg	Price
LC129912	SimLabSolutions Loaded BLS Jump Bag™	1 ea	\$325.99
LC129922	SimLabSolutions Loaded BLS Jump Bag™ Refill Kit	1 ea	\$270.99

TUESDAY TEACHINGS  
SEE THIS PRODUCT  
IN ACTION ▶▶▶▶▶



### PPE/Infection Control

- Two (2) Pairs Sterile Exam Gloves, Vinyl
- Nine (9) Alcohol Gel Hand Wipes



### Wound Care/Dressings

- One Hundred (100) Plastic Adhesive Strips - 0.75" x 3"
- One (1) Elastic Bandage - 3" x 5 yd
- Two (2) Sterile Trauma Dressings - 12" x 30"
- Two (2) Sterile Burn Dressings
- Two (2) Triangular Bandages
- Five (5) Sterile ABD Pads - 5" x 9"
- Two (2) Sterile Transparent Dressings - 4" x 4"
- Fifteen (15) Sterile Tongue Depressors
- Two (2) Sterile Sof-Form Gauze Rolls - 1"
- Three (3) Sterile Sof-Form Gauze Rolls - 3"
- Three (3) Sterile Gauze - 4" x 4"
- Four (4) Sterile Eye Pads
- Two (2) Transparent Adhesive Tape Rolls - 1"



### Diagnostic

- One (1) Glucose Meter
- Fifty (50) Test Strips
- Fifty (50) Safety Lancets
- One (1) Aneroid BP Cuff, Adult
- One (1) Stethoscope
- One (1) Penlight with Pupil Gauge
- One (1) Digital Thermometer
- One (1) Pulse Oximeter, Fingertip



### Respiratory

- One (1) BVM, Adult
- One (1) Disposable CPR Face Shield
- One (1) Each Disposable Guedel Airway - 50 mm, 60 mm, 70 mm, 80 mm, 90 mm & 100 mm
- One (1) Each Rusch Nasal Airway - 16 FR, 26 FR & 28 FR
- Six (6) Sterile Lubricating Jelly Packets
- One (1) Bite Stick



### Additional Supplies

- One (1) Splinter Forceps
- One (1) Utility Scissors - 7.5"
- One (1) Kelly Forceps, Straight - 5.5"
- One (1) Instant Cold Pack - 5" x 9"
- One (1) Obstetrical Kit
- One (1) Emergency Survival Blanket
- Twenty-Five (25) Alcohol Prep Pads
- One (1) Each Cervical Collars, Adult & Pediatric



simlabsolutions

# SimLabSolutions Loaded ALS Jump Bag™

Stay life support ready with the SimLabSolutions Loaded ALS Jump Bag™! This loaded jump bag arrives stocked with airway supplies, diagnostic equipment, IV start supplies and first aid supplies required for advanced life support.



For Training Use Only

SKU	Description	Pkg	Price
LC129911	SimLabSolutions Loaded ALS Jump Bag™	1 ea	\$745.99
LC129921	SimLabSolutions Loaded ALS Jump Bag™ Refill Kit	1 ea	\$595.99

TUESDAY TEACHINGS  
SEE THIS PRODUCT  
IN ACTION ▶▶▶▶



### Diagnostic

- One (1) Glucose Meter
- Fifty (50) Test Strips
- Fifty (50) Safety Lancets
- One (1) Aneroid BP Cuff, Adult
- One (1) Stethoscope
- One (1) Penlight with Pupil Gauge
- One (1) Digital Thermometer
- One (1) Pulse Oximeter, Fingertip



### PPE/Infection Control

- Two (2) Pairs Sterile Exam Gloves, Vinyl
- Nine (9) Alcohol Gel Hand Wipes
- One (1) Procedure Mask with Face Shield



### Infusion/Injection

- One (1) SimRx 0.9% NaCl Simulated IV Bag - 1,000 mL
- Twenty-Five (25) Alcohol Prep Pads
- Three (3) Latex-Free Tourniquets, Rolled - 1" x 18"
- One (1) IV Armboard, Adult - 3" x 9"
- Two (2) Each Luer Lock Syringes - 3mL, 5mL & 10 mL
- Two (2) Tuberculin Safety Syringes
- Four (4) Each Safety Hypodermic Needles - 25G & 18G
- Four (4) IV Start Kits with Chloraprep
- One (1) IV Extension Set - 7"
- One (1) 3-Way Stopcock
- Two (2) Primary IV Sets - Microdrip & Macro drip
- Twenty (20) Alcohol Swabsticks
- One (1) Each Safety IV Catheters - 24G, 22G, 18G & 14G



### Respiratory

- One (1) BVM, Adult
- One (1) Intubation Bag
- One (1) Emergency Cricothyrotomy Kit
- One (1) Each Nasal Cannulas - Adult & Pediatric
- One (1) Each Non-Rebreather Oxygen Masks - Adult & Pediatric
- Six (6) Sterile Lubricating Jelly Packets
- One (1) Bite Stick
- Five (5) Sterile Tongue Depressors
- One (1) Each Disposable Guedel Airway - 50 mm, 60 mm, 70 mm, 80 mm, 90 mm & 100 mm
- One (1) Each Rusch Nasal Airway - 16 FR, 26 FR, 28 FR & 30 FR
- One (1) Disposable Laryngoscope Handle
- One (1) Each Disposable Laryngoscope Blades, Miller - #0, #1 & #2
- One (1) Each Disposable Laryngoscope Blades, Macintosh - #2, #3 & #4
- One (1) Each Cuffed Endotracheal Tubes - 3.0 mm, 4.0 mm, 5.0 mm, 6.0 mm, 7.0 mm, 8.0 mm & 9.0 mm
- One (1) Each Intubation Stylets - 10 FR & 14 FR
- One (1) CO2 Indicator, Colorimetric
- One (1) ET Tube Holder
- One (1) Magill Catheter Forceps, Adult - 9.75"



### Wound Care/Dressings

- One Hundred (100) Plastic Adhesive Strips - 0.75" x 3"
- Two (2) Sterile Trauma Dressings - 12" x 30"
- Two (2) Sterile Sof-Form Gauze Rolls - 1"
- Three (3) Sterile Sof-Form Gauze Rolls - 3"
- Three (3) Sterile Gauze - 4" x 4"
- Five (5) Sterile Gauze - 2" x 2"
- Two (2) Sterile Burn Dressings
- Five (5) Sterile ABD Pads - 5" x 9"
- Four (4) Sterile Eye Pads
- Two (2) Sterile Transparent Dressings - 4" x 4"
- One (1) Elastic Bandage - 3" x 5 yd
- Two (2) Triangular Bandages



### Additional Supplies

- Two (2) Transparent Adhesive Tape Rolls - 1"
- One (1) Splinter Forceps
- One (1) Utility Scissors - 7.5"
- One (1) Kelly Forceps, Straight - 5.5"
- One (1) Instant Cold Pack - 5" x 9"
- One (1) Obstetrical Kit
- One (1) Emergency Survival Blanket
- One (1) Each Cervical Collars - Adult & Pediatric

simlabsolutions

# SimLabSolutions Loaded Pediatric Jump Bag™

Are your students trained on how to respond to a pediatric emergency? By incorporating SimLabSolutions' Loaded Pediatric Jump Bag™ into your curricula, instructors can simulate pediatric emergencies that incorporate the Broselow System, pediatric respiratory care and medication dosage calculations.



For Training Use Only

SKU	Description	Pkg	Price
LC129913	SimLabSolutions Loaded Pediatric Jump Bag™	1 ea	\$755.99
LC129923	SimLabSolutions Loaded Pediatric Jump Bag™ Refill Kit	1 ea	\$605.99

TUESDAY TEACHINGS  
SEE THIS PRODUCT  
IN ACTION ▶▶▶▶



## Color-Coded Weight Based Modules

- One (1) Suction Catheter Kit
- One (1) Cuffed Murphy Endotracheal Tube
- One (1) Disposable Laryngoscope Blade, Miller and/or Macintosh
- One (1) Intubating Stylet
- One (1) Disposable Guedel Airway
- One (1) Rusch Nasal Airway
- One (1) Safety IV Catheter
- One (1) Butterfly Infusion Set
- Two (2) IV Start Kits
- Ten (10) Alcohol Prep Pads
- Three (3) Sterile Lubricating Jelly Packets



## Additional Supplies

- One (1) SimRx 0.9% NaCl Simulated IV Bag - 250 mL
- One (1) SimRx 0.9% NaCl Simulated IV Bag - 500 mL
- One (1) Each Disposable IV Armboards - Neonatal, Infant & Adult
- Two (2) Each Luer Lock Syringes - 10 mL & 3 mL
- Four (4) Pediatric/Neonatal ECG Electrodes
- One (1) Each BVMs - Child & Infant
- One (1) Each Aneroid BP Cuffs - Pediatric & Infant
- One (1) Pediatric Stethoscope
- Two (2) Sterile Bulb Syringes - 60 mL
- One (1) Obstetrical Kit
- One (1) PediaQuick Tape
- One (1) Each Primary IV Sets - Microdrip & Macrodrip
- One (1) IV Extension Set - 7"
- One (1) Disposable Laryngoscope Handle
- One (1) Whole Blood Tube Plastic Vacutainer
- One (1) 3-Way Stopcock
- One (1) Magill Catheter Forceps, Child
- One (1) Each Nasal Cannulas - Pediatric & Infant
- One (1) Non-Rebreather Oxygen Mask - Pediatric
- One (1) CO2 Indicator, Colorimetric

## TUESDAY TEACHINGS

### Bring Our Products Into Your Practice

Learn why our EMS Loaded Jump Bags™ are an essential tool for life-saving procedures & training scenarios!

[WATCH NOW](#) ▶▶▶▶



simlabsolutions

# SimLabSolutions Loaded Simulated Medication Box™

Prepare your soon-to-be paramedics for pre-hospital medication administration with SimLabSolutions' Loaded Simulated Medication Box™. This medication box arrives loaded with the most commonly used pre-hospital simulated medications and IV fluids.



For Training Use Only

SKU	Description	Pkg	Price
LC129910	SimLabSolutions Loaded Simulated Medication Box™	1 ea	\$820.99
LC129920	SimLabSolutions Loaded Simulated Medication Box Refill Kit™	1 ea	\$535.99



### Oral Medication

- One (1) Simulated Acetaminophen Bottle - 600 mg
- One (1) Simulated Aspirin Bottle - 325 mg
- One (1) Simulated Prednisone Bottle
- One (1) Simulated Nitroglycerin Bottle



### Respiratory Medications

- Six (6) Simulated Albuterol Nebules
- One (1) Nebulizer



### Pre-Filled Syringes & IV Solutions

- Three (3) SimRx Atropine Sulfate Injection USP - 10 mL
- Four (4) SimRx Epinephrine Injection - 10 mL
- One (1) SimRx Adenosine 6 mg Injection USP - 2 mL
- One (1) SimRx Adenosine 12 mg Injection USP - 4 mL
- One (1) SimRx Sodium Bicarbonate 8.4% Injection USP - 50 mL
- Two (2) SimRx 50% Dextrose Injection USP - 50 mL
- One (1) SimRx Lidocaine HCl 2% Injection USP - 5 mL
- One (1) SimRx Calcium Chloride 10% Injection USP - 10 mL
- Two (2) SimRx Naloxone Hydrochloride Injection USP - 2 mL
- Ten (10) SimRx 0.9% Saline Flush Syringes - 10 mL



### Medication Vials

- Three (3) Simulated Fentanyl Ampules - 5 mL
- Two (2) Simulated Midazolam Vials
- Two (2) Simulated Diphenhydramine Vials
- Two (2) Simulated Ondansetron Vials
- Two (2) Simulated Furosemide Vials
- Two (2) Simulated Vials - 20 mL



### Label/Vial Sets

- Three (3) Simulated Amiodarone HCl Peel & Stick Labels - 3 mL
- Three (3) Simulated Tint Vials - 3 mL
- Two (2) Simulated Dopamine HCl Peel & Stick Labels - 5 mL
- Two (2) Simulated Vials - 5 mL
- Four (4) Simulated Magnesium Sulfate Peel & Stick Labels - 2 mL
- Four (4) Simulated Tint Vials - 2 mL
- Five (5) Simulated Morphine Peel & Stick Labels - 10 mL
- Five (5) Simulated Mini Ampules - 1 mL
- Four (4) Simulated Epinephrine Peel & Stick Labels - 1 mL
- Four (4) Simulated Mini Ampules - 1 mL



### Injection Supplies

- Thirty (30) Alcohol Prep Pads
- Five (5) Simulated IV Solution Additive Labels
- Two (2) Primary IV Sets - 10 mL
- Twenty (20) Nasal Cannulas
- Two (2) 3-Way Stopcock
- Two (2) Tuberculin Safety Syringe
- Two (2) Each Luer Lock Syringes - 3 mL, 5 mL, 10 mL, 20 mL
- Three (3) Each Safety Hypodermic Needles - 25G x 1", 22G x 1", 18G x 1.5"
- One (1) Transparent Adhesive Tape Roll - 1"



### Additional Medications & Supplies

- Two (2) SimRx 0.9% NaCl Simulated IV Bags - 250 mL
- Six (6) Simulated Albuterol Sulfate Nebules
- One (1) Sharps Container
- One (1) M.A.D. Intranasal Atomizer

simlabsolutions

# SimLabSolutions Emergency Opioid Overdose Training Kit™

Are you prepared for an opioid overdose? The SimLabSolutions Emergency Opioid Overdose Training Kit™ provides users with simulated Naloxone and other basic PPE to practice proper intranasal or injection administration of Naloxone.



For Training Use Only



SKU	Description	Pkg	Price
LC139910	Syringe, Single Use	1 ea	\$35.99
LC139911	Syringe, 10 Uses	10/pk	\$95.99
LC139912	Syringe, 20 Uses	20/pk	\$160.99

SKU	Description	Pkg	Price
LC139920	Vial, Single Use	1 ea	\$35.99
LC139921	Vial, 10 Uses	10/pk	\$95.99
LC139922	Vial, 20 Uses	20/pk	\$160.99



## Simulated Naloxone Syringe Opioid Overdose Training Kit™

- One (1), Ten (10), or Twenty (20) Pre-Filled SimRx Simulated Naloxone Syringe
- One (1), Ten (10), or Twenty (20) Pair Non-Latex Gloves
- One (1) or Two (2) M.A.D. Intranasal Atomizer
- One (1) or Two (2) Pull-Through Breakaway Seal
- One (1) CPR Face Shield
- One (1) Reusable Bag



## Simulated Naloxone Vial Opioid Overdose Training Kit™

- One (1), Ten (10), or Twenty (20) SimRx Simulated Naloxone Vial
- One (1), Ten (10), or Twenty (20) Syringe with Needle - 3 mL
- One (1), Ten (10), or Twenty (20) Pair Non-Latex Gloves
- One (1) or Two (2) M.A.D. Intranasal Atomizer
- One (1) or Two (2) Pull-Through Breakaway Seal
- One (1) CPR Face Shield
- One (1) Reusable Bag



SimRig® The Ambulance Trainer



SimRig® The Ambulance Trainer provides an unparalleled level of realism and precision to EMS training, allowing students to enter the workforce with a complete understanding of how to operate within an ambulance.

FEATURES FOR BOTH CLASSROOM & MOBILE SIMRIG®

Functioning action control station

Flip-Up Bumper built to the specifications of an actual ambulance

Integrated cot mounting apparatus

Functioning suction and simulated oxygen

Over thirty (30) additional standard features

*\*Lifetime structural warranty and five (5) year limited warranty on all other components*



**Established Brand**

DiaMedical is at the forefront of simulation training, providing tools for training and simulation.



**Facilitates Group Learning**

SimRig's® open design and HD camera system facilitates instructor and student observation.



**Unparalleled Realism**

SimRig® is built to precisely mimic a real ambulance. All aspects, including the back door and bumper, are built to actual specifications of a road ready ambulance.



**Scientifically Backed**

Studies have shown that hands-on simulation training can help enhance performance and reduce errors in critical life-saving emergencies.



SRM00220



## SimRig® The Classroom Ambulance Trainer

SRM00220

Teach with an ambulance in your classroom! The unique modular design allows for the SimRig® to be installed in any classroom or training center! The SimRig® Classroom Ambulance Trainer is built to precisely mimic a road-ready ambulance to help prepare the students for the hardships of working within a confined space.

**CALL (877) 593-6011 FOR PRICING**



SRM00261



## SimRig® The Mobile Ambulance Trainer

SRM00261

Simulate the entire call process anywhere! SimRig® The Mobile Ambulance Trainer is designed to be transported behind any trailer compatible vehicle to different training sites and arrives with standard trailer hitch and functioning electrical components. By being completely mobile, instructors can simulate every aspect of an emergency call from a parking lot, shopping center, simulation site or even streetside.

**CALL (877) 593-6011 FOR PRICING**

### PARTS & ACCESSORIES FOR BOTH CLASSROOM & MOBILE SIMRIG®

SKU	Description
SRM00231	HD 4-Color Camera System In Ceiling
SRM00230	46" HD Color Monitor
SRM00237	SimRig Rollok Rear Security Door
SRM00238	SimRig Rollok Side Security Door
SRM00223	Observation Window
SRM00232	Siren Package
MS099232	2-Way Radio
MS099211	Akro Storage Bin Containers
MS099221	SimRig® Disposable Supply Package



**CONTACT YOUR ACCOUNT MANAGER TO CUSTOMIZE YOUR SIMRIG®**

(877) 593-6011 | [SimRig@DiaMedicalUSA.com](mailto:SimRig@DiaMedicalUSA.com)



# SIM RESCUE™

THE AUTO EXTRICATION SYSTEM



## Why Educators Love SimRescue™ The Auto Extrication System

- **Simulates Dozens of Rollover Positions:**

SimRescue™ is constructed with strength-rated steel construction so the system can be safely setup in numerous positions for unique extrication scenarios.

- **Easy Transportation:**

SimRescue™ easily separates in half for transportation and storage. When disassembled, SimRescue™ can fit through any standard size doorway.

## SimRescue™ The Auto Extrication System



Retractable steering wheel with full dash and airbags



Tempered glass breakage simulator, pry box simulator and seat belt cutting simulator



Removable Lexan windows, adjustable seats and foot pedals with foot pinning simulator



SimRescue™ Starter Package:

- Twenty (20) seat belts
- Twenty (20) tempered glass pieces
- Twenty (20) pry plate pieces
- One (1) crowbar
- One (1) window punch
- Tools for assembly





**SimRescue™ Facilitates  
Dozens of Simple &  
Complex Rollover  
Positions.**



## SimRescue™ The Auto Extrication System

EM012818

Emergency Medical Service providers and paramedics provide patient care in a variety of settings from calm to chaotic. However, when it comes to responding to an emergent motor vehicle accident, every second matters. With SimRescue™ The Auto Extrication System, students and professionals can safely practice patient removal scenarios! SimRescue™ is the only auto extrication system on the market that was specifically created for patient extrication. From basic seatbelt cutting simulations to lethal roll-over accidents, SimRescue™ allows EMS professionals to feel confident when responding to an auto extrication call.

**CALL (877) 593-6011 FOR PRICING**



EM012818

### PARTS & ACCESSORIES

EM090148	10' Replacement Seatbelt
EM090149	Pry Plate
EM090150	5" x 5" Tempered Glass



Scan to Watch  
SimRescue™ in Action!




Patient Simulator Manikins



# SIMBODIES® Trauma Care EMS-Man Manikin

- Anatomically correct airway, chest and pelvis
- Blood-filled humerus and tibia sites for IO access
- Adult male approximately 30-40 years of age
- Skin tone: White Caucasian
- Clothing size: Medium/Large
- Shoe size: US 11, UK 10, ER 45
- Splashproof, washable skin
- Lifelike look and feel with flexible limbs
- Wig is made from real human hair



MS150294

SKU	Description	Pkg	Price
MS150294	Trauma Care Manikin; Full Body	1 ea	\$32,000.00

## EMS-Man Manikin Scenario Training



### Airway Management

- Head tilt, chin lift, jaw thrust
- Oral and nasopharyngeal airway suctioning
- Insertion of oral, supraglottic and nasopharyngeal airways
- Laryngoscopy and endotracheal intubation
- BVM ventilation
- Needle thoracentesis with chest seal application



### Immobilization

- Application of spinal collar, pelvic binder and traction devices
- Manual in-line stabilization



### Hemorrhage Management

- Thigh/groin wound packing with application of direct and indirect pressure
- Tourniquet application for hemorrhage control
- Removable lower leg amputation with bleeding capability
- Simulated circulation assessment



### Other Features

- IO access in bilateral humeral heads and tibias
- Chest compressions
- Application of AED/defibrillator pads
- Patient transfers, dressing, bed baths



## SIMBODIES® Emergency Medical Skills Torso Trainer

- Skin safe, prosthetic grade platinum silicone
- Anatomically correct airway for intubation wand ventilation
- Surgical airway (cricothyrotomy) capability
- Nasal passages for nasopharyngeal airway insertion and suctioning
- Humeral heads for IO access and infusion
- Needle decompression in the 2nd intercostal space
- Chest drain in the 6th intercostal space
- Wig is made from real human hair
- Weight: 55.1 lbs
- Dimensions: 31" x 23.6" x 11"



MS150296-LGT

SKU	Description	Pkg	Price
MS150296-LGT	Skills Trainer; Torso - Light Skin	1 ea	\$21,000.00
MS150296-MED	Skills Trainer; Torso - Medium Skin	1 ea	\$21,000.00
MS150296-DRK	Skills Trainer; Torso - Dark Skin	1 ea	\$21,000.00


CAE

# CAE Ares Emergency Care Manikin

CAE Ares is mid-fidelity manikin designed to fulfill the requirements for ALS, ACLS and emergency care team training. Ares is portable, lightweight and offers a variety of key features needed for emergency care simulation including realistic articulation, bilateral and unilateral chest rise and ECG monitoring.

**Standard Equipment:**

- CAE Ares wireless and tetherless manikin
- Male and female gender configuration
- Physiological Modeling for CAE Maestro
- One wireless StethoSym
- Instructor tablet
- CAE Maestro software suite—instructor-driven (manual mode) only
- One CAE Maestro Standalone license (manual mode)
- Six (6) CAE simulated clinical experiences (SCEs)
- One (1) year CAE assurance warranty



SKU	Description	Pkg	Price
MS140087	Advanced Manikin; Medium Skin Tone	1 ea	\$25,390.00
MS140087-D	Advanced Manikin; Dark Skin Tone	1 ea	\$25,390.00
MS140088	Complete Manikin; Medium Skin Tone	1 ea	\$28,940.00
MS140088-D	Complete Manikin; Dark Skin Tone	1 ea	\$28,940.00

+ customer favorite

## Key Features

**Neuro:** Tristate pupils and optional SymEyes

**Articulation:** Realistic articulation, range of motion in hips, knees, ankles and shoulders

**Airway:** Bag-valve-mask ventilation, orotracheal and nasotracheal intubation, placement of various airway adjuncts, endotracheal tube intubation, retrograde and fiberoptic intubation, transtracheal jet ventilation, needle and surgical cricothyrotomy, tracheostomy, right main stem intubation with detection, stomach distention with esophageal intubation and laryngospasm (manual)

**Breathing:** Spontaneous breathing, bilateral and unilateral chest rise and fall, lung auscultation, upper airway sounds and bilateral needle decompression

**IM Medication Administration:** Four (4) intramuscular sites

**Urinary:** Urinary catheterization without fluids and interchangeable genitalia

**Cardiac:** 4-Lead ECG monitoring with real equipment, 12-lead dynamic ECG display, defibrillation, cardioversion and pacing using live equipment via external defib box (optional)

**Circulation:** Unilateral blood pressure measurement by auscultation and palpation, bilateral carotid and femoral pulses, unilateral radial and brachial pulses with variable pulse strength

**CPR:** CPR compressions generate palpable pulse, blood pressure, waveform and ECG artifacts, realistic chest compression depth and resistance, software metrics detect hand placement, rate and depth of compressions, recoil, ventilators and ventilation

	Base	Advanced	Complete
Manikin	✓	✓	✓
StethoSym	✓	✓	✓
SymEyes	✗	✓	✓
SymDefib	Add-On	✓	✓
One Maestro Standalone License	✓	✓	✓
Modeled Physiology Software	Add-On	Add-On	✓
IO Legs	Add-On	Add-On	✓
Accessories Available	✓	✓	✓



### Simulaids STAT Manikin

- Automatic cardioversion upon defibrillation
- Seventeen (17) rhythms and pacing capability
- Bilateral chest tube insertion and decompression sites
- Twelve (12) pulse point with three (3) zones used to demonstrate BP loss
- Venous access from dorsal hand to upper arm with all major veins
- Upper teeth break out with laryngoscope pressure
- Tongue swelling and laryngospasm
- Removable cricoid membrane for cricothyrotomy
- Insertion of ET tubes, Combitube®, LMA, LT, OPA, NPA and ET devices



MS154010

SKU	Description	Pkg	Price
MS154010	Manikin with Deluxe Airway Head	1 ea	\$9,117.95



### Simulaids Stat Baby Simulator

- Bilateral chest tube insertion
- Defibrillation with shock level of 2-4 J/kg
- ET tube insertion
- Feeding tube care
- Foley catheter insertion (female)
- Intramuscular and subcutaneous injection sites
- Intraosseous (IO) leg
- LMA insertion
- NG tube placement fluid reservoir
- Peripheral IV site on hand and foot
- PICC line
- Pneumothorax
- Rectal medication
- Six (6) Pulse locations - two (2) brachial, two (2) femoral and two (2) carotid
- Tongue edema and laryngospasm
- Three (3) year warranty



MS153976

SKU	Description	Pkg	Price
MS153976	Infant Simulator; 9 Month Old	1 ea	\$6,699.95



### Life/form® Complete Infant CRiSis™ Manikin

The Life/form® Complete Infant CRiSis™ Manikin features a complete resuscitation system designed for teaching life-saving techniques for infants. This manikin allows practice of 4-lead ECG monitoring, IV therapy, umbilical catheterization, airway management, intraosseous (IO) infusion and CPR.

- Designed to offer realistic training and practice of infant airway management and CPR education
- Practice 4-lead ECG monitoring, IV therapy, umbilical catheterization, airway management, intraosseous infusion and CPR
- Airway suctioning and intubation using up to 4.0 mm ID endotracheal tubes
- Chest skin enables 4-lead ECG monitoring with simulated patient monitoring systems
- Five (5) year warranty



MS200034

SKU	Description	Pkg	Price
MS200034	Infant PALS Manikin; Full Body	1 ea	\$2,403.95

**Life/form.**

**Life/form® Deluxe Child CRiSis™ Manikin with ECG Simulator**

The Life/form® Deluxe Child CRiSis™ Manikin is designed for teaching individuals or groups life-saving techniques for children. This manikin provides complete PALS training capabilities it combines all of the skill stations on one manikin, making training more realistic and giving students an appreciation for the smaller scale of the patient.

- Enables BVM ventilation, oral/nasal intubation and 4-lead ECG Simulator recreates seventeen (17) adult heart rhythms and seventeen (17) pediatric rhythms
- Oral, nasal, and digital intubation with ET, EOA, PTL, LMA, EGTA, Combitube and King Video
- Includes electronic BP arm, femoral access, IO infusion and injection arm
- Five (5) year warranty



MS140034

SKU	Description	Pkg	Price
MS140034	Pediatric PALS Manikin, Full Body; 5 Year Old	1 ea	\$6,252.95

**Rescue Manikins**



**Simulaids Rescue Randy**

- Rescue Randy was developed for lifelike adult victim handling, transportation and extrication training
- Can be used in situations too hazardous or uncomfortable for human volunteers
- Made of durable vinyl with zinc-plated steel frame for longevity and rugged use
- Fully articulated for positioning in widely diversified training environments
- Three (3) year limited warranty



MS153977

SKU	Description	Pkg	Price
MS153977	55 lbs; 5' 5"	1 ea	\$1,239.95
MS153974	105 lbs; 5' 5"	1 ea	\$1,361.95
MS150183	125 lbs; 5' 5"	1 ea	\$1,736.95
MS150109	135 lbs; 5' 5"	1 ea	\$2,166.95
MS200002	145 lbs; 5' 5"	1 ea	\$1,736.95
MS150111	160 lbs; 5' 5"	1 ea	\$2,102.95
MS153975	165 lbs; 5' 5"	1 ea	\$1,864.95



**Simulaids Trauma Randy Upgrade Flash Moulage Overlay**

- Four (4) bleeding moulages
- One (1) deep bleeding packing wound
- Three (3) shallow bleeding wounds
- One (1) third degree burn
- One (1) flushing syringe to clear lines
- One (1) gallon of TrueClot® Blood Simulant Concentrate
- TrueClot® Simulant Hemostatic Gauze
- 500 mL TrueClot® Simulant Clotting Solution
- Foot pump with gravity tank



MS153999

SKU	Description	Pkg	Price
MS153999	Moulage Overlay for Manikin	1 ea	\$1,519.95





### Simulaids Rescue Rudy

- Rescue Rudy was developed for lifelike adult victim handling, transportation and extrication training
- Can be used in situations too hazardous or uncomfortable for human volunteers
- Made of durable vinyl with zinc-plated steel frame for longevity and rugged use
- Fully articulated for positioning in widely diversified training environments
- Three (3) year limited warranty

SKU	Description	Pkg	Price
MS150107	55 lbs; 5' 5"	1 ea	\$1,239.95
MS150108	105 lbs; 5' 5"	1 ea	\$1,361.95
MS150110	145 lbs; 5' 5"	1 ea	\$1,736.95
MS150112	165 lbs; 5' 5"	1 ea	\$1,864.95



MS150107



### Simulaids Combat Challenge Rescue Randy

- Combat Challenge Rescue Randy from Simulaids is the Official Manikin of the Firefighter Combat Challenge with articulated joints made from 4,200-lb. test plastic-coated cables
- Durable vinyl with a zinc-plated steel frame for rugged use in all types of environments
- Useful for simulating the removal of victims from utility poles, conduits, mine cave-ins, tunnels, and limitless other training procedures
- Three (3) year limited warranty

SKU	Description	Pkg	Price
MS150115	Dark Skin; 145 lbs	1 ea	\$2,152.95
MS150116	Dark Skin; 165 lbs	1 ea	\$2,010.95
MS150114	Light Skin; 145 lbs	1 ea	\$1,738.95
MS150188	Light Skin; 165 lbs	1 ea	\$1,867.95



MS150115



### Laerdal® Crash Kelly

- Crash Kelly is a durable, rugged training manikin with an intubation head for advanced airway management training and realistic articulation allowing the manikin to be placed in various settings for extrication or rescue
- Educationally effective for training in airway management and adult extrication
- Extremely durable, rugged and lifelike, made to withstand years of use

SKU	Description	Pkg	Price
MS139918	Laerdal® Crash Kelly	1 ea	\$2,999.00



MS139918



### Simulaids Pediatric Rescue Water Rescue Manikin

- Rescue Cathy newborn water filled weight: 7 lbs, Dimensions: 26" L x 8" W x 8" H
- Rescue Billy 6-9 month old water filled weight: 12 lbs, Dimensions: 26" L x 8" W x 8" H
- Rescue Timmy 3 year old water filled weight: 20 lbs, Dimensions: 34" L x 11" W x 7" H
- Three (3) year warranty


SKU	Description	Pkg	Price
MS140052	Newborn; Water Rescue Manikin	1 ea	\$290.95
MS140054	Infant; Water Rescue Manikin	1 ea	\$306.95
MS140053	Child; Water Rescue Manikin	1 ea	\$298.95





## Simulaids Water Rescue Manikins

- Durable vinyl plastic with rust-resistant stainless steel skeletons and articulated joints
- Manikin will submerge to the neck when filled with water
- Three (3) year warranty



MS140051

SKU	Description	Pkg	Price
MS140051	Adolescent	1 ea	\$1,021.95
MS140050	Adult	1 ea	\$1,130.95

## Airway Trainers



### Life/form® Airway Larry Trainer with Stand

- Simulates a non-anesthetized patient for practicing intubation, ventilation, suction and CPR techniques
- Offers tongue swelling and laryngospasm in addition to all the features on the standard models
- Features durable, rugged, one-piece construction and bifurcated lungs
- Five (5) year warranty



MS200016

SKU	Description	Pkg	Price
MS200016	Airway Trainer	1 ea	\$1,292.95



### Life/form® Airway Management Trainer with Stand

- Simulates a non-anesthetized patient for practicing intubation, ventilation, suction and CPR techniques
- A variety of anatomic features and airway skills trainer
- Oral and nasal intubation
- Includes: Airway management trainer mounted on sturdy base, 240 mL bottle of airway lubricant and a hard carry case
- Five (5) year warranty



MS024507

SKU	Description	Pkg	Price
MS024506	Infant	1 ea	\$944.95
MS024505	Child	1 ea	\$1,088.95
MS024507	Adult	1 ea	\$1,266.95



### Life/form® S.A.L.A.D. Simulator

- Suction Assisted Laryngoscopy and Airway Decontamination (S.A.L.A.D)
- Simulates a non-anesthetized patient for practicing intubation, ventilation and suction
- Anatomy and landmarks include teeth, tongue, oral and nasal pharynx, larynx, epiglottis, arytenoids, false cords and more
- Rigid suction catheter distracts lower mandible and tongue to permit laryngoscope insertion
- Offers room to maneuver and a slightly anterior position, making the simulator great for introductory classes as well as advanced students
- Five (5) year warranty



MS153991

SKU	Description	Pkg	Price
MS153992	Child	1 ea	\$1,858.95
MS153991	Adult	1 ea	\$1,918.95



### Simulaids 3 Year Old Airway Management Trainer

- Designed to teach oral and nasal endotracheal intubation procedures on children
- Detailed upper airway anatomy with highlighted vocal cords for laryngoscopy
- Includes sturdy nylon bag
- Three (3) year warranty

SKU	Description	Pkg	Price
MS140040	Pediatric Airway Trainer	1 ea	\$1,025.95



MS140040



### Simulaids Airway Management Trainer with Board

- Detailed respiratory anatomy with breakout teeth, tongue edema and laryngospasm
- Accepts up to size 8.0 oral and nasal endotracheal tubes
- Anatomical stomach swells with esophageal intubation or excessive BVM pressure during rescue breathing
- Vomiting capability requires students to clear the airway and prevent aspiration
- Includes hard carry case
- Three (3) year warranty

SKU	Description	Pkg	Price
MS029922	Adult Airway Trainer with Board	1 ea	\$2,341.95



MS029922



### Laerdal® Airway Trainer

- Simulates real-world complications when practicing a variety of intubation, ventilation and suction techniques
- Provides visual inspection of lung expansion and auscultation of breath sounds
- Simulates stomach inflation and vomiting situations
- One (1) year warranty

SKU	Description	Pkg	Price
MS024203	Adult Airway Trainer	1 ea	\$2,799.00



MS024203



### Simulaids Tension Pneumothorax Simulator

- Practice proper needle insertion to enable built-up air pressure to escape the thorax
- Foot pump, replacement two (2) piece pleural chamber, rib cage and lifelike skin
- Soft carry bag included
- Procedural needle not included

SKU	Description	Pkg	Price
MS046901	Tension Pneumothorax Simulator	1 ea	\$836.95



MS046901



### Life/form® Chest Tube Manikin

- Teach the theory, anatomy and skills needed to manage pre-hospital chest trauma and ongoing chest tube maintenance
- Teach the concepts and mechanics of closed water-seal drainage systems
- Three (3) year warranty

SKU	Description	Pkg	Price
MS044501	Chest Tube Simulator	1 ea	\$2,366.95



MS044501

**Life/form.**

**Life/form® Pericardiocentesis Simulator with Chest Tube and Pneumothorax**

- Pressurized tension pneumothorax site and a site for the surgical placement of a functional chest tube
- Fluid color, volume and viscosity can be controlled by the instructor
- Includes fluid reservoir bag, foot pump, surgical skin pads, subcutaneous surgical pads, nurse training pads, pneumothorax pads, blood powder, hydroxyethyl cellulose thickener, simulated pericardium, IV bag and a hard carry case
- Five (5) year warranty



EM090322

SKU	Description	Pkg	Price
EM090322	Pericardiocentesis Simulator	1 ea	\$2,268.95



**MST™ Needle Decompression Kit**

- One (1) protective vest
- One (1) rib plate with realistic muscle
- One (1) chest skin
- 0.125 lbs of glass shrapnel
- Twenty (20) needle catheters
- Four (4) simulated pleural space
- Four (4) air seals
- Two (2) battery operated air pumps
- Two (2) hand operated air pumps
- 4' of tubing
- Two (2) straight tubing connectors & clamps
- One (1) 4-Pack of D batteries
- One (1) 1 oz. thick blood
- One (1) 2 oz. simulated blood
- One (1) silicone adhesive
- One (1) mixing palette
- One (1) palette Knife
- One (1) scissors
- Fifty (50) sanitizing towelettes

SKU	Description	Pkg	Price
MG119138-LT	Light Skin Tone	1 kit	\$6,449.99
MG119138-MD	Medium Skin Tone	1 kit	\$6,499.99
MG119138-DK	Dark Skin Tone	1 kit	\$6,499.99



MG119138-LT

**Cricothyrotomy Trainers**



**MST™ Cricothyrotomy Kit**

- One (1) protective neck wrap with trachea bracket
- Thirty (30) simulated cricothyroid membranes
- Two (2) trachea structure
- Four (4) neck skins
- Thirty (30) neck skin inserts with blood pouch
- One (1) emergency cricothyrotomy kit
- One (1) silicone adhesive
- One (1) mixing palette
- One (1) palette knife
- One (1) scissors
- Fifty (50) sanitizing towelettes

SKU	Description	Pkg	Price
MG119137-LT	Light Skin Tone	1 kit	\$2,150.00
MG119137-MD	Medium Skin Tone	1 kit	\$2,150.00
MG119137-DK	Dark Skin Tone	1 kit	\$2,150.00



MG119137-LT





### Laerdal® Cricoid Stick Trainer

- Anatomically accurate landmarks for site training
- Interchangeable tracheas facilitate realistic simulation of needle and surgical cricothyrotomy procedures
- Replaceable neck skin with self-repairing memory allows repeated practice
- Mounted on base for stability and safety
- Soft trachea is suitable for surgical cricothyrotomy with incision
- Rigid trachea is suitable for 16G or 18G needle cricothyrotomy
- Includes one (1) reusable neck skin
- One (1) year warranty

SKU	Description	Pkg	Price
MS200040	Cricoid Stick Trainer	1 kit	\$729.00
MS150129	Replacement Skin	1 ea	\$35.00



MS200040



### Techline Trauma® Cric in a Box™ Task Trainer

- Includes simulated skin and membranes that allow the trainee to practice surgical incisions required to establish an emergency airway
- Simulated trachea has interior rings to provide tactile feedback when probed

SKU	Description	Pkg	Price
EM090323	Cricothyrotomy Task Trainer	1 kit	\$620.00
EM090324	Replacement Cric Membrane	25/pk	\$40.00
EM090325	Replacement Non-Bleed Neck Skins	10/pk	\$450.00



EM090323


## CPR Manikins



### Aero Brayden CPR Manikin

- Students see the effectiveness of CPR with visual LED blood flow to the brain and other organs
- LED lights indicate the quality of chest compression blood circulation CPR
- Features audible feedback, head tilt and chin lift and realistic chest movements with correct ventilations
- Compliant with new AHA feedback guidelines for CPR training

SKU	Description	Pkg	Price
MS150230	Brayden STATIC CPR Manikin	1 ea	\$270.00
MS150231	Brayden LED CPR Manikin with Blood Flow Circulation	1 ea	\$495.00
MS150232	Brayden PRO CPR Manikin with Bluetooth	1 ea	\$1,195.00



MS150231



### Laerdal® Little Baby QCPR

- Anatomically correct oral and nasal passages for nose-pinch, head tilt, chin lift and jaw thrust for chest rise
- Features an audible click when the correct compression depth is achieved
- Instructors can remotely monitor up to six (6) manikins at a time with the QCPR app
- Meets the 2019 AHA Feedback Requirements for CPR Manikins

SKU	Description	Pkg	Price
MS150218	Light Skin Tone; Infant CPR Manikin	1 ea	\$299.00
MS150219	Dark Skin Tone; Infant CPR Manikin	1 ea	\$299.00
MS150220	Light Skin Tone; Infant CPR Manikin	4/pk	\$1,149.00
MS150221	Dark Skin Tone; Infant CPR Manikin	4/pk	\$1,149.00



MS150219



### Laerdal® Little Anne QCPR

- Anatomically correct oral and nasal passages for nose-pinch, head tilt, chin lift and jaw thrust for chest rise
- Features an audible click when the correct compression depth is achieved
- Real-time CPR feedback allows instructors to review compression depth, release and rate, and ventilations
- Meets the 2019 AHA Feedback Requirements for CPR Manikins



MS139938

SKU	Description	Pkg	Price
MS139937	Light Skin Tone; Pediatric CPR Manikin	1 ea	\$319.00
MS139938	Dark Skin Tone; Pediatric CPR Manikin	1 ea	\$319.00
MS139939	Light Skin Tone; Pediatric CPR Manikin	4/pk	\$1,276.00
MS139940	Dark Skin Tone; Pediatric CPR Manikin	4/pk	\$1,276.00



### Laerdal® Resusci Anne QCPR

- Resusci Anne QCPR from Laerdal features realistic anatomical landmarks for adult chest compressions, head tilt/jaw thrust and chest rise with assisted ventilation
- Disposable lung system ensures students do not come into contact with exhaled air
- Laerdal QCPR, SimPad PLUS, and Skill Guide meet AHA requirements for CPR training and certification
- Includes 1 year manufacturer limited warranty



MS139974

SKU	Description	Pkg	Price
MS139974	Resusci Anne QCPR Torso with Carry Bag	1 ea	\$2,349.00
MS139975	Resusci Anne QCPR Full Body with Trolley Suitcase	1 ea	\$2,949.00
MS139976	Resusci Anne QCPR AW Torso with Carry Bag	1 ea	\$3,299.00
MS139977	Resusci Anne QCPR AW Full Body with Trolley Case	1 ea	\$3,899.00
MS139978	Resusci Anne QCPR AED Torso with Carry Bag	1 ea	\$2,949.00
MS139979	Resusci Anne QCPR AED Full Body with Trolley Suitcase	1 ea	\$3,649.00
MS139981	Resusci Anne QCPR AED AW Torso with Carry Bag	1 ea	\$3,899.00
MS139982	Resusci Anne QCPR AED AW Full Body with Trolley Suitcase	1 ea	\$4,599.00



### Prestan® Professional Infant Manikin with CPR Rate Monitor

- Revolutionary CPR Rate Monitor that allows for instant rate and depth feedback to both the instructor and student
- Allows for realistic infant CPR training
- Sculpted body contours and realistic skin feel help prepare students for real life rescues
- Face/head tilt enhances the CPR learning experience and teaches students how to open the airway during CPR



EM090395

SKU	Description	Pkg	Price
EM090395	Medium Skin Tone; Infant CPR Manikin	1 ea	\$156.00
EM090395-D	Dark Skin Tone; Infant CPR Manikin	1 ea	\$156.00
EM090395-DIVERSE	Diverse Tone; Infant CPR Manikin	1 ea	\$156.00
EM090396	Medium Skin Tone; Infant CPR Manikin	4/pk	\$599.00
EM090396-D	Dark Skin Tone; Infant CPR Manikin	4/pk	\$599.00
EM090396-DIVERSE	Diverse Tone; Infant CPR Manikin	4/pk	\$599.00



### Prestan® Professional Adult Manikin with CPR Rate Monitor

- Revolutionary CPR Rate Monitor that allows for instant rate and depth feedback to both the instructor and student
- Realistic to the eye and the touch, unlike any other CPR Training Manikin on the market
- Lightweight with fast and simple setup
- Fully compliant with current industry guidelines, including the most recent AHA Integrated Feedback Directive



EM090397

SKU	Description	Pkg	Price
EM090397	Medium Skin Tone; Infant CPR Manikin	1 ea	\$177.00
EM090397-D	Dark Skin Tone; Infant CPR Manikin	1 ea	\$177.00
EM090398	Medium Skin Tone; Infant CPR Manikin	4/pk	\$668.00
EM090398-D	Dark Skin Tone; Infant CPR Manikin	4/pk	\$668.00



### Simulaids Choking Manikins

- Functions on the principle of clearing a blocked airway with forced air generated with the proper procedure
- Visual confirmation of expelled object when appropriate force is applied
- Allow abdominal thrusts and back blows
- Includes three (3) anatomical landmarks including: Rib Cage, Xiphoid Process and Suprasternal Notch



MS026902

SKU	Description	Pkg	Price
MS140055	Infant	1 ea	\$371.95
MS026904	Child	1 ea	\$285.95
MS026903	Adolescent	1 ea	\$342.95
MS026902	Adult	1 ea	\$529.95
MS150128	Obese Adult	1 ea	\$586.95

## IO Trainers



### Simulaids Humerus Intraosseous (IO) Trainer

- Train the access procedure for fluid management
- Humerus access is now recognized as the best alternative when IV access cannot be obtained
- Presents three (3) orientations that can be demonstrated to provide the correct simulation for application
- Realistic trainer includes the correct hallmarks for accurate and safe access for fluid administration



MS154005

SKU	Description	Pkg	Price
MS154005	Humerus IO Trainer	1 ea	\$522.95



### Life/form® Adult Intraosseous Infusion Simulator

- Features palpable landmarks with soft skin covering and a pressurized system that allows both fluid aspiration and infusion
- Compatible with almost all IO infusion devices including B.I.G.™ and EZ-IO™
- Lower leg trainer comes with ten (10) replacement bones and four (4) replacement skin pads
- Five (5) year warranty



MS154002

SKU	Description	Pkg	Price
MS154002	Leg IO Simulator; Adult	1 ea	\$846.95

# SimVS EMS Series



SKU: DG073218  
SKU: DG073215

## Dynamic 12-Leads



Dynamic 12-leads available on both X-Series & LP-15 Defibrillator Interfaces

## Mobile Connectivity



Mobile router and battery pack allow for simulations in any environment

## Flexible Checklists



All students actions are recorded and available to be exported in a variety of formats

## Fit for Any Environment



Easily project instructor/student tablets for full classroom instruction



## SimVS Service Commitment

SUPPORT - 24/7 Support by Phone, Email, Webinar

TRAINING - Included With Packages

NO HIDDEN COSTS - No Annual Fees or Extra Screen Costs

(877) 593-6011 | [Info@DiaMedical.com](mailto:Info@DiaMedical.com) | [DiaMedical.com](http://DiaMedical.com)



**Life/form.**

**Life/form® Infant Intraosseous Infusion Simulator, Single Leg**

- Can be attached to Life/form® Infant CRiSis™ Manikin
- Fluid administration kit allows flashback and infusion of fluids into the leg
- Includes 15G IO needle, ten (10) replaceable simulated bones and two (2) sets of replacement skin
- Five (5) year warranty

SKU	Description	Pkg	Price
MS200037	Leg IO Simulator; Infant	1 ea	\$718.95



MS200037

**Life/form.**

**Life/form® Intraosseous Infusion Simulator**

- Designed to demonstrate and simulate the intraosseous infusion procedure
- Both legs have also been designed so that aspiration is possible
- Includes ten (10) bone inserts, two (2) sets of leg skins, two (2) towelettes, syringe with tubing, blood powder, sealing wax, 0.5 oz lubricant, mixing bottle, needle and hard carrying case
- Five (5) year warranty

SKU	Description	Pkg	Price
MS139916	Leg IO Simulator; Infant	1 ea	\$997.95



MS139916

**Life/form.**

**Life/form® Adult Sternal Intraosseous Infusion Simulator**

- Soft skin and accurate landmarks allow students to learn sternal IO access with skill and accuracy
- Replaceable simulated bones can be rotated and reused up to 10-12 times
- Sternum may be filled with simulated blood to better visualize flashback
- Includes ten (10) simulated bones, simulated blood packet and hard carry case
- Five (5) year warranty

SKU	Description	Pkg	Price
MS150138	Sternal IO Simulator; Adult	1 ea	\$1,961.95



MS150138

**SAM MEDICAL**

**SAM Medical® IO Training Kit**

- One (1) Training Case
- One (1) Driver
- One (1) Stabilizer
- One (1) Adapter
- One (1) 15 mm Needle Set
- One (1) 25 mm Needle Set
- One (1) 45 mm Needle Set
- One (1) Proximal Humerus Sawbone
- One (1) Proximal Tibia Sawbone
- One (1) Distal Tibia Sawbone

SKU	Description	Pkg	Price
EM090328	Kit with Proximal Humerus, Proximal Tibia & Distal Tibia	1 kit	\$785.99
EM090329	Kit with Proximal Humerus	1 kit	\$579.99
EM090330	Kit with Proximal Tibia	1 kit	\$579.99
EM090331	Kit with Distal Tibia	1 kit	\$579.99



EM090328

Wound Trainers



**TacMed Solutions™ Tourniquet Task Trainer Arm**

- Hyper-realistic arm looks exactly like a real human extremity
- Integrated gunshot entry wound can bleed
- Self-healing skin is designed for repeated use increased longevity
- Includes tourniquet trainer arm, hard carrying case, connector tubing and clamp and simulated IV bag



EM090333

SKU	Description	Pkg	Price
EM090333	Full-Length Arm; Tourniquet Trainer with Gunshot Wound	1 kit	\$490.50



**SIMBODIES® Shoulder Wound Packing Trainer**

- Contains one (1) large internal cavity and two (2) entry/exit gunshot through and through wound with internal bleeding tubes for simulated major hemorrhage
- Hemorrhage management using tourniquet, bandages and hemostatic agents
- Handmade using lifelike prosthetic grade platinum silicone that is durable and washable
- One (1) year warranty



MS150299-LT

SKU	Description	Pkg	Price
MS150299-LT	Light Skin Tone	1 ea	\$990.00
MS150299-MD	Medium Skin Tone	1 ea	\$990.00
MS150299-DK	Dark Skin Tone	1 ea	\$990.00

TUESDAY TEACHINGS  
SEE THIS PRODUCT  
IN ACTION >>>>



**SIMBODIES® Thigh Wound Packing Trainer**

- Contains one (1) superficial deep oval abrasion and one (1) entry/exit through and through gunshot wound with internal bleeding tubes for simulated major hemorrhage
- Hemorrhage management using tourniquet, bandages and hemostatic agents
- Handmade using lifelike prosthetic grade platinum silicone that is durable and washable
- One (1) year warranty



MS150300-LT

SKU	Description	Pkg	Price
MS150300-LT	Light Skin Tone	1 ea	\$1,210.00
MS150300-MD	Medium Skin Tone	1 ea	\$1,210.00
MS150300-DK	Dark Skin Tone	1 ea	\$1,210.00

TUESDAY TEACHINGS  
SEE THIS PRODUCT  
IN ACTION >>>>

**SIMBODIES**

**SIMBODIES® Full Leg Major Injury Task Trainer**

- Features an integrated major injury to the lower leg, exposing internal bone and muscle
- Accurate weight for added realism when applying a tourniquet to an immobile limb
- Enables students to practice hemorrhage control using tourniquets, bandages and haemostatic agents

SKU	Description	Pkg	Price
MS150295-LT	Light Skin Tone	1 ea	\$3,850.00
MS150295-MD	Medium Skin Tone	1 ea	\$3,850.00
MS150295-DK	Dark Skin Tone	1 ea	\$3,850.00



MS150295-LT

TUESDAY TEACHINGS  
SEE THIS PRODUCT  
IN ACTION ▶▶▶

**TECHLINE TRAUMA**  
TRAIN MORE. BLEED LESS!

**Techline Trauma® Gunshot in a Box™**

- Simulates a gunshot to soft-tissue, used to instruct students in proper wound packing technique
- Skin and tissue is made from high quality silicone rubber with layers of coloration for added realism
- Blood tubing is integrated into the construction to easily create simulated blood flow
- Bleeding systems sold separately

SKU	Description	Pkg	Price
EM090334	Light Skin Tone	1 ea	\$281.25
EM090335	Medium Skin Tone	1 ea	\$281.25
EM090336	Dark Skin Tone	1 ea	\$281.25



EM090334

**TECHLINE TRAUMA**  
TRAIN MORE. BLEED LESS!

**Techline Trauma® Laceration in a Box™**

- Simulates a soft-tissue wound, used to instruct students in proper wound packing technique
- Features a hard bone (fractured) inside of the wound cavity to elevate skill level
- Skin and tissue is made from high quality silicone rubber with layers of coloration for added realism
- Blood tubing is integrated into the construction to easily create simulated blood flow
- Bleeding systems sold separately

SKU	Description	Pkg	Price
EM090337	Light Skin Tone	1 ea	\$322.50
EM090338	Medium Skin Tone	1 ea	\$322.50
EM090339	Dark Skin Tone	1 ea	\$322.50



EM090337

**Pocket Nurse**  
Simulation & Education Supplies

**Pocket Nurse Suturing Skills Pad**

- Dimensions: 5" x 4.13" x 1.75"
- Realistic skin-like surface
- Illustrates epidermis, dermis and muscle layers
- Anti-skid base
- Excellent starter pad for incision and wound closure techniques

SKU	Description	Pkg	Price
DM70000	Suturing Skills Pad	1 ea	\$86.99



DM70000





## Take Emergency Medical Training to the Next Level in a Confined Ambulance Environment



### Why Educators Love SimRig® The Ambulance Trainer

SimRig® helps students understand the confines of an ambulance and enables them to provide life saving care within the spatial limitations of a rig. The Mobile SimRig® arrives at your facility road ready and can be trailered behind your vehicle to any training setting.

Contact Your Account Manager to Customize Your SimRig®

(877) 593-6011

| [DiaMedical.com](http://DiaMedical.com)

| [Info@DiaMedical.com](mailto:Info@DiaMedical.com)



Moulage



**Simulaids Advanced Military Casualty Simulation Kit**

Victims can appear to be in shock or pain, traumatized, lacerated, burned, bleeding or dirty. Bleeding wounds have reservoir pump assemblies to assure a continuous lifelike flow of simulated blood.

- Two (2) compound fractures of tibia, lower leg
- Two (2) compound fractures of humerus, upper arm
- Five (5) dozen packs of assorted adhesive wounds
- Shrapnel and debris
- Simulated body fluids
- Makeup and application supplies



MG126904

SKU	Description	Pkg	Price
MG126904	Military Casualty Moulage Kit	1 kit	\$937.95



**Simulaids Deluxe Casualty Simulation Kit**

This new kit contains everything needed for mass casualty and emergency medicine simulations. The full line of accessories allows creation of seventy (70) individual wounds.

- Fifteen (15) moulage stick-on wounds
- Seven (7) bleeding strap-on wounds with blood reservoir bags and pumps
- Makeup and application supplies
- Carry case



MG126905

SKU	Description	Pkg	Price
MG126905	Military Casualty Moulage Kit	1 kit	\$1,161.95



**Simulaids Xtreme Trauma Moulage Kit**

The Xtreme Trauma Moulage Kit represents wounds that result from explosions and chemical disasters. These simulated trauma wounds will assist in training exercises for military and government responders and anyone providing treatment to victims injured in disasters and war.

- Four (4) fractures
- Three (3) chemical burns
- One (1) avulsion
- One (1) burned face
- Three (3) lacerations
- Makeup and application supplies
- Shrapnel and debris



MG126975

SKU	Description	Pkg	Price
MG126975	Trauma Moulage Kit, Explosions & Chemical Disasters	1 kit	\$2,081.95



**Moulage Concepts Deluxe Nursing Moulage Kit**

- Features a combination of both basic and advanced moulage kits in the perfect way to get started in building medical moulage injuries in the simulation suite
- Includes the essentials needed to create basic, intermediate and advanced (medical, trauma, obstetrics and triage) moulage designs
- Contains moulage make-up supplies to create diabetic foot ulcers, surgical incisions, infiltrated IV's, rashes, cyanosis, blisters and pus
- Makeup and application supplies
- Carry case



MG126927

SKU	Description	Pkg	Price
MG126927	Military Casualty Moulage Kit	1 kit	\$899.99

## Simulaids Deluxe Xtreme Trauma Moulage Kit

The Deluxe Xtreme Trauma Moulage Kit includes all of the wounds from the Xtreme Trauma Moulage Kit and Xtreme Trauma Moulage Kit 2, and has sufficient makeup and accessories to create even more realism.

- One (1) Bleeding exit gunshot palm
- One (1) Burned face, 3rd degree
- One (1) Bleeding eviscerated intestine
- One (1) Bleeding forearm avulsion with skin flap
- One (1) Bleeding impaled stick
- One (1) Bleeding exit wound
- One (1) Bleeding laceration, 12 cm
- One (1) Electrical burn palm
- One (1) Intermediate shotgun entry
- One (1) 45-caliber close entry
- One (1) Bleeding open fracture, tibia
- One (1) Bleeding open fracture, femur
- One (1) Bleeding open fracture, clavicle
- One (1) Bleeding open wrist fracture
- One (1) Bleeding open rib fracture
- One (1) Half-face chemical burn
- One (1) Right hand 4th degree
- One (1) Chemical burn
- One (1) Dorsal hand 4th degree
- One (1) Chemical burn
- One (1) Bleeding exposed denture avulsion
- One (1) Bleeding open skull fracture
- One (1) Burned face
- One (1) Bleeding skin laceration
- One (1) Bleeding major muscle laceration
- One (1) Bleeding open ankle dislocation
- Makeup and application supplies



MG126902

SKU	Description	Pkg	Price
MG126902	Trauma Moulage Kit	1 kit	\$3,125.95



## Moulage Concepts EMS Trauma Kit

- Contains moulage make-up supplies to create bruises, cuts, bloody contusions, blisters, broken bones, pale pallor, near death demeanor and all manner of accident related theatrical training effects on live participants
- Ideal for mass casualty and emergency medicine simulation training

SKU	Description	Pkg	Price
MG126933	Military Casualty Moulage Kit	1 kit	\$1,208.99



MG126933



## Moulage Concepts Active Shooter Moulage & Training Kit

- This kit enables coordinated practice of rapid assessment, airway management and hemorrhagic response
- Utilizing interactive wounds, this kit increases visual and olfactory realism needed for an effective Emergency Operations Exercise

SKU	Description	Pkg	Price
MG126931	Active Shooter/Mass Casualty Moulage Kit	1 kit	\$1,854.00



MG126931



## Moulage Concepts Trauma/All Hazards Moulage Kit

- This moulage kit is appropriate for Active Shooter, Industrial Accident and Explosion simulations
- Contains necessary moulage make-up, supplies and accessories to create three-dimensional, interactive wounds that can be sutured, triaged, decontaminated and debrided
- Includes supplies for over two-hundred fifty (250) types of wounds

SKU	Description	Pkg	Price
MG126932	EMS/Nursing Trauma Moulage Kit	1 kit	\$2,015.99



MG126932



### MST™ Full Leg Amputation

- Simulates a complete traumatic leg amputation
- Made from realistic, durable and reusable silicone rubber
- Includes tubing and blood ports
- No adhesive necessary

SKU	Description	Pkg	Price
MG119101-LT-LT	Light Skin Tone; Left Leg	1 kit	\$1,299.99
MG119101-LT-MD	Medium Skin Tone; Left Leg	1 kit	\$1,299.99
MG119101-LT-DK	Dark Skin Tone; Left Leg	1 kit	\$1,299.99
MG119101-RT-LT	Light Skin Tone; Right Leg	1 kit	\$1,299.99
MG119101-RT-MD	Medium Skin Tone; Right Leg	1 kit	\$1,299.99
MG119101-RT-DK	Dark Skin Tone; Right Leg	1 kit	\$1,299.99



MG119101-LT-LT



### MST™ Partial Leg Amputation

- Simulates a partial traumatic leg amputation with foreleg still attached
- Made from realistic, durable and reusable silicone rubber
- Comes with tubing and blood ports
- Designed for the Meti CAE MetiMan
- No adhesive necessary
- Available in various skin tones

SKU	Description	Pkg	Price
MG119102-LT-LT	Light Skin Tone; Left Leg	1 kit	\$1,999.99
MG119102-LT-MD	Medium Skin Tone; Left Leg	1 kit	\$1,999.99
MG119102-LT-DK	Dark Skin Tone; Left Leg	1 kit	\$1,999.99
MG119102-RT-LT	Light Skin Tone; Right Leg	1 kit	\$1,999.99
MG119102-RT-MD	Medium Skin Tone; Right Leg	1 kit	\$1,999.99
MG119102-RT-DK	Dark Skin Tone; Right Leg	1 kit	\$1,999.99



MG119102-LT-LT



### MST™ Full Arm Amputation

- Simulates a complete traumatic arm amputation
- Made from realistic, durable and reusable silicone rubber
- Includes tubing and blood ports
- No adhesive necessary

SKU	Description	Pkg	Price
MG119103-LT-LT	Light Skin Tone; Left Arm	1 kit	\$999.99
MG119103-LT-MD	Medium Skin Tone; Left Arm	1 kit	\$999.99
MG119103-LT-DK	Dark Skin Tone; Left Arm	1 kit	\$999.99
MG119103-RT-LT	Light Skin Tone; Right Arm	1 kit	\$999.99
MG119103-RT-MD	Medium Skin Tone; Right Arm	1 kit	\$999.99
MG119103-RT-DK	Dark Skin Tone; Right Arm	1 kit	\$999.99



MG119103-LT-LT



### MST™ Partial Arm Amputation

- Simulates a partial traumatic arm amputation with forearm still attached
- Made from realistic, durable and reusable silicone rubber
- Includes tubing and blood ports
- No adhesive necessary

SKU	Description	Pkg	Price
MG119104-LT-LT	Light Skin Tone; Left Arm	1 kit	\$1,799.99
MG119104-LT-MD	Medium Skin Tone; Left Arm	1 kit	\$1,799.99
MG119104-LT-DK	Dark Skin Tone; Left Arm	1 kit	\$1,799.99
MG119104-RT-LT	Light Skin Tone; Right Arm	1 kit	\$1,799.99
MG119104-RT-MD	Medium Skin Tone; Right Arm	1 kit	\$1,799.99
MG119104-RT-DK	Dark Skin Tone; Right Arm	1 kit	\$1,799.99



MG119104-LT-LT



### MST™ Full Limb Amputation Kit

- Includes full left arm, full right arm, full left leg and full right leg
- Made from realistic, durable and reusable silicone rubber
- Comes complete with tubing and blood ports, allowing trainees to practice treating various types of amputations

SKU	Description	Pkg	Price
MG119105-LT	Light Skin Tone; 4-Piece	1 kit	\$5,999.99
MG119105-MD	Medium Skin Tone; 4-Piece	1 kit	\$5,999.99
MG119105-DK	Dark Skin Tone; 4-Piece	1 kit	\$5,999.99



MG119105-MD



### MST™ Multi-Amputation Kit

- Includes full left arm, full right arm, partial left arm, partial right arm, full left leg, full right leg, partial left leg and partial right leg
- Made from realistic, durable and reusable silicone rubber
- Includes tubing and blood ports
- No adhesive necessary

SKU	Description	Pkg	Price
MG119106-LT	Light Skin Tone; 8-Piece	1 kit	\$10,999.99
MG119106-MD	Medium Skin Tone; 8-Piece	1 kit	\$10,999.99
MG119106-DK	Dark Skin Tone; 8-Piece	1 kit	\$10,999.99



MG119106-LT



### MST™ Gunshot Wound

- Made from silicone rubber
- Highly realistic, durable, reusable and hypoallergenic
- Apply with silicone or manikin adhesive

SKU	Description	Pkg	Price
MG119119-1-LT	Light Skin Tone, Exit Wound	1 ea	\$10.99
MG119119-1-MD	Medium Skin Tone, Exit Wound	1 ea	\$10.99
MG119119-1-DK	Dark Skin Tone, Exit Wound	1 ea	\$10.99
MG119120-1-LT	Light Skin Tone, Entry Wound	1 ea	\$5.99
MG119120-1-MD	Medium Skin Tone, Entry Wound	1 ea	\$5.99
MG119120-1-DK	Dark Skin Tone, Entry Wound	1 ea	\$5.99



MG119119-1-LT





### MST™ Bleeding Wounds Skills Trainer Kit

- Silicone simulated glass
- Resin simulated shrapnel
- Resin simulated nuts, bolts and screws
- Twelve (12) simulated exit wounds
- Eighteen (18) simulated entrance wounds
- Ten (10) simulated basic lacerations
- One (1) simulated perforated bowel
- Ten (10) simulated dynamic bleeding wounds
- One (1) simulated dynamic compound fracture

SKU	Description	Pkg	Price
MG119135-LT	Light Skin Tone	1 kit	\$1,449.99
MG119135-MD	Medium Skin Tone	1 kit	\$1,449.99
MG119135-DK	Dark Skin Tone	1 kit	\$1,449.99



MG119135



### MST™ Simulated Urine

- Allows students to identify medical conditions based on urine color
- Will not evaporate during training

SKU	Description	Pkg	Price
MG119153-B	Odor; Blood	1 ea	\$59.99
MG119153-D	Odor; Dehydrated	1 ea	\$59.99
MG119153-N	Odor; Normal	1 ea	\$59.99
MG119153-C	Odor; Clear	1 ea	\$59.99
MG119152-B	Odorless; Blood	1 ea	\$25.99
MG119152-D	Odorless; Dehydrated	1 ea	\$25.99
MG119152-N	Odorless; Normal	1 ea	\$25.99
MG119152-C	Odorless; Clear	1 ea	\$25.99



MG119135



### MST™ Simulated Odors

- Designed to provide realistic olfactory simulation in addition to the visual simulation of MST prosthetics
- Scent oils are designed to simulate unpleasant scents of scenarios
- Do not use in eyes, in mouth or on human skin

SKU	Description	Pkg	Price
MG119150-BF	Burnt Flesh; 2 oz Bottle	1 ea	\$42.99
MG119150-CB	Car Bomb; 2 oz Bottle	1 ea	\$42.99
MG119150-DF	Decaying Flesh; 2 oz Bottle	1 ea	\$42.99
MG119150-FE	Feces; 2 oz Bottle	1 ea	\$42.99
MG119150-GP	Gunpowder; 2 oz Bottle	1 ea	\$42.99
MG119150-UR	Urine; 2 oz Bottle	1 ea	\$42.99
MG119150-VO	Vomit; 2 oz Bottle	1 ea	\$42.99



MG119150-BF



### Life/form® Venous Simulated Blood

- Life/form® simulation blood provides a sense of reality to any activity
- For use with Life/form® products



SKU	Description	Pkg	Price
EM090343	1 Gallon	1 ea	\$47.32



### Life/form® Arterial Blood

- Life/form® simulation blood provides a sense of reality to any activity
- For use with Life/form® products



SKU	Description	Pkg	Price
EM090341	1 Quart	1 ea	\$34.95

IV Skills Trainers



**Demo Dose® Inject-Ed™ Injection Trainer**

- Intramuscular (IM), subcutaneous (SQ), intradermal (ID) and z-track simulation
- Diameter: 3"
- Resilient and durable
- Lifelike skin surface
- Safety inject air or liquids


SKU	Description	Pkg	Price
DM50808-LGT	Light Skin Tone	1 ea	\$24.99
DM50808-MED	Medium Skin Tone	1 ea	\$24.99
DM50808-DRK	Dark Skin	1 ea	\$24.99



**BT Inc. Adult IV Injection Training Arm**

- Features a blood circulation pump for realistic flashback when venipuncture has been achieved
- Skin and veins have realistic texture and resistance to venipuncture
- Practice on basilic, cephalic, median cubital and metacarpal veins
- Skin pad includes epidermis, vein and subcutaneous fat

SKU	Description	Pkg	Price
MS010607	Large; 3.5"	1 ea	\$2,500.00



MS010607



**Life/Form® Advanced Injection IV Simulator - Arm**

- Provides complete venous access for IV therapy and phlebotomy, plus sites for intramuscular and intradermal injections
- Eight (8) line vascular system
- Simplified setup with one external fluid bag supplying artificial blood to all veins

SKU	Description	Pkg	Price
MS074503	Light Skin Tone	1 ea	\$949.95
MS074512	Medium Skin Tone	1 ea	\$949.95
MS074593	Dark Skin Tone	1 ea	\$949.95



MS074593

Simulated IV Fluids & Medications

TUESDAY TEACHINGS  
SEE THESE PRODUCTS  
IN ACTION ▶▶▶



**Demo Dose® Prefilled Syringe with Water**

- Available as 3 or 10 mL
- Packaging: Pre-filled, pre-assembled plastic syringe

SKU	Pkg	Price
PM300272	3 mL	\$2.99
PM300271	10 mL	\$3.09



**Demo Dose® Heparin Lock Flush Syringe**

- Dosage: 1000 units/mL
- Packaging: Pre-filled, pre-assembled plastic syringe

SKU	Pkg	Price
PM208050	1 ea	\$2.39



**Demo Dose® NaCl 0.9% Normal Saline Prefilled Syringe**

- Dosage: 0.9% Sodium Chloride
- Packaging: Pre-filled, pre-assembled plastic syringe

SKU	Pkg	Price
PM300207	3 mL	\$2.99
PM300190	10 mL	\$3.09



DEMO DOSE  
Simulated Medication for Education



**Demo Dose® Adenocrd 6 mg/2 mL**

SKU	Pkg	Price
PM207991	1 ea	\$4.49



DEMO DOSE  
Simulated Medication for Education



**Demo Dose® Adenocrd 12 mg/4 mL**

SKU	Pkg	Price
PM207992	1 ea	\$4.49



DEMO DOSE  
Simulated Medication for Education



**Demo Dose® Atropin 1 mg/10 mL**

SKU	Pkg	Price
PM207998	1 ea	\$5.19



DEMO DOSE  
Simulated Medication for Education



**Demo Dose® Calcium Chlorid 10% 1 gram/10 mL**

SKU	Pkg	Price
PM207999	1 ea	\$5.19



DEMO DOSE  
Simulated Medication for Education



**Demo Dose® Dextros 50% 25 mg/50 mL**

SKU	Pkg	Price
PM208004	1 ea	\$8.29



DEMO DOSE  
Simulated Medication for Education



**Demo Dose® EPINEPHrin 1 mg/10 mL**

SKU	Pkg	Price
PM207996	1 ea	\$5.19




DEMO DOSE  
Simulated Medication for Education



**Demo Dose® Lidocain 2% 100 mg/5 mL**

SKU	Pkg	Price
PM207997	1 ea	\$5.19



DEMO DOSE  
Simulated Medication for Education



**Demo Dose® Naloxon 2 mg/2mL**

SKU	Pkg	Price
PM300172	1 ea	\$4.09



DEMO DOSE  
Simulated Medication for Education



**Demo Dose® Sodim Bicarbonat 8.4% 50 mEq/50 mL**

SKU	Pkg	Price
PM208003	1 ea	\$8.29

**Demo Dose® Code Med Set #1**

- One (1) Epinephrine
- One (1) Atropine
- One (1) Calcium Chloride



- One (1) Lidocaine
- One (1) Dextrose
- One (1) Sodium Bicarb

SKU	Description	Pkg	Price
PM300415	6 Pack	6/pkg	\$38.99

DEMO DOSE  
Simulated Medication for Education

**Demo Dose® Code Med Set #2**

- One (1) Epinephrine
- One (1) Atropine
- One (1) Calcium Chloride
- One (1) Lidocaine
- One (1) Dextrose



- One (1) Sodium Bicarb
- One (1) Vasopressin 10mg
- One (1) Amiodaron 150mg
- One (1) Furosemid 10mg

SKU	Description	Pkg	Price
PM300417	9 Pack	9/pkg	\$44.99



TUESDAY TEACHINGS  
SEE THESE PRODUCTS  
IN ACTION



# SIMULATED BLOOD, PLASMA & PLATELETS

SimRx Simulated Blood, Plasma and Platelet Bags can be used to teach blood transfusion procedures. Bags are available in all eight (8) blood types, have barcoded labels and are filled with fluid that is realistic in color and texture.

Blood Types	A+	A-	B+	B-	O+	O-	AB+	AB-	Variety Pack (All 8 Blood Types)	Price
Stainless Blood Bags	IV158640	IV158641	IV158648	IV158642	IV158645	IV158646	IV158643	IV158644	IV158647	\$16.99   Filled \$5.99   Empty \$102.99   Variety Pack
Simulated Blood Bags	IV058640	IV058641	IV058642	IV058643	IV058646	IV058647	IV058644	IV058645	IV058648	\$12.99   Filled \$5.99   Empty \$77.99   Variety Pack
Simulated Plasma Bags	IV058650	IV058651	IV058652	IV058653	IV058656	IV058657	IV058654	IV058655	IV058658	\$12.99   Filled \$5.99   Empty \$77.99   Variety Pack
Simulated Platelet Bags	IV058660	IV058661	IV058662	IV058663	IV058666	IV058667	IV058664	IV058665	IV058668	\$12.99   Filled \$5.99   Empty \$77.99   Variety Pack



TUESDAY TEACHINGS  
SEE THESE PRODUCTS  
IN ACTION



# SIMULATED IV FLUIDS

Now you can teach IV administration with this cost-effective solution from SimRx. SimRx Simulated IV Fluids were specifically designed for healthcare educators to provide a safe and realistic learning environment. SimRx Simulated IV Fluids contain distilled water so they are safe to use with your task trainers, high fidelity manikins and infusion pumps.

IV Fluid	0.9% Sod Chlor-ide	0.9% Sod Chlor-ide Blue Capped Port	5% Dex-trose	5% Dex-trose & 0.9% Sod Chlor-ide	Lac-tated Ringr's	Individual	Quantity & Price
50 mL	IV058601	IV058601-B	IV058610	IV058620	-	\$2.99 Filled	100   \$263.99
100 mL	IV058602	IV058602-B	IV058611	IV058621	-	\$2.99 Filled	100   \$274.99
250 mL	IV058603	IV058603-B	IV058612	IV058622	-	\$2.99 Filled	50   \$134.99
500 mL	IV058604	IV058604-B	IV058613	IV058623	IV058630	\$2.99 Filled	50   \$145.99
1,000 mL	IV058605	IV058605-B	IV058614	IV058624	IV058631	\$2.99 Filled	25   \$77.99





**EMERGENCY**  
BAGS







## Field Ready Loaded BLS Jump Bag™

Save valuable time with the Field Ready Loaded BLS Jump Bag™. This jump bag comes loaded with first aid supplies and diagnostic equipment required to perform basic life support.



SKU	Description	Pkg	Price
LC129912-FR	DiaMedical Field Ready Loaded BLS Jump Bag™	1 ea	\$325.99
LC129922-FR	DiaMedical Field Ready Loaded BLS Jump Bag™ Refill Kit	1 ea	\$269.99
LC129912	SimLabSolutions Loaded BLS Jump Bag™; For Training Use Only	1 ea	\$325.99



### PPE/Infection Control

- Two (2) Pairs Sterile Exam Gloves, Vinyl
- Nine (9) Alcohol Gel Hand Wipes



### Diagnostic

- One (1) Glucose Meter
- Fifty (50) Test Strips
- Fifty (50) Safety Lancets
- One (1) Aneroid BP Cuff, Adult
- One (1) Stethoscope
- One (1) Penlight with Pupil Gauge
- One (1) Digital Thermometer
- One (1) Pulse Oximeter, Fingertip



### Respiratory

- One (1) BVM, Adult
- One (1) Disposable CPR Face Shield
- One (1) Each Disposable Guedel Airway - 50 mm, 60 mm, 70 mm, 80 mm, 90 mm & 100 mm
- One (1) Each Rusch Nasal Airway - 16 FR, 26 FR & 28 FR
- Six (6) Sterile Lubricating Jelly Packets
- One (1) Bite Stick



### Wound Care/Dressings

- One Hundred (100) Plastic Adhesive Strips - 0.75" x 3"
- One (1) Elastic Bandage - 3" x 5 yd
- Two (2) Sterile Trauma Dressings - 12" x 30"
- Two (2) Sterile Burn Dressings
- Two (2) Triangular Bandages
- Five (5) Sterile ABD Pads - 5" x 9"
- Two (2) Sterile Transparent Dressings - 4" x 4"
- Fifteen (15) Sterile Tongue Depressors
- Two (2) Sterile Sof-Form Gauze Rolls - 1"
- Three (3) Sterile Sof-Form Gauze Rolls - 3"
- Three (3) Sterile Gauze - 4" x 4"
- Four (4) Sterile Eye Pads
- Two (2) Transparent Adhesive Tape Rolls - 1"



### Additional Supplies

- One (1) Splinter Forceps
- One (1) Utility Scissors - 7.5"
- One (1) Kelly Forceps, Straight - 5.5"
- One (1) Instant Cold Pack - 5" x 9"
- One (1) Obstetrical Kit
- One (1) Emergency Survival Blanket
- Twenty-Five (25) Alcohol Prep Pads
- One (1) Each Cervical Collar - Adult & Pediatric



## Field Ready Loaded ALS Jump Bag™

Stay life support ready with the Field Ready Loaded ALS Jump Bag™! This loaded jump bags arrives stocked with airway supplies, diagnostic equipment, IV supplies and first aid supplies required for advanced life support.



SKU	Description	Pkg	Price
LC129911-FR	DiaMedical Field Ready Loaded ALS Jump Bag™	1 ea	\$745.99
LC129921-FR	DiaMedical Field Ready Loaded ALS Jump Bag™ Refill Kit	1 ea	\$595.99
LC129911	SimLabSolutions Loaded ALS Jump Bag™; For Training Use Only	1 ea	\$745.99



### PPE/Infection Control

- Two (2) Pairs Sterile Exam Gloves, Vinyl
- Nine (9) Alcohol Gel Hand Wipes
- One (1) Procedure Mask with Face Shield



### Diagnostic

- One (1) Glucose Meter
- Fifty (50) Test Strips
- Fifty (50) Safety Lancets
- One (1) Aneroid BP Cuff, Adult
- One (1) Stethoscope
- One (1) Penlight with Pupil Gauge
- One (1) Digital Thermometer
- One (1) Pulse Oximeter, Fingertip



### Infusion/Injection

- Twenty-Five (25) Alcohol Prep Pads
- Three (3) Latex-Free Tourniquets, Rolled - 1" x 18"
- One (1) IV Armboard, Adult - 3" x 9"
- Two (2) Each Luer Lock Syringes - 3mL, 5mL & 10 mL
- Two (2) Tuberculin Safety Syringes
- Four (4) Each Safety Hypodermic Needles - 25G & 18G
- Four (4) IV Start Kits with Chloraprep
- One (1) IV Extension Set - 7"
- One (1) 3-Way Stopcock
- Two (2) Primary IV Sets - Microdrip & Macro drip
- Twenty (20) Alcohol Swabsticks
- One (1) Safety IV Catheter - 24G, 22G, 18G & 14G



### Additional Supplies

- Two (2) Transparent Adhesive Tape Rolls - 1"
- One (1) Splinter Forceps
- One (1) Utility Scissors - 7.5"
- One (1) Kelly Forceps, Straight - 5.5"
- One (1) Instant Cold Pack - 5" x 9"
- One (1) Obstetrical Kit
- One (1) Emergency Survival Blanket
- One (1) Each Cervical Collars - Adult & Pediatric



### Respiratory

- One (1) BVM, Adult
- One (1) Intubation Bag
- One (1) Emergency Cricothyrotomy Kit
- One (1) Each Nasal Cannulas - Adult & Pediatric
- One (1) Each Non-Rebreather Oxygen Masks - Adult & Pediatric
- Six (6) Sterile Lubricating Jelly Packets
- One (1) Bite Stick
- Five (5) Sterile Tongue Depressors
- One (1) Each Disposable Guedel Airway - 50 mm, 60 mm, 70 mm, 80 mm, 90 mm & 100 mm
- One (1) Each Rusch Nasal Airway - 16 FR, 26 FR, 28 FR & 30 FR
- One (1) Disposable Laryngoscope Handle
- One (1) Each Disposable Laryngoscope Blades, Miller - #0, #1 & #2
- One (1) Each Disposable Laryngoscope Blades, Macintosh - #2, #3 & #4
- One (1) Each Cuffed Endotracheal Tubes - 3.0 mm, 4.0 mm, 5.0 mm, 6.0 mm, 7.0 mm, 8.0 mm & 9.0 mm
- One (1) Each Intubation Stylets - 10 FR & 14 FR
- One (1) CO2 Indicator, Colorimetric
- One (1) ET Tube Holder
- One (1) Magill Catheter Forceps, Adult - 9.75"



### Wound Care/Dressings

- One Hundred (100) Plastic Adhesive Strips - 3/4" x 3"
- Two (2) Sterile Trauma Dressings - 12" x 30"
- Two (2) Sterile Sof-Form Gauze Rolls - 1"
- Three (3) Sterile Sof-Form Gauze Rolls - 3"
- Three (3) Sterile Gauze - 4" x 4"
- Five (5) Sterile Gauze - 2" x 2"
- Two (2) Sterile Burn Dressings
- Five (5) Sterile ABD Pads - 5" x 9"
- Four (4) Sterile Eye Pads
- Two (2) Sterile Transparent Dressings - 4" x 4"
- One (1) Elastic Bandage - 3" x 5 yd
- Two (2) Triangular Bandages



## Field Ready Loaded Pediatric Jump Bag™

Pediatric emergencies are stressful enough, let DiaMedical make sure your jump bag is ready! The Field Ready Loaded Pediatric Jump Bag™ arrives with pediatric emergency equipment and supplies organized in Broselow Color-Coded modules for simple identification.



SKU	Description	Pkg	Price
LC129913-FR	DiaMedical Field Ready Loaded Pediatric Jump Bag™	1 ea	\$755.99
LC129923-FR	DiaMedical Field Ready Loaded Pediatric Jump Bag™ Refill Kit	1 ea	\$605.99
LC129913	SimLabSolutions Loaded Pediatric Jump Bag™; For Training Use Only	1 ea	\$755.99



### Color-Coded Weight Based Modules

- One (1) Suction Catheter Kit
- One (1) Cuffed Murphy Endotracheal Tube
- One (1) Disposable Laryngoscope Blade, Miller and/or Macintosh
- One (1) Intubating Stylet
- One (1) Disposable Guedel Airway
- One (1) Rusch Nasal Airway
- One (1) Safety IV Catheter
- One (1) Butterfly Infusion Set
- Two (2) IV Start Kits
- Ten (10) Alcohol Prep Pads
- Three (3) Sterile Lubricating Jelly Packets



### Additional Supplies

- One (1) Each Disposable IV Armboards - Neonatal, Infant & Adult
- Two (2) Each Luer Lock Syringes - 10 mL & 3 mL
- Four (4) Pediatric/Neonatal ECG Electrodes
- One (1) Each BVMs - Child & Infant
- One (1) Each Aneroid BP Cuffs - Pediatric & Infant
- One (1) Pediatric Stethoscope
- Two (2) Sterile Bulb Syringes - 60 mL
- One (1) Obstetrical Kit
- One (1) PediaQuick Tape
- One (1) Each Primary IV Sets - Microdrip & Macro drip
- One (1) IV Extension Set - 7"
- One (1) Disposable Laryngoscope Handle
- One (1) Whole Blood Tube Plastic Vacutainer
- One (1) 3-Way Stopcock
- One (1) Magill Catheter Forceps, Child
- One (1) Each Nasal Cannulas - Pediatric & Infant
- One (1) Non-Rebreather Oxygen Mask - Pediatric
- One (1) CO2 Indicator, Colorimetric

## TUESDAY TEACHINGS

### Bring Our Products Into Your Practice

Learn why our EMS Loaded Jump Bags™ are an essential tool for life-saving procedures & training scenarios!

[WATCH NOW](#) ▶▶▶▶▶





## Equipment Bags

**Kemp USA Ultimate EMS Backpack**

- Features two (2) large frontal pockets, two (2) side pockets, a large main compartment with four (4) utility pouches and two (2) hook and loop straps for "D" tank cylinder
- Material: 600D Tear-Resistant Nylon
- Dimensions: 18" L x 13" W x 20" D



SKU	Description	Pkg	Price
EM090348	Red	1 ea	\$264.00

**Kemp USA Premium Ultimate EMS Backpack**

- Features two (2) large frontal pockets, two (2) side pockets, a large main compartment with four (4) utility pouches and two (2) hook and loop straps for "D" tank cylinder
- Material: 600D Tear-Resistant Nylon
- Dimensions: 18" L x 13" W x 20" D



SKU	Description	Pkg	Price
EM090349-RED	Red	1 ea	\$297.00
EM090349-NVY	Navy	1 ea	\$297.00
EM090349-RTP	Red Tarpaulin; Fluid Proof	1 ea	\$306.90

**Kemp USA Large Professional Trauma Bag**

- Includes a front pocket with loops and pouches, an easy access top pocket and internal compartments
- Features replaceable bottoms, tuck-away backpack harness and padded shoulder strap
- Material: 600D Nylon
- Dimensions: 17.5" L x 13" W x 10" D



SKU	Description	Pkg	Price
EM090350	Red	1 ea	\$118.13
EM090351	Royal Blue	1 ea	\$118.13

**Kemp USA Maxi Trauma Bag**

- Features an adjustable should strap, elastic loops on top flap for oral airways/ small equipment, customizable padded center compartment and two (2) large zippered end pockets
- Material: 600D Nylon
- Dimensions: 18" L x 13" W x 10" D



SKU	Description	Pkg	Price
EM090353	Orange	1 ea	\$49.50
EM090354	Red	1 ea	\$49.50
EM090355	Navy	1 ea	\$49.50

**Kemp USA Premium Large Professional Trauma Bag**

- Includes a front pocket with elastic loops and pouches, an easy-access top pocket with two (2) external side pockets and a large internal compartment with removable dividers and zippered see-through pouch
- Features tuck-away backpack harness, gel-foam shoulder strap and pull-put kneeling pad
- Material: 1680D Denier Nylon
- Dimensions: 17.5" L x 13" W x 10" D



SKU	Description	Pkg	Price
EM090352-RED	Red	1 ea	\$155.93
EM090352-NVY	Navy	1 ea	\$148.05
EM090352-RTP	Red Tarpaulin; Fluid Proof	1 ea	\$118.13

**Kemp Ultra EMS Bag**

- Features eight (8) compartments, two (2) pockets on each side and two (2) removable zippered pockets
- Includes bottom feet, padded nylon duffel handles and an adjustable shoulder strap
- Material: 600D Nylon
- Dimensions: 23" L x 12" W x 12" D



SKU	Description	Pkg	Price
EM039310	Red	1 ea	\$130.50



### StatPacks G3 Load N' Go Backpack

- Features transparent sleeves, mesh pockets and elastic holsters for organization and foam molded shoulder straps
- Includes side pockets and quick stash pockets
- Material: 4200 Denier Nylon
- Dimensions: 21" L x 16" W x 7" D



SKU	Description	Pkg	Price
EM090128	Red	1 ea	\$345.45
EM090129	Black	1 ea	\$300.00
EM090423	Green	1 ea	\$300.00



### StatPacks G3 Clinician Jump Bag

- Features three (3) zippered exterior pockets and one (1) large main zippered compartment
- Unique foam lined compartments and foam molded shoulder straps
- Material: Tarpaulin
- Dimensions: 22" L x 16" W x 9.5" D



SKU	Description	Pkg	Price
EM090132	Red	1 ea	\$315.00



### Iron Duck® BLS Backpack

- Top front pocket features multi-sized elastic loops
- Main compartment includes four (4) color-coded removable pouches with clear windows
- Oxygen compartment equipped with quick release tie-downs to secure a "D" or Jumbo "D" cylinder
- Dimensions: 21" L x 14" W x 8" D



SKU	Description	Pkg	Price
EM090358-RED	Red	1 ea	\$259.95
EM090358-BLK	Black	1 ea	\$259.95
EM090358-GRE	Green	1 ea	\$259.95
EM090358-RYB	Royal Blue	1 ea	\$259.95



### StatPacks G3 Medslinger

- Unique single shoulder strap design specifically built for medics
- Features clear, durable urethane windows on inside pouches and exterior quick stash pocket
- Material: Tarpaulin
- Dimensions: 18.5" L x 13" W x 3" D



SKU	Description	Pkg	Price
EM090158	Red	1 ea	\$142.80
EM090159	Tactical Black	1 ea	\$142.80
EM090424	Green	1 ea	\$142.80
EM090157	Blue	1 ea	\$142.80



### Iron Duck® Backpack Plus

- Includes two (2) exterior pockets and one (1) main compartment with three (3) removable color-coded pouches with clear windows
- O2 compartment equipped with quick-release tie-downs to secure a "D" or Jumbo "D" cylinder
- Dimensions: 23" L x 14" W x 12" D
- Additional colors available online



SKU	Description	Pkg	Price
EM090357-BLK	Black	1 ea	\$408.45
EM090357-RED	Red	1 ea	\$408.45
EM090357-RYB	Royal Blue	1 ea	\$408.45
EM090357-BUP	Black; UP Material	1 ea	\$449.30
EM090357-RUP	Red; UP Material	1 ea	\$449.30
EM090357-RBUP	Royal Blue; UP Material	1 ea	\$449.30



### Iron Duck® Pack Case Triple

- Features main compartment with removable padded dividers, three (3) exterior pockets and top lip has flat zippered UP pocket with elastic tool holders
- Material: 1000D Nylon
- Dimensions: 21" L x 15" W x 10" H



SKU	Description	Pkg	Price
EM090356-BLK	Black	1 ea	\$164.95
EM090356-ORG	Orange	1 ea	\$164.95
EM090356-RED	Red	1 ea	\$164.95
EM090356-RBL	Royal Blue	1 ea	\$164.95
EM090356-OUP	Orange; UP Material	1 ea	\$182.95
EM090356-RUP	Red; UP Material	1 ea	\$182.95
EM090356-BUP	Royal Blue; UP Material	1 ea	\$182.95

**MERET OMNI™ PRO X**

- Internal main compartment can be customized by utilizing the slider dividers system
- Hideaway, mountaineering style backpack harness
- OneView™ technology allows the bag to open completely facing one direction
- Dimensions: 15" L x 22" W x 9.5" D
- Additional colors available online



SKU	Description	Pkg	Price
EM039316-BLK	Black; Standard	1 ea	\$399.95
EM039316-BBL	Black; Ballistics	1 ea	\$599.95

**MERET DEFENDER™ PRO**

- Features slim rear zippered pocket, large middle zippered pocket, full-length outer zippered pocket, small front zippered pocket, multifunction bottom padded zippered pocket and two (2) external mesh pockets
- Dimensions: 18" x 8" x 11.5"



SKU	Description	Pkg	Price
EM039315	Black; Standard	1 ea	\$119.95
EM039315-AR	Black; Ballistics	1 ea	\$319.95

**Oxygen Bags****Iron Duck® Oxygen Bag; Size "E"**

- Designed to carry a "E" oxygen tank, regulator and accessories
- Dimensions: 32" L x 8.5" W x 8.5" D



SKU	Description	Pkg	Price
EM090429	Black	1 ea	\$177.63
EM090430	Green	1 ea	\$177.63
EM090361	Red; UP Material	1 ea	\$195.39

**Iron Duck® Oxygen Bag; Size "D"**

- Designed to carry a "D" oxygen tank, regulator and accessories
- Dimensions: 21" L x 8.5" W x 8.5" D



SKU	Description	Pkg	Price
EM090359	Black	1 ea	\$147.00
EM090425	Green	1 ea	\$147.00
EM090426	Red; UP Material	1 ea	\$161.70

**Iron Duck® Oxygen Bag with Pockets**

- Designed to carry a "D" or class Jumbo "D" oxygen tank, regulator and accessories
- Features two (2) full length side pockets
- Dimensions: 21" L x 11.5" W x 8.5" D



SKU	Description	Pkg	Price
EM090360	Black	1 ea	\$182.00
EM090427	Green	1 ea	\$182.00
EM090428	Red; UP Material	1 ea	\$200.19

**StatPacks G3 Breather**

- Designed to carry a "D" or Jumbo "D" oxygen tank
- Includes main compartment, side pockets and stash pocket with clear urethane windows, mesh pockets and elastic holsters
- Dimensions: 20" L x 19" W x 7" D




SKU	Description	Pkg	Price
EM090376-RED	Red	1 ea	\$383.25
EM090376-GRN	Green	1 ea	\$383.25
EM090376-RYB	Royal Blue	1 ea	\$383.25



## Iron Duck® Breathsaver Plus

- Designed to carry a “D” or Jumbo “D” oxygen tank
- Fully-opening end pocket offers quick access to oxygen tank regulator and shut off valve
- Large main compartment offers a removable liner over closed-cell foam and plastic stiffener
- Features three (3) large exterior pockets and interior walls are lined with mesh pouches and elastic holders
- Dimensions: 27” L x 13” W x 11” D



EM090362-GRE

SKU	Description	Pkg	Price
EM090362-BLK	Black	1 ea	\$299.25
EM090362-GRE	Green	1 ea	\$299.25
EM090362-RED	Red	1 ea	\$299.25
EM090362-RYB	Royal Blue	1 ea	\$299.25
EM090362-BUP	Black; UP Material	1 ea	\$329.18
EM090362-GUP	Green; UP Material	1 ea	\$329.18
EM090362-RUP	Red; UP Material	1 ea	\$329.18
EM090362-RBUP	Royal Blue; UP Material	1 ea	\$329.18

## Hard Cases



### Pelican™ 1500EMS Protector EMS Case

- Watertight, crushproof and dustproof
- Multi-layer lid organizer with pockets and Velcro® brand fasteners
- Includes foam inserts
- Designed with an automatic purge valve to equalize air pressure, a watertight silicone O-ring lid, over-molded rubber handles and stainless steel hardware
- Dimensions: 18.5” L x 14” W x 7” D



SKU	Description	Pkg	Price
EM039336-ORG	Orange	1 ea	\$224.95
EM039336-BLK	Black	1 ea	\$224.95



### Pelican™ 1600EMS Protector EMS Case

- Watertight, crushproof and dustproof
- Multi-layer lid organizer with pockets and Velcro® brand fasteners
- Includes dividers with adjustable walls for custom sized compartments
- Designed with an automatic purge valve to equalize air pressure, a watertight silicone O-ring lid, over-molded rubber handles and stainless steel hardware
- Dimensions: 25” L x 20” W x 9” D



SKU	Description	Pkg	Price
EM039338-ORG	Orange	1 ea	\$334.95
EM039338-BLK	Black	1 ea	\$334.95



### Pelican™ 1460EMS Protector EMS Case

- Top lid opens up to 180° for on-the-field anti-roll stability when open
- Includes sharps container storage bracket, adjustable tray dividers and lockable compartment for controlled substances
- Designed with an automatic purge valve to equalize air pressure, a watertight silicone O-ring lid, over-molded rubber handles and stainless steel hardware
- Dimensions: 21” L x 13” W x 13” D



EM039314-ORG

SKU	Description	Pkg	Price
EM039314-ORG	Orange	1 ea	\$399.99




## Supply Modules

**StatPacks G3 Airway Cell**

- Unique semi-rigid, quad-fold design helps protect airway supplies
- Color-coded for easy identification from other EMT airway bags
- Two (2) handles with "AIRWAY" written to identify the difference in EMT Gear
- Dimensions: 12" L x 7" W x 6" D



SKU	Description	Pkg	Price
EM090377	 Green	1 ea	\$154.88

**StatPacks G3 Medicine Cell**

- Unique semi-rigid, quad-fold design helps protect and organize medication
- Color-coded for easy identification from other EMS medication bags
- Two (2) handles with "MEDICINE" written to identify the difference in EMT Gear
- Dimensions: 12" L x 7" W x 6" D




SKU	Description	Pkg	Price
EM090378	 Red	1 ea	\$166.95

**StatPacks G3 Intravenous Cell**

- Color-coded for easy identification from other EMS I.V. Bags
- Two (2) handles with "INTRAVENOUS" written to identify the difference in EMT Gear
- Dimensions: 12" L x 7" W x 5" D



SKU	Description	Pkg	Price
EM090379	 Blue	1 ea	\$131.25

**MERET IV MED™ PRO X Medication/IV Start Kit**

- Designed to carry essential supplies for advanced life support medics
- Constructed from non-absorbing, wipe-clean infection control material
- OneView™ technology completely opens all main compartments
- Dimensions: 8.5" L x 16" W x 1.5" D







SKU	Description	Pkg	Price
EM090368	 Black	1 ea	\$64.95

**MERET AIRWAY™ PRO X Intubation Module**

- Roll-out, tri-fold easy access system organizes airway supplies
- Features clear internal pockets for quick identification
- Constructed from wipe-clean infection control materials internally and externally
- Dimensions: 6.25" L x 10.25" W x 3.95" D







SKU	Description	Pkg	Price
EM090367	 Yellow	1 ea	\$99.95
EM090367-RB	 Navy	1 ea	\$99.95
EM090367-RED	 Red	1 ea	\$99.95
EM090376-BLK	 Black	1 ea	\$99.95

**MedKit Pro X Medications Module**

- Easy access medication filing system
- Four (4) storage options for first line meds, IV's, pre-filled meds, ampules or boxes
- Constructed from non-absorbing, wipe-clean infection control material
- Lockable zippers for pharmacy seal
- Dimensions: 6.25" L x 10.25" W x 5.7" D



SKU	Description	Pkg	Price
EM090369	 Yellow	1 ea	\$109.95
EM090369-RB	 Royal Blue	1 ea	\$109.95
EM090369-RED	 Red	1 ea	\$109.95
EM090369-BLK	 Black	1 ea	\$109.95



### Iron Duck® IDX Expandable Drug Module

- Holds up to eight (8) removable and customizable pages with multiple pockets and elastic tacking for supplies
- Expandable zipper increases the storage capacity by 175.5 cubic inches
- Front clear vinyl window for easy identification of medications and expiration dates
- Dimensions: 9" L x 6.5" W x 3.5" D



SKU	Description	Pkg	Price
EM090370	Red	1 ea	\$111.88



### Iron Duck® Opioid Overdose Case

- Secure storage for two (2) pre-filled syringes
- Interior mesh pouch for storing multiple MAD devices
- Room for a CPR Barrier Device and a pair of exam gloves
- Dimensions: 7.5" L x 3.5" H x 2" D



SKU	Description	Pkg	Price
EM090120	Yellow; Empty Case	1 ea	\$80.99

### Bag Accessories



### StatPacks G3 Vial Strand

- Holds multi-dose meds in place with yellow pull-tabs for easy removal
- Each strand may hold five (5) vials of various sizes



SKU	Description	Pkg	Price
EM090371	Red	4/pk	\$31.50



### Iron Duck® Ergonomic Shoulder Strap

- Ergonomic padded shoulder strap with speed clip ends
- Can be added to many Iron Duck bags to make carrying more comfortable and convenient



SKU	Description	Pkg	Price
EM090373	Black	1 ea	\$42.50



### MERET X 12™ Divider

- Durable infection control materials with inner foam construction maintains shape
- 1" Velcro® strips provides secure connections
- Sliders can be attached to create multiple organizing options



SKU	Description	Pkg	Price
EM090431	Grey; Standard	1 ea	\$9.95



### MERET X SLIDERS™

- Durable infection control materials with inner foam construction maintains shape
- 1" Velcro® strips provides secure connections
- Sliders interlock to create adjustable dividers from 3" to 5.5"



SKU	Description	Pkg	Price
EM090374	Grey; Infection Control	1 ea	\$7.95



### MERET X STRAP™ Shoulder Strap

- Constructed from infection control materials and a 1.5" webbing strap
- Strap includes all metal hardware

SKU	Description	Pkg	Price
EM090375	Yellow	1 ea	\$19.95
EM090375-RB	Navy	1 ea	\$19.95
EM090375-RED	Red	1 ea	\$19.95
EM090375-BLK	Black	1 ea	\$19.95



EM090375-RED



SimRig® The Ambulance Trainer provides an unparalleled level of realism and precision to EMS training, allowing students to enter the workforce with a complete understanding of how to operate within an ambulance.

FEATURES FOR BOTH CLASSROOM & MOBILE SIMRIG®

Functioning action control station

Flip-Up Bumper built to the specifications of an actual ambulance

Integrated cot mounting apparatus

Functioning suction and simulated oxygen

Over thirty (30) additional standard features

*\*Lifetime structural warranty and five (5) year limited warranty on all other components*



**Established Brand**

DiaMedical is at the forefront of simulation training, providing tools for training and simulation.



**Facilitates Group Learning**

SimRig's® open design and HD camera system facilitates instructor and student observation.



**Unparalleled Realism**

SimRig® is built to precisely mimic a real ambulance. All aspects, including the back door and bumper, are built to actual specifications of a road ready ambulance.



**Scientifically Backed**

Studies have shown that hands-on simulation training can help enhance performance and reduce errors in critical life-saving emergencies.





SRM00220



## SimRig® The Classroom Ambulance Trainer

SRM00220

Teach with an ambulance in your classroom! The unique modular design allows for the SimRig® to be installed in any classroom or training center! The SimRig® Classroom Ambulance Trainer is built to precisely mimic a road-ready ambulance to help prepare the students for the hardships of working within a confined space.

**CALL (877) 593-6011 FOR PRICING**



SRM00261



## SimRig® The Mobile Ambulance Trainer

SRM00261

Simulate the entire call process anywhere! SimRig® The Mobile Ambulance Trainer is designed to be transported behind any trailer compatible vehicle to different training sites and arrives with standard trailer hitch and functioning electrical components. By being completely mobile, instructors can simulate every aspect of an emergency call from a parking lot, shopping center, simulation site or even streetside.

**CALL (877) 593-6011 FOR PRICING**

### PARTS & ACCESSORIES FOR BOTH CLASSROOM & MOBILE SIMRIG®

SKU	Description
SRM00231	HD 4-Color Camera System In Ceiling
SRM00230	46" HD Color Monitor
SRM00237	SimRig Rollok Rear Security Door
SRM00238	SimRig Rollok Side Security Door
SRM00223	Observation Window
SRM00232	Siren Package
MS099232	2-Way Radio
MS099211	Akro Storage Bin Containers
MS099221	SimRig® Disposable Supply Package



**CONTACT YOUR ACCOUNT MANAGER TO CUSTOMIZE YOUR SIMRIG®**

(877) 593-6011 | [SimRig@DiaMedicalUSA.com](mailto:SimRig@DiaMedicalUSA.com)





# SIMRESCUE™

THE AUTO EXTRICATION SYSTEM



## Why Educators Love SimRescue™ The Auto Extrication System

- **Simulates Dozens of Rollover Positions:**

SimRescue™ is constructed with strength-rated steel construction so the system can be safely setup in numerous positions for unique extrication scenarios.

- **Easy Transportation:**

SimRescue™ easily separates in half for transportation and storage. When disassembled, SimRescue™ can fit through any standard size doorway.

## SimRescue™ The Auto Extrication System



Retractable steering wheel with full dash and airbags



Tempered glass breakage simulator, pry box simulator and seat belt cutting simulator



Removable Lexan windows, adjustable seats and foot pedals with foot pinning simulator



SimRescue™ Starter Package:

- Twenty (20) seat belts
- Twenty (20) tempered glass pieces
- Twenty (20) pry plate pieces
- One (1) crowbar
- One (1) window punch
- Tools for assembly





**SimRescue™ Facilitates  
Dozens of Simple &  
Complex Rollover  
Positions.**



## SimRescue™ The Auto Extrication System

EM012818

Emergency Medical Service providers and paramedics provide patient care in a variety of settings from calm to chaotic. However, when it comes to responding to an emergent motor vehicle accident, every second matters. With SimRescue™ The Auto Extrication System, students and professionals can safely practice patient removal scenarios! SimRescue™ is the only auto extrication system on the market that was specifically created for patient extrication. From basic seatbelt cutting simulations to lethal roll-over accidents, SimRescue™ allows EMS professionals to feel confident when responding to an auto extrication call.

**CALL (877) 593-6011 FOR PRICING**

### PARTS & ACCESSORIES

EM090148	10' Replacement Seatbelt
EM090149	Pry Plate
EM090150	5" x 5" Tempered Glass



EM012818



Scan to Watch  
SimRescue™ in Action!

# Custom Student Supply Kits

Choose from **over 500 items** to pack your kit!

STEP

01

Pick Your Bag

STEP

02

Select Your Products

Tote Bag



Oxygen Delivery



Infection Control Supplies



IV Start Supplies



Wound Care Supplies




Duffel Bag



Simulated Medications & IV Fluids




CPR Masks



Suction



ID & Organization



Backpack



Pocket Supplies



Diagnostics



Syringes & Safety Needles



Injection Supplies



\* Contact your rep for a complete list of bag options





STEP

**03**

**Personalize Your Kit**

STEP

**04**

**Contact Your  
DiaMedical Account  
Manager to Design  
Your Custom Kit**

[\(877\) 593-6011](tel:8775936011)

[Info@DiaMedical.com](mailto:Info@DiaMedical.com)

IV Tubing Sets



Medication Administration Supplies



Respiratory Supplies



Print Your Logo



Embroider Your Logo









# TERMS & CONDITIONS

## HOW TO ORDER

DiaMedical accepts orders online or by phone, fax, email or mail. Contact information and mailing address are available at the bottom of this page.

## PAYMENT TERMS

Net 30 day payment terms or terms acceptable to your facility. DiaMedical accepts valid purchase order numbers or payment by credit card (VISA, MasterCard, American Express). Payment must be made in U.S.D.

## PRICES & PRICE CHANGES

Prices listed may not include freight or handling fees and are subject to correction or change without notice. In the case of substantial increases, customers will be notified on an as needed basis. However, we always try to keep our prices as low as possible.

## TAX

DiaMedical collects sales tax on all of its products. If your facility qualifies for a tax exemption you must email a completed tax exemption certificate to [accounting@diamedical.com](mailto:accounting@diamedical.com) prior to your order being invoiced.

## WARRANTY

DiaMedical provides the warranty of the manufacturer or a minimum one (1) year limited warranty on all of its products.

## SHIPPING & HANDLING

Shipping charges may vary depending on distance, value, weight and cargo. All orders are shipped UPS Ground or common carrier unless otherwise specified. Shipping charges will be added to your order total. Should you need to know the shipping charges in advance, please contact DiaMedical for a shipping cost to be added to your original quote. Items may ship separately.

## RETURNS

In most cases we can accept returns and provide credit for all purchases returned within 30 days. In the case where we operate as a distributor, we extend the return policies as set forth by our manufacturer or supplier.

## INTERNATIONAL ORDERS

DiaMedical can ship internationally. Canadian orders will be shipped UPS Standard unless otherwise specified. All other international orders will be shipped by DHL or USPS unless otherwise specified. Please note that orders shipping through the USPS are only traceable through the United States.

## USE OF PRODUCTS

All products purchased from DiaMedical are intended for educational purposes only and to be used under the supervision of a qualified instructor. All instructional products are not intended for human or animal use, injection, application and/or ingestion. Purchaser further agrees that all products will be used for their intended instructional purposes only, in a safe and legal manner consistent with all applicable state and federal laws and Purchaser's own internal health and safety policies.

## PRODUCT SPECIFICATIONS

DiaMedical reserves the right to change product specifications on all items without notice. All efforts are made to display accurate product information, however, due to uncontrollable manufacturing conditions, product specifications and colors are subject to change without notice. Actual product color may differ from catalog product color.

## BID LIST NOTICE

DiaMedical wishes to be placed on any bid list. We assure prompt replies on all bids. Please contact us if our W-9 and/or vendor registration information is needed.

## HOW TO ORDER



**PHONE**  
(877) 593-6011




**EMAIL**  
[Orders@DiaMedical.com](mailto:Orders@DiaMedical.com)



**WEBSITE**  
[DiaMedical.com](http://DiaMedical.com)





Resources Contact Us Quick Quote (877) 593-6011

**DiaMedical** Shop Products Medical Mattresses Service & Parts Who We Serve About Us

Login | Register Your Cart

---

[Contact My Rep](#) [Start Shopping](#)

## Hi There, Welcome Back!

[Not You? Log Out](#)

- Overview
- Account History
- Favorites Lists
- Billing
- Settings
- Cases

**23**  
[View All](#)

**\$4,050**  
[Manage Invoices](#)

**13**  
[View All](#)

**5**  
[View All](#)

**Shipping Address** [Edit](#)

DiaMedical  
31440 Northwestern Hwy #150  
Farmington Hills, MI 48334

**Billing Address** [Edit](#)

DiaMedical  
73140 Northwestern Hwy #150  
Farmington Hills, MI 48334

**National Service Network**  
Delivery & Installation, Service, Rentals & Contracts

Call Service Department  
(877) 593-6011, Option 2

**A Turnkey Solution for Healthcare Management**

DiaMedical USA improves the quality of healthcare by creating customized solutions for medical facilities and educational institutions to streamline operational efficiencies allowing clinicians and educators to prioritize patient care.

[Start Shopping](#)



SCAN TO  
[Create an Account >](#)



Meet Your Account Manager

(877) 593-6011

[Info@DiaMedical.com](mailto:Info@DiaMedical.com)

## Why Create an Online DiaMedical Account?



### Request Quotes

Easily request **new quotes** or **convert a quote** to a sales order with the click of a button!



### Account History

Manage open orders, view past orders, find tracking updates, **simplify re-ordering and return** requests.



### Billing

Access all of your **invoices in one location** and stay up to date on payments and statements.



### Favorites List

Create lists of your favorite products to **stay organized** and easily find the items you need.



### Customization

Update your account profile with **full control** over saved credit cards, default addresses and email preferences.



31440 Northwestern Hwy, Suite 150  
Farmington Hills, MI 48334

# Your Single Source for Advancing Pre-Hospital Care & Realistic Emergency Medical Training

Infection Control | Airway & Oxygen Delivery | Trauma & Wound Care | Patient Handling & Immobilization  
Monitoring & Diagnostics | IV & Drug Delivery | Training & Simulation | Emergency Bags



## Quick Ship Medical Supplies

Over 3,000 Warehoused  
Medical Products



(877) 593-6011 | [DiaMedical.com](http://DiaMedical.com) | [Info@DiaMedical.com](mailto:Info@DiaMedical.com)